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# Sterisart Universal Pump 16420

## Sterisart® Universal-Pumpe 16420

The enclosed CD contains the operating instructions for Sterisart Universal Pump 16420 as a PDF file in French, Italian and Spanish.

If the CD is missing, you can obtain a copy from us by specifying the order number:

Order number: 85032-537-64  
Publication number: SLD6026-I

Sartorius Stedim Biotech GmbH  
Technical Editorial Department  
August-Spindler-Strasse 11  
37079 Goettingen, Germany

[www.sartorius-stedim.com](http://www.sartorius-stedim.com)

System requirements:  
Windows® XP or higher  
Adobe® Reader® 5.0 or higher

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# Contents

Manufacturer   Distributor .....	4
EC Declaration of Conformity .....	4
General Safety Information .....	5
<b>1 Intended Use of the Sterisart Universal Pump .....</b>	<b>7</b>
<b>2 Important Notes on the Safe and Proper Usage of the Sterisart Universal Pump .....</b>	<b>7</b>
<b>3 Unpacking the Pump and Checking the Equipment Supplied .....</b>	<b>8</b>
<b>4 Description of the Sterisart Universal Pump .....</b>	<b>9</b>
<b>5 Startup .....</b>	<b>11</b>
<b>6 Application .....</b>	<b>13</b>
6.1 Operating the Sterisart Universal Pump 16420... ..	14
<b>7 Configuring the Sterisart Pump .....</b>	<b>18</b>
<b>8 Emergency Operation Mode .....</b>	<b>24</b>
<b>9 Maintenance .....</b>	<b>24</b>
<b>10 Cleaning and Sanitization .....</b>	<b>25</b>
<b>11 Installation in an Isolator .....</b>	<b>27</b>
<b>12 Important Notes Concerning the Warranty .....</b>	<b>27</b>
<b>13 Troubleshooting .....</b>	<b>28</b>
<b>14 Technical Specifications .....</b>	<b>29</b>
<b>15 Ordering Information for the Sterisart Universal Pump .....</b>	<b>30</b>
<b>16 Ordering Information for Consumables .....</b>	<b>31</b>
<b>17 Literature .....</b>	<b>32</b>
<b>18 Information and Instructions on Disposal .....</b>	<b>32</b>
Non-hazardous Declaration .....	33

# Manufacturer | Distributor

**Manufacturer:**

**Möller Medical GmbH**

Wasserkuppenstr. 29-31  
D-36043 Fulda, Germany

**Distributor:**

**Sartorius Stedim Biotech GmbH**

August-Spindler-Strasse 11  
D-37079 Goettingen, Germany

Tel.: +49.551.308.0  
Fax. +49.551.308.3289  
[www.sartorius-stedim.com](http://www.sartorius-stedim.com)

## EC Declaration of Conformity

**We hereby declare that the Sterisart Universal pump 16420 has been designed and manufactured according to the following directives:**

- 2006/95/EC – Low voltage
- 2006/42/EC – Machinery
- 2004/108/EC – Electromagnetic compatibility
- 2002/95/EC – RoHS Directive

**The Sterisart Universal pump 16420 is labeled with the CE mark.**

# General Safety Information

## Explanation of Safety Symbols Used

These operating instructions contain important information which is visibly marked. This information must be observed in order to avoid hazards to users as well as to avoid damage and | or malfunctions to the device.

### Symbols in the operating instructions:

**Important Note!** Important information  
**Danger!** General danger with risk of injury

**Laser warning**

**General information and | or help**

---

### Symbols on the manufacturer's ID label

**AC voltage**

**Return and disposal**

**Observe instructions for use**

**Manufacturer**

---

# 1 Intended Use of the Sterisart Universal Pump

The new Sterisart Universal pump is a peristaltic pump equipped with appropriate devices for attaching a closed sterility test system and a sample container.

The pump draws liquid from the sample container into the specially designed double tubing, distributing the liquid equally to both sample containers of the sterility test system. The sample is rinsed in a similar way. After the nutrient media is added, the sample containers are incubated for a prescribed time and evaluated.

If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired!

## 2 Important Notes on the Safe and Proper Usage of the Sterisart Universal Pump

- The Sterisart Universal pump must be switched off before any electrical connections are made or broken (e.g. when plugging in or unplugging the power cable, foot switch or control unit with display).
- The Sterisart Universal pump must only be used with the original current-carrying parts supplied with it (power cable with ferrite core, D-Sub connecting cable and foot switch). The D-Sub protective cap supplied must cover the port when the D-Sub connecting cable is not connected.
- The Sterisart Universal pump is designed for use with Sterisart NF filtration systems. The suitability of any other tubing must be approved by Sartorius.
- Only authorized Sartorius service technicians are allowed to open the pump housing.
- Do not expose the Sterisart Universal pump to extreme temperatures.
- Improper use or handling will invalidate the manufacturer's guarantee.
- Sterisart Universal pump may only be used with a grounded, readily accessible and fuse-protected electrical outlet that complies with national standards.
- When using device components that do not correspond to their original version, performance, safety and EMC behavior can be negatively affected.

## 3 Unpacking the Pump and Checking the Equipment Supplied

When unpacking the Sterisart Universal pump 16420 make sure not to leave any parts in the packaging. The equipment supplied includes the following accessories (1 piece each):

- Sterisart Universal pump 16420
- Control unit
- Clamp for control unit
- Power cable (country-specific)
- Draining container
- Draining container cover (Sterisart units)
- Drain tubing
- Fastening arm
- Bottle holder
- Support rod
- Bow
- Clamping lever with sleeve
- Operating instructions in German | English

## 4 Description of the Sterisart Universal Pump


Fig. 1

The Sterisart Universal pump (Model 16420, see Fig. 1) consists of pump housing (1) with pump head (2), lever (3), bow (4) and handles (5). On the right handle, there is also a bolt with a tilting lock for the draining container (14). The selecting knob (6) is located on the front right. The support rod (7), including bottle holder (8), fastening arm (9) and optional ampoule breaker (10), is located on the left side of the pump. Additionally, the control unit (11) is also fastened to the support rod via the control unit clamp (12) and connected to the pump via the connecting cable plug (13) (see Fig. 1).

**Laser warning**  
The internal barcode scanner emits a laser beam.  
Do not look directly into the beam.

Fig. 2

The display and control unit consists of (Fig. 2): glass display (a), 2 soft keys (b), 1 Esc key (c), window for barcode scanner (d), pump connecting cable with plug (e).

The draining container (Fig. 1; 14) for the sterility test containers is fastened to the right of the pump. Here, a draining tube can be connected to the designated flange connection (Fig. 1; 15).

Fig. 3

Located on the rear panel (see Figure 3): On | Off switch (1), AC power socket (2), socket for the foot switch (3), an RS-232 interface port for connecting a PC/MMC card reader (4), an RS-232 interface port for connecting an external barcode scanner (5), the syringe actuator (6, optional), the RS-485 port (7), a laser warning label (8) and the manufacturer's ID label (9).

## 5 Startup

Set up the Sterisart universal pump ergonomically in a suitable area (cleanroom, laminar flow work bench, isolator).

Before initial start-up, the accessories supplied must be assembled on the Sterisart Universal pump. Please proceed in the following order:

At the back left, push the support rod (Fig. 1, 7) onto the pins of the upper part of the pump until it is firmly in place. The holes in the foot of the support rod can be used to change the position of the support rod.

- First, fasten the bottle holder to the support rod (Fig. 1, 8), next fasten the control unit clamp (Fig. 1, 12) with the display and control unit (Fig. 1, 11; only for version 16420) and the fastening arm (Fig. 1, 9). Then make the electrical connection to the pump via the control unit's cable and plug contact (Fig. 1, 13). If necessary, attach the ampoule breaker (Fig. 1, 10) beforehand.
  - Insert the bow (Fig. 1, 4) and the clamping lever (Fig. 1, 3) on the corresponding bolts on the top of the pump. The lever is jacketed in a plastic sleeve that can be removed for cleaning. If signs of wear appear on this sleeve after it has been in use for some time, it can be re-ordered as a spare part (see Chapter 14).
  - After attaching the drain tubing to the flange connection (Fig. 1, 15) on the draining container (Fig. 1, 14), attach the draining container to the appropriate bolt with tilting lock on the right handle (Fig. 1, 5). The draining container and the cover plate are attached by moving the tilting lock.
- Fasten the end of the drain tubing to the appropriate drain.

### Very Important Note!

**Always make sure that there is a sufficient liquid level differential (> 5 cm) between the draining container and the drain. This makes safe drainage possible and prevents backflow into the draining container. Make sure that the drain tubing is kept lower than the pump to prevent the drain pan from overflowing.**

**Treat the waste as a biohazard. Only dispose of it in accordance with current regulations.**

---

- Plug the optional foot switch into the designated socket located on the rear panel. (see Fig 3, 3). Fasten the plug to the socket by turning it clockwise.
- Plug the power cable into the male power socket on the back of the Sterisart Universal pump and into a grounded electrical outlet (Fig. 3, 2).
- Use the On | Off switch (Fig. 3, 1) at the front of the Sterisart Universal pump to turn it on.

After connecting the power, a boot sequence will appear for 6 seconds with relevant information (Fig. 4). This displays the current software version of the main board (main firmware), the software version of the display and control unit (UI firmware) as well as the device bus address (Fig. 4).

During startup, the control unit carries out a self test whereby the internal barcode scanner is turned on for a brief time.

After a successful switch-on test, the system is factory set to switch into the manual operating mode (see Fig. 6).

Fig. 4

## 6 Application

**Caution!**  
General danger with risk of injury

**Do not operate the Sterisart Universal pump until the tubing has been properly installed. The pump may only be operated with the tubing inserted and only by one person at a time. Never pump any flammable, corrosive or explosive fluids. Users should always wear tight-fitting clothing (no ties, etc.).**

To be able to operate the pump, Sterisart sterility test systems must be installed beforehand. To do this, follow the steps given below:

Remove the sterility test systems from their respective packaging and place the sterility test sample containers in the container for draining on the right side of the pump.

- Now open the lever and the bow.
- Keep the tubing around the rotor of the pump so that no kinks occur and the tubing is not under tension (the arrow on the fixed part surrounding the pump head indicates the later direction of rotation, Fig.5). The tubing is attached by simply inserting it (with it lying on the ground) and fastening the lever (3). In this way, the tubing is automatically correctly positioned.

**Important Note!**  
When threading the tubing, make sure that it is not kinked or twisted.

Fig. 5

- To start the pumping process, the tubing must be fully inserted and the lever must be engaged. A safety function ensures that the pump does not start if the lever or the bow are open.

- After threading the tubing and closing the lever, insert the removal spike | needle of the Sterisart NF system as described in the instructions for the Sterisart NF sterility test system. Then attach the test liquid container (e.g., bag or bottle) by hanging it on the support rod (bag) or by placing it in the holding ring on the support rod (bottle). The holding ring accepts bottles of various sizes. The holder must be turned so that the sterility test system tubing is not under tension after the spike | needle has been inserted and so that the spike | needle cannot slip out of the container.

For testing ampoules, an ampoule breaker is also attached to the support rod. Use the integrated knife to scratch the neck of the ampoule, then snap off the top. The broken-off top of the ampoule falls into the collection container (which can be emptied later).

- For the following pumping process, make sure that the Sterisart NF sample containers are correctly positioned in the intended holders, and that the tubing has been inserted without tension or kinking.

**Important Note!**

**When the lever is opened, the Sterisart Universal pump is turned off immediately by a safety function. In this event, the Sterisart Universal pump can be restarted after adequately closing all components.**

---

## 6.1 Operating the Sterisart Universal Pump 16420

All operating parameters of the Sterisart Universal pump 16420 are saved in programs. Once you have switched on the pump, it will enter the first program saved (default = "Manual"). Press the right-hand soft key Menü (German for "menu") to access the list of programs (program menu).

The following programs are predefined:

1. **Manual**  
see section 6.1.1
2. **Manual Record**  
see section 6.1.2

### 6.1.1 Manual Program

We will use the Manual program here as an example to explain all the information that can be displayed (Fig. 6):


Fig. 6

The following information is shown in the display (Fig. 6):

- 1: Operating mode (Manual)
- 2: Current program step (here: Prg. 1/1)
- 3: Date (dd:mm:yy)
- 4: Time 24h
- 5: Speed, factory-set
- 6: Level of memory capacity already used (Mem: indicates the amount of memory used in %. Please note: 50 data records (see 6.1.2) take up approx. 1% of the internal memory capacity).
- 7: Pump time (current) (shows how long the pump has been running)
- 8: Preprogrammed pump time (in the Manual operating mode, this setting is 0; i.e., inactive so that this information is not shown on the display .
- 9: Soft key: current function (in this case: Speed). The white background indicates that the value (Speed) can be changed (by turning the selecting knob). If the soft key does not have a white background (Esc key pressed), the value cannot be changed. Pressing the soft key (Speed) again returns it to a white background. Pressing the soft key with a white background (current function: Speed) switches the display to Duration. Here the preprogrammed pump time can be changed (turning the selecting knob). If the soft key does not have a white background (Esc key pressed), the value cannot be changed. Pressing the soft key (Duration) again, returns it to a white background. Here you can return to the Speed display by re-pressing the soft key (Duration).
- Soft key: current function (in this case: Speed). The black background indicates that you can change this value (Speed) (by turning the selecting knob). If the soft key does not have a black background (Esc key pressed), the value displayed cannot be changed. Pressing the soft key (Speed) again returns it to black background. When you press the soft key with a black background (current function: Speed), the display will change to Duration. Here you can change the preprogrammed time by turning the selecting knob. If the soft key does not have a black background (Esc key pressed), the value cannot be changed. Pressing the soft key (Duration) again, returns it to a black background. Here you can press the soft key (Duration) one more time to return to the Speed display.
- 10: Soft key: Menu (for selecting the individual programs: Program Menu)

### Program Sequence

The pump is in the basic operating state (Manual program). To start the pumping process, press the selecting knob or the foot switch. To stop the pump, press the selecting knob or the foot switch again. In this case, you can change the pump rate (Speed) before and during a pumping procedure in a range between 10 – 100 (infinitely; does not correspond to revolutions per minute).

After re-pressing the left Speed soft key, the Duration soft key is available. A specific pump time can be entered via Duration.

**Note:** When you enter a specific pump time in Duration, the current pump time is set to zero.

The Esc key has two different functions, depending on the state of the program:

- The left soft key (Speed | Duration) has a white background (active). Pressing the Esc key deactivates the current soft key function (Speed | Duration) so that the respective value (Speed | Duration) can no longer be set (no white background).
- The left soft key (Speed | Duration) does not have a white background (inactive). By pressing the Esc key, the current pump time is reset.

To switch off the pump, use the On | Off switch on the back panel of the pump housing.

### 6.1.2 Manual Record Program (Manual Operating Mode with Data Collection)

For preparing the pump for initial start-up, see Chapter 6.

This mode is reached by pressing the Menu key once. The Program Menu appears on the display and lists the following options (Fig. 7):

- Manual
- Manual Record

The different operating modes are selected by turning the selecting knob. For final selection of an operating mode, press the selecting knob (in this case: Manual Record).

The difference between the Manual program and this Manual Record program is that when the latter begins, a data record can be collected in the manual mode.

The display will now prompt a barcode scan (Scan Operator, Fig. 8).

The control unit now waits until either data recording is started by pressing the left Scan soft key (barcode is scanned via the internal barcode scanner), the current data set recording (current scan process) is bypassed by pressing the Esc key or the program is stopped and returned to the Program menu by pressing the Menu key.

**Note:** The Scan soft key does not have to be pressed for an external scanner. The control unit automatically reads the data from the external scanner and assigns it to the appropriate data records. In the example in Fig. 8, this data is assigned to the "Operator" parameter. The control unit automatically switches to the next parameter as soon as the data have been scanned.


Fig. 7


Fig. 8

Fig. 9

To read in a barcode, hold it under the laser beam within 10 seconds after you have pressed the Scan key (see Fig. 9). If no barcode is recognized or read in during this time, press the Scan key again to repeat this procedure. A successful scanning process is confirmed by an acoustic signal. At the same time, the barcode recognized is shown on the display for approx. 1 second.

Follow the same procedure to scan in the 5 additional barcodes in the order given below:

- Scan sterility test unit
- Scan sample
- Scan rinsing fluid
- Scan tryptic soy (TSB)
- Scan thioglycollate (FTM)

You can press the Esc key at any time during scanning to skip an item.

After scanning, the pump will be in the basic operating state as in Manual Program (see section 6.1.1)

At the end of the sterility test, the data records collected can be output using a PC or an optionally available MMC/SD card reader.

For more on the topic of data evaluation, see separate instructions for use.

## 7 Configuring the Sterisart Pump

Press the Menu soft key twice to go to the pump's configuration menu (Fig. 10). Here, the following parameters can be set.

If you do not press any key or switch within 60 seconds (default can be changed in menu item "Operation Timeout Menu"), you will automatically exit from the system menu and the pump will return to the last program step it was performing.

In general, a changed value is first saved by pressing the selecting knob. You can also exit the System menu by pressing the Esc key. In this way, the set value is not applied.

### Storage Card (Memory Card)

To access this menu, you will first need to connect an external card reader and insert an MMC/SD. Interface configuration is not required.

#### Save data to card (Fig. 11):

Default setting: No  
Selectable options: No | Yes

By confirming this menu prompt (selectable option: Yes), all data records logged up to this point will be transferred to the MMC/SD card in the card reader and deleted from the pump memory to provide free memory space. Active data transfer is indicated by a series of bars. While this is being displayed, do NOT remove the card.

**Note: If there is no MMC/SD card in the card reader, the "File error" error message is displayed.**

#### Read programs (Fig. 12):

Default setting: No  
Selection options: No | Yes

By pressing Okay to confirm this prompt (selectable option: YES), the program structure currently on the storage card will be transferred to the pump. Active data transfer is indicated by a series of bars. While this is being displayed, do NOT remove the card.

**Here, it is important that the current program be in the root directory (main directory) of the card and exactly matches the file name "Programm.smm."**


Fig. 10


Fig. 11


Fig. 12


Fig. 13

### Password

Inputting a password (Fig. 13):

Default setting: 0  
 Minimum value: 0  
 Maximum value: 29999

The password can be a 5-digit number and it is entered using the selecting knob. Password entry must be confirmed by re-pressing the selecting knob.

Entries visible for the user in the System menu are released by entering two freely-selectable passwords. These levels are then accessible to all users without the further need of entering a password. The following relationship exists between the respective password and the entries in the System menu.

#### Password 1:

Password factory setting 101  
 Levels: Memory card, password, Programs, display, Operation, audio signals

#### Password 2:

Password factory setting 202  
 Levels: Memory card, password, Programs, display, Operation, audio signals, Date | Time, data recording, interfaces

**Note: If there is no agreement with the saved password and the password just entered, the system automatically returns to the System menu.**

After entering the password, the user can change the respective default password. The following relationship exists between the password entered at the beginning and the changes that can be made to the default passwords:


Fig. 14

### Changing password 1 (Fig. 14):

Default setting: 101  
 Minimum value: 0  
 Maximum value: 29999

The 101 default password can be changed in this menu.

### Changing password 2:

Default setting: 202  
 Minimum value: 0  
 Maximum value: 29999

The 202 default password can be changed in this menu.

**Note: When entering password 2, the user also has the option of changing password 1.**


Fig. 15

### Programs (Fig. 15):

Using this menu item, adjustments can be made to the program step that is currently activated. In detail, these are:

#### Speed Minimum

Default setting: 10  
 Minimum value: 10  
 Maximum value: 100

This parameter displays the lower limit value of the speed of the current program step down to which the later user may change the speed.

#### Standard Speed

Default setting: Current speed  
 Minimum value: Speed Minimum  
 Maximum value: Speed Maximum

This value displays the speed of the current program step. The operator can change the speed within the limit values Speed Minimum and Speed Maximum.

#### Speed Maximum

Default setting: 100  
 Minimum value: 10  
 Maximum value: 100

This parameter displays the upper limit value of the speed of the current program step up to which the later user may change the speed

#### Speed Steps (Increments)

Default setting: 10  
 Minimum value: 1  
 Maximum value: 10

This value represents the amount the current speed is increased or reduced by turning the wheel by one position.

#### Speed Adjustment Timeout [s]

Default setting: 0 s (Speed soft key always with white background)  
 Minimum value: 0 s (Speed soft key always with white background)  
 Maximum value: 120 s

The Speed soft key loses the white background (becomes inactive) if the key is not pressed before this time expires. The soft key would then have to be pressed again to make a new change to the Speed value.

The value 0 means that the time period does not end, i.e. the soft key no longer deactivates automatically.

#### Duration Minimum [s]

Default setting: 0 s (infinite pump time)  
 Minimum value: 0 s (infinite pump time)  
 Maximum value: 1800 s

This parameter represents the lower limit value of the current program step's duration down to which the later user may change the duration

### **Duration [s]**

Default setting: Current duration  
Minimum value: Duration minimum  
Maximum value: Duration maximum

This value displays the duration of the current program step. The operator can change any value between Duration Minimum and Duration Maximum.

### **Duration Maximum [s]**

Default setting: 300 s  
Minimum value: 0 s  
Maximum value: 1800 s

This parameter displays the upper limit value for the pump time up to which the later user may change the duration.

### **Increment Duration**

Default setting: 10 s  
Minimum value: 1 s  
Maximum value: 300 s

This value represents the amount that the current duration will be increased or reduced by turning the selecting knob around one position.

### **Duration Adjustment Timeout [s]**

Default setting: 0 s (Duration soft key always with white background)  
Minimum value: 0 s (Duration soft key always with white background)  
Maximum value: 120 s

The Duration soft key loses the white background (becomes inactive) if the key is not pressed before this time expires. The soft key would then have to be pressed again to make a new change to the Duration value.

The value 0 means that the time period does not end, i.e. the soft key no longer deactivates automatically.

### **Behavior of the foot switch**

Standard setting: Tactile  
Minimum: Switch  
Maximum: Tactile

Use this setting to adapt the behavior of the foot switch. In the "switch" setting, you can start the rotor and put it on hold by operating the foot switch. In the "tactile" setting, the rotor will keep turning as long as the foot switch is pressed. Releasing the foot switch stops the rotor.


Fig. 16

### **Display (Fig. 16)**

#### **Contrast Display**

Default setting: 50  
Minimum value: 10  
Maximum value: 100

Use this parameter to set contrast of the display.

### Contrast [%]

Default setting: 50  
Minimum value: 10  
Maximum value: 100

This parameter is used to set the display contrast.

### Display Time

Default setting: Yes  
Selectable options: Yes | No

Use this parameter to define whether the time is shown or not shown when the Main Menu is displayed.

### Operation (Fig. 17)

#### Language

Default setting: English  
Selectable options: German | English

The current language can be changed between English and German using this menu item.

#### Timeout Menu [s]

Default setting: 60 s  
Minimum value: Off (current menu is not exited)  
Maximum value: 600 s

After a certain time has elapsed, the Setup menu closes automatically, if no key is pressed within this time.

#### Summertime Automatic (Automatic Daylight Savings Time)

Default setting: Yes  
Selectable options: Yes | No

Use this parameter to activate or de-activate automatic switching of the clock to daylight savings time ("summer time") or to winter time.

### Audio Signals (Fig. 18)

#### Switching on the Audio Signal

Default setting: 2  
Minimum value: Off  
Maximum value: 4

This parameter can be used to adjust the volume of the confirmation signal between OUT and level 4 by pressing a key.

#### Signal for program end

Default setting: 2  
Minimum value: OFF  
Maximum value: 4

This parameter is used to set the volume of the sound signaling the program end (range can be set between off to level 4).


Fig. 17


Fig. 18


Fig. 19

### Date | Time (Fig. 19)

#### Year

Default setting: Current year  
 Minimum value: 2010  
 Maximum value: 2099

In this menu item, you can set the current year:

#### Month

Default setting: Current month  
 Minimum value: 1  
 Maximum value: 12

Use this menu item to set the current month.

#### Day

Default setting: Current day  
 Minimum value: 1  
 Maximum value: 31

Use this menu item to set the current day.

#### Hour

Default setting: Current hour  
 Minimum value: 0  
 Maximum value: 23

Use this menu item to set the current hour.

#### Minutes

Default setting: Current minute(s)  
 Minimum value: 0  
 Maximum value: 59

Use this menu item to set the current minutes.

#### Seconds

Default setting: Current second(s)  
 Minimum value: 0  
 Maximum value: 59

Use this menu item to set the current seconds.


Fig. 20

### Data format (Fig. 20)

#### Data recording

Default setting: Yes  
 Selectable options: Yes | No

Using this parameter, data recording via the barcode scanner can basically be turned on and off. This change affects all programs in equal measure (Manual Record).


Fig. 21

## Interfaces (Fig. 21)

### Barcode Scanner

Default setting: Internal  
 Selectable options: External | Internal

This parameter is used to select between the internal and the optional external barcode scanner.

### Baud rate PC RS-232

Default setting: 38.400  
 Values: 9600 | 19200 | 38400

Use this parameter to set the baud rate for communication with a computer (PC).

## 8 Emergency Operation Mode

The pump is equipped with an emergency operation mode (see Fig. 22) that enables the user to continue to operate the pump in case the display and control unit malfunctions. However, the user must make sure that the supplied protective cover is attached to the connector plug.

In emergency operation mode, the pump can only be operated at two different speeds (speed 20 and 50). The emergency operation mode starts automatically when the electrical connection between the display and control unit and the pump is interrupted. The start of the emergency operation mode is signaled by the switch-on signal. If the display and control unit is connected in this mode, the emergency operation mode is displayed as per Fig. 22. The user also has the option of restarting the pump via the soft key to return to normal operation.

### Speed

Default setting: 20  
 Values: 20 | 50

The pumping process is started by pressing the selecting knob | foot switch. The current pumping process is stopped by re-pressing the selecting knob | foot switch. The speed can only be changed during the pumping process. Turning the selecting knob clockwise changes the speed from 20 to 50 (counterclockwise from 50 to 20).

## 9 Maintenance

We recommend that you have a trained and qualified Sartorius service technician perform annual maintenance on your Sterisart Universal-Pump. However, your Sterisart Universal-Pump may need maintenance at more frequent intervals, depending on its particular use and on your specific process conditions. If you have any further questions, please contact your local representative who will be glad to assist you.

Fig. 22

# 10 Cleaning and Sanitization

**Caution!**  
Ports on the back panel are not tight.

The Sterisart Universal pump is not suited for sanitization (gaseous or atomized) with peracetic acid, hydrogen peroxide, or other decontamination procedures commonly used in isolators, unless it is built into an isolator countertop. These types of disinfectant can be used when the pump has been professionally installed in the isolator.

Make sure that no liquids penetrate the socket of the external display and control unit. If the external display and control unit is unplugged during the cleaning process, the supplied protection cover must be inserted into the socket.

- Before cleaning the pump rotor, turn off the main power switch and unplug the power cable.

The following parts of the pump have to be removed for cleaning:

- Bow, lever, support rod and draining container
- Disconnect the pump head:
  - The rotor must be turned by hand until the marking on the upper side (black dot) is pointing backwards.
  - The marking indicates the position of the fastening screw that is located in the rotor mid-section.
  - Loosen this screw with a 4-mm Allen wrench. Do not turn the screw out completely.
  - Pull the rotor upwards.
- Pull the speed-selecting knob upwards.

## General Instructions:

**Comply with the IP protection rating! No liquid may enter the Sterisart Universal pump. Disconnect the pump from line voltage before any cleaning, maintenance, or repair work.**

## Cleaning stainless steel surfaces:

The stainless steel housing can be cleaned at regular intervals. Clean stainless steel parts with a wet cloth. Customary cleaning agents suitable for stainless steel can be used without risk. For surface sanitization, use commercially available disinfectants that are suitable for stainless steel.

Do not use any cleaning agents on the stainless steel housing that contain sodium hydroxide or acetic, sulfuric or citric acid. The use of steel wool cleaning sponges is prohibited.

All stainless steel surfaces can be cleaned by wiping with a moist cloth and with the following products:

- Peracetic acid, H<sub>2</sub>O<sub>2</sub> and acetic acid (2.5 %)
- Quaternary ammonium
- Formaldehyde and glutaraldehyde
- Isopropyl alcohol
- Ethanol 70%

Autoclavable at 121°C for 20 min. are bow, lever, speed selecting knob, support rod, bottle holder, counter-bottle holder and draining container (complete).

**Never autoclave the pump rotor and never submerge it into any liquid.  
Only use a slightly moistened cloth | tissue for cleaning | disinfecting.**

---

After cleaning, all parts should be secured as described in Chapter „Startup“ (5).

The rotor is assembled as follows:

- The rotor shaft has a groove on the side turned towards the back.
- Push the rotor upwards onto the shaft.
- Turn the marking on the rotor towards the back, so that the groove on the rotor shaft aligns with the marking on the rotor
- Tighten the screw in the rotor to guarantee that the rotor is firmly seated.

# 11 Installation in an Isolator

**The installation must be carried out by the isolator manufacturer.**

To install the Sterisart Universal pump in an isolator, it must be opened up first. Make sure to unplug it from AC current (mains supply) before opening it!

During assembly, please follow the installation instructions precisely (part of the installation kit for isolators, order number 1ZE---0050). The instructions are available separately from the manufacturer.

# 12 Important Notes Concerning the Warranty

The correct functioning and safe use of the Sterisart Universal pump is only guaranteed when it is operated according to the information given in these instructions.

The equipment manufacturer is responsible only for effects on safety, dependability and performance of the instrument if assembly, extensions, new settings, changes or repairs are carried out by persons authorized by him to do so, if the electrical installation in the room is in conformance with IEC requirements and if the equipment is used in accordance with these operating instructions.

Any changes made to the Sterisart Universal pump will invalidate the manufacturer's warranty unless they are previously agreed to by authorized Sartorius Service staff, or carried out by them.

## 13 Troubleshooting

This section lists a few problems that could arise in connection with the Sterisart pump. Several solutions are given for each problem. The first solution given is usually the most likely one. If this does not solve the problem, go down the list of suggestions until the problem has been solved.

**The pump must be switched off before electrical connections are unplugged or plugged in!**

Should it not be possible to locate and solve the problem in this way, contact the appropriate service center (see the back cover of this manual).

Problem	Solution
No function; no display	<p>The pump is not switched on, or is not correctly connected to the electric supply. Plug the power cable into the electric outlet correctly and to the Sterisart Universal pump.</p> <p>Switch on the on-off switch.</p> <p>Check the electric supply; switch on multiple outlet sockets; check the wiring.</p> <p>Check the plug connection of the display</p> <p>If these possibilities are not feasible, the primary fuse in the main line filter may be defective.</p> <p>Changing this, with the power cable unplugged, must only be performed by authorized personnel!</p>
Emergency operation mode activated	<p>The display was disconnected during the pumping process.</p> <p>Check the plug connection of the display.</p> <p>You may have to replace the complete display and control unit.</p>
The pump rotor does not rotate.	<p>The locking lever is open or the foot switch is not properly connected.</p> <p>The bow has been removed from the pump.</p> <p>Check the connection of the foot switch to the pump.</p>
Data transmission to a PC (using the RS-232 cable) does not work.	<p>The serial cable is not plugged into the correct pump connection.</p> <p>Attach the cable to the upper RS-232 port.</p> <p>Check whether the right COM Port has been selected.</p> <p>The interface settings between the pump and the PC are not compatible.</p> <p>Compare the interface settings in the operating system with the interface settings in the configuration menu</p> <p>The default settings of the interface port are: Baud rate: 38400 Handshake: None Data bits: 8 Stop bits: 1</p>
Barcodes are not recognized	<p>Check whether System menu   Communication for barcode scanners is set to Internal. Otherwise, the controls can address an optionally available external barcode scanner. Another cause might be an unknown barcode format. Default settings support the following barcodes:</p> <p>CODE 39 CODABAR CODE128</p>

# 14 Technical Specifications

Dimensions pump in mm:	(W × D × H) approx. 336 × 260 × 210 (with lever)
Dimensions pump with bottle holder, draining container in mm:	(W × D × H) 440 × 365 × 485
Weight	
Version 16420:	Approx. 14.6 kg
Maximum operating temperature:	+10 – +40° C

## Electrical connection

Voltage:	100 – 240 VAC (alternating voltage)
Frequency:	50 – 60 Hz
Power consumption:	100 VA

## Pump

Pump flow rate:	70 – 650 ml/min (depending on the tubing)
-----------------	---

## Tubing

Correct functioning of the Sterisart Universal pump is only guaranteed when used with tubing that has been approved by Sartorius.

## Materials

Base plate, lever, bow, rotor:	Stainless steel 316L (1.4435)
Housing, connection port, housing control unit, support rod, bottle holder:	Stainless steel 316 (1.4401)
Connecting cable:	PVC
Connector grommet:	TPE
Display connecting plug:	Mat chrome
Draining container, plastic sleeve for lever, hydrodynamic bearing for bow:	PET
Hydrodynamic bearing for lever:	Iglidur X
Seals, gaskets, drain tubing:	Silicone

## Environmental Conditions

Protection class:	ⓘ
Humidity:	less than 80% relative humidity
Air pressure:	700 – 1050 hPa
Operating altitude:	max. 2000 m MSL
Overvoltage category:	Ⓜ
Pollution degree:	2"

## 15 Ordering Information for the Sterisart Universal Pump

16420	Sterisart Universal pump, extended version with display   upgraded user software
16973	Ampoule tongs
1ZG---0009	Bottle holder
1ZF---0004	Fastening arm
1ZA---0002	Drain tubing
1ZW---0002	Ampoule breaker
1ZE---0033	Foot switch
1ZF---0013	Support rod
1ZG---0022	Draining container
1ZB---0001	Plastic sleeve for lever
1ZF---0014	Clamp for control unit
1ZF---0007	Wing nut
1ZS---0005	Bow
1ZS---0002	Lever
1ZL---0003	Rotor
1ZE---0034	RJ45-cable
1ZE---0035	Serial cable
1ZG---0023	Draining container cover (Sterisart sterility test units)
1ZG---0024	Draining container cover (Millipore sterility test units)
1ZG---0014	Adapter for Sterisart draining container for use with a Millipore Equinox pump
1ZE---0039	Transport trolley
1ZE---0036	MMC/SD card reader for 16420
1ZE---0037	MMC card for 16420
1ZE---0040	Communication kit
1ZE---0050	Installation kit for isolators

Additional accessories are available on request.

## 16 Ordering Information for Consumables

Recommended disposable sterility test units for use with pump

Order No.	Description
16466-ACD	Sterisart NF alpha, dual-needle metal spike for closed containers (box of 10, individually sterilized with ETO; single-packed).
16467-ACD	Sterisart NF alpha, metal needle for open containers (box of 10, individually sterilized with ETO; single-packed)
16468-ACD	Sterisart NF alpha, system for medical devices with Luer lock connectors (box of 10, individually sterilized with ETO; single-packed)
16466-GBD	Sterisart NF gamma, dual-needle metal spike for closed containers (box of 10, individually gamma sterilized, double-packed). Optimal for use in isolators.
16467-GBD	Sterisart NF gamma, metal needle for open containers (box of 10, individually gamma sterilized, double-packed). Optimal for use in isolators.
16468-GBD	Sterisart NF gamma, system for medical devices with Luer or Luer lock connectors (box of 10, individually gamma sterilized, double-packed). Optimal for use in isolators.
16469-GBD	Sterisart NF gamma, system with adapter for prefilled syringes (box of 10, individually gamma sterilized, double-packed). Optimal for use in isolators.
16470-GBD	Sterisart NF gamma, system for slightly soluble powders in unvented vials (box of 10, individually gamma sterilized, double-packed). Optimal for use in isolators.
16475-GBD	Sterisart NF gamma, system for lyophilized or freely soluble powders in unvented vials (box of 10, individually gamma sterilized, double-packed). Optimal for use in isolators.
16476-GBD	Sterisart NF gamma, system with short dual-needle metal spike for closed containers (box of 10, individually gamma sterilized, double-packed). Optimal for use in isolators.
16477-GBD	Sterisart NF gamma, system for liquids in plastic containers (box of 10, individually gamma sterilized, double-packed). Optimal for use in isolators.
16478-GBD	Sterisart NF gamma, adapter for medical containers with Luer Lock male connectors (box of 10, sterilized, double-packed). Optimal for use in isolators.

Additional Sterisart units are available on request.

## 17 Literature

Data sheet Sterisart Universal Pump      English  
Publ. No.: SLD2010-e      Order No.: 85030531-17

Brochure Sterisart Universal Pump      English  
Publ. No.: SLD1003-e      Order No.: 85030-531-14

## 18 Information and Instructions on Disposal

Attention must be paid to any applicable laws, rules and regulations, when disposing of contaminated components.

The packaging is to be taken to a local waste disposal site if no longer required. The packaging is made of environmentally-friendly materials that can be used as secondary raw materials.

The equipment, including accessories and batteries, does not belong in your regular household waste. EU legislation requires its Member States to collect electrical and electronic equipment and dispose of it separately from other unsorted municipal waste with the aim of recycling it.

In Germany and several other countries, Sartorius Stedim Biotech GmbH itself assumes responsibility for the return and conformant disposal of its electronic and electrical products. These products may not be placed with household waste or brought to collection centers run by local public disposal operations – not even by small commercial operators.

For disposal in Germany and in the other member nations of the European Economic Area (EEA), please contact our local service technicians or our Service Center in Goettingen, Germany:

Sartorius Stedim Biotech GmbH  
Servicezentrum Biotechnologie  
August-Spindler-Str. 11  
37079 Goettingen, Germany

WEEE-Reg.- No. DE 89907997

In countries that are not members of the European Economic Area (EEA) or where no Sartorius Stedim subsidiaries or dealerships are located, please contact your local authorities or a commercial disposal operator.

Prior to disposal and/or scrapping of the equipment, any batteries should be removed and disposed of in local collection boxes.

Sartorius Stedim Biotech will not take back equipment contaminated with hazardous materials (ABC contamination) – either for repair or disposal. Please refer to the package leaflet | the included CD ROM or to our web site ([www.sartorius-stedim.com](http://www.sartorius-stedim.com)) for more detailed information and addresses for repair service or disposal of your device.

### Non-hazardous Declaration

Declaration of Safety and Acceptability as Non-Hazardous to Health

The safety and health of our employees, the legal regulations governing the handling of materials that pose hazards to health and the environment and the regulations governing occupational safety mandate and require that this Non-hazardous Declaration be filled out completely for all equipment returned to the respective Sartorius affiliate, subsidiary, dealer or distributor listed for your country. No repair or return is possible unless a completed Non-hazardous Declaration form has been submitted.

Fax or mail a copy of a completed Non-hazardous Declaration form in advance to the respective Sartorius affiliate, subsidiary, dealer or distributor listed for your country (for the fax number or address, please refer to our list of countries). We require this information before the equipment | component arrives. An additional copy must accompany the equipment | component. If necessary, the shipping company must be notified.

Incomplete entries or non-compliance with this procedure will automatically lead to considerable delays in processing.

In the event of return, we do not differentiate between

- still utilizable and | or repairable used electrical and electronic equipment and no longer utilizable and | or repairable used electrical and electronic equipment or
- between waste that is subject to monitoring by the supervisory authorities and waste that is not subject to any special monitoring by supervisory authorities (used electrical and electronic equipment).

Model:	Serial No.:
<p>I   we guarantee that the following measures have been taken:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> The equipment   component has been freed of hazardous materials to ensure that persons involved with the handling   repair are not exposed to any hazard or danger whatsoever.</li> <li><input type="checkbox"/> The equipment has been packaged for safety and the complete labeling has been affixed to it.</li> <li><input type="checkbox"/> The shipping agent has been notified (as prescribed) about the hazardous nature of the shipment.</li> </ul>	
<p>The person sending in the equipment   component is hereby notified that said person shall be liable to Sartorius Stedim Biotech GmbH, its affiliates, subsidiaries, dealers and distributors as well as to any third parties – in particular to any of employees of Sartorius Stedim Biotech GmbH, its affiliates, its subsidiaries, dealers or distributors involved with the handling   repair of the equipment   component – for any damage, caused by entries that are incomplete or incorrect on account of negligence, gross negligence or willful intent.</p>	
<p>Name: _____</p> <p>Position: _____</p> <p>Date: _____</p> <p>Signature: _____</p> <p>Company stamp:</p>	

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Phone: +49.551.308.0  
www.sartorius.com

The information and figures contained in these instructions correspond to the version date specified below.

Sartorius reserves the right to make changes to the technology, features, specifications and design of the equipment without notice.

Masculine or feminine forms are used to facilitate legibility in these instructions and always simultaneously denote the other gender as well.

**Copyright notice:**

This instruction manual, including all of its components, is protected by copyright. Any use beyond the limits of the copyright law is not permitted without our approval. This applies in particular to reprinting, translation and editing irrespective of the type of media used.

© Sartorius Germany

Last updated:  
03 | 2018

**Software Version:**

--

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Sartorius Lab Instruments GmbH & Co. KG