Virosart® HC Minisart®

Directions for use

1. Introduction

The purpose of this document is to provide the appropriate method for identifying, installing and using Sartorius Stedim Biotech Virosart® HC Minisart®. The procedure should be strictly followed to ensure operator safety and to maintain optimum performance of Virosart® HC Minisart®.

This document pertains to Sartorius Stedim Biotech Virosart® HC filter elements with the following part number:

Virosart® HC Minisart®:

- 539VM-----B (box of 5)
- 539VM-----ITA (box of 4)

Virosart® HC Minisart® has a female luer lock inlet & male luer slip outlet connection.

2. Labeling of Virosart® HC Minisart®

The device consists of a hydrophilic surface modified polyethersulphone (PES) double layer membrane with 5 cm² effective filtration surface within a polycarbonate housing.

The type, lot number and serial number are imprinted on the outer support of the housing of the Virosart® HC Minisart®.


The Virosart® HC Minisart® tested for integrity by diffusion testing in house at Sartorius Stedim Biotech are single packed into a blister package which is labeled with type, order code, lot number and use before date. In addition, these specifications are printed on the labels of the cardboard box.

3. Installing of Virosart® HC Minisart®

△ The Virosart® HC Minisart® tested for integrity in house at SSB (539VM-----ITA) is sterilized in a validated production step by ETO sterilization according to EN ISO 11135. Sterility, however, is only guaranteed when this primary packaging of the Minisart® is intact before using the Virosart® HC Minisart® as a sterile product. If the packaging is damaged in any way the Virosart® HC Minisart® cannot be handled as a sterile product!

If the Virosart® HC Minisart® is to be used as a sterile product, care must be taken to prevent contamination when connecting the Minisart® to an appropriate filtration system (e.g. connection within a clean bench).

The Virosart® HC Minisart® must be connected in the direction of filtration at the female luer lock (printed side) to an appropriate filtration set-up (see picture below). An application specialist from Sartorius Stedim Biotech can provide assistance in installing the system at your facility to make sure that requirements for a correct scale down trial are met. Virosart® HC Minisart® can be connected to the liquid reservoir by using a sterilized disposable 2-way valve. The order code for such valves is 7ZV--K0112, the order code for the stainless steel liquid reservoir is 16249-----6.

Picture 1: Typical filtration set-up for a Virosart® HC Minisart®

4. Wetting of Virosart® HC Minisart® (In the direction of filtration)

Before operation the Virosart® HC Minisart® as well as the entire system needs to be rinsed with buffer or WFI to ensure that it is particle-free as well as for complete wetting of the membrane and proper use of the whole filter capacity. The Virosart® HC Minisart® must be wetted with buffer or WFI prior use. Please connect the Minisart® to the filtration system. Sartorius Stedim Biotech strongly recommends using constant pressure mode for operating the Virosart® HC Minisart®. Carefully increase the inlet pressure to 2 bar | 29 psi and flush the filter approx. 15 minutes with buffer or WFI. Please make sure that no air is entrapped within the Minisart®.

Air bubbles within the Minisart® may decrease the effective filtration area of the device. Typical WFI flow rates may vary between 0.2 – 0.6 ml/min at 2 bar pressure. If the Virosart® HC Minisart® has been flushed according to the instructions it is thoroughly wetted. Now carefully remove the buffer or WFI from the reservoir and replace it with product. Please ensure there is no air trapped upstream of the Minisart® and start the filtration at the desired filtration pressure.

5. Autoclaving Virosart® HC Minisart®

Sterilization temperature: up to 121 °C, 30 min.

**The Virosart® HC Minisart® can only be autoclaved dry!
Please do NOT wet the Minisart® prior to autoclaving.**

Please use a liquid cycle for autoclaving to minimize pressure peaks. Once the sterilization temperature has been attained, the Minisart® can be autoclaved for 30 minutes at a maximum sterilization temperature of 121 °C, 1 bar | 14.5 psi. The inlet and outlet of the Virosart® HC Minisart® should be wrapped with steam permeable autoclave paper or aluminum foil. It is important not to close the Minisart® with plugs (or a pre-filter). During the autoclave cycle, please make sure that any connecting piece on the Minisart® (e.g. made of stainless steel) are not attached too tightly and do not have a tractive effect on the Minisart® housing which could cause deformation.

The cooling phase has to be conducted without deliberately decreasing the pressure (opening the autoclave at temperatures above 85°C), because pressure drops at elevated temperature may lead to bumping and device damages.

Autoclaved Minisart® must be flushed with buffer for 15 min @ 2 bar | 29 psi and used for filtration within 4 h after cooling down to room temperature. With longer storage times after autoclaving there is a possibility of trapping air in the housing. Air in the housing would lead to less effective filtration area and with that misleading filtration results.

Virosart® HC Minisart® can be sterilized by autoclaving up to 2 cycles.

6. Maximal Allowable Differential Pressure

Virosart® HC Minisart®

In the direction of filtration:
At 20 °C max. 5.0 bar | 73 psi

7. Filtration with Virosart® HC Minisart®

Virus filtration with Virosart® HC should be realized at the end of the purification process prior to the last buffer exchange | form and fill. The protein solution which is subject to the virus filtration must be as much as possible free of aggregates, host cell protein & DNA to avoid premature blockage of Virosart® HC filter. Even if these contaminants should be removed either during the initial recovery process or latest during the polishing process of your target molecule, small amounts might be sufficient to cause premature blockage of the final virus filter. To prevent this, an efficient adsorptive pre-filtration step, such as the Virosart® Max, might be required as protection for the Virosart® HC. The optimum pre-filter – final filter ratio has to be identified during development of this process step. It is recommended to use fresh and representative material for these studies.

To make full use of the filtration area, make sure that no air is entrapped within the Virosart® HC Minisart®. The pre-filter should be used in-line with the final virus filter under the same conditions. The filtration can be realized at any pressure up to 5 bar | 73 psi pressure. The filter qualification was performed at an operating pressure of 2 bar | 29 psi to meet customer standard requirements. The flow rate achieved through the Virosart® HC Minisart® during the filtration is a function of the purity of the feed stream, the target molecule size and concentration. The flow rate will decline with the progress of filtration and indicate the respective blockage of the filter.

After filtration the Virosart® HC Minisart® can be flushed with buffer to maximize protein transmission. Therefore please flush the Minisart® with >3 ml of buffer. Make sure that no air is entrapped within the Virosart® HC Minisart® at this step.

8. Post Use Integrity Test for Virosart® HC Minisart®

After filtration, an integrity test can be performed. Therefore, please submerge the Minisart® still connected to the filtration system, into a beaker of WFI or WFI equivalent. Now please pressurize the Minisart® carefully up to 4.5 bar | 65 psi.

The specification for a passed integrity test is that no air bubbles may be detected at the Minisart® outlet for 30 seconds.

9. Changing out the Virosart® HC Minisart®

The Virosart® HC Minisart® must be changed, at the latest, when the maximum allowable differential pressure is reached.

10. Applications Support by Sartorius Stedim Biotech

If required, please contact your local Sartorius Stedim Biotech representative to obtain further information on technical data, product integrity testing or general information concerning specific applications. Our application specialists from Purification Technologies will gladly assist in your filtration trial work if needed. Moreover our validation service and support team are always available to assist with questions concerning product and process validation, as well as training of staff and optimization of filtration processes.

11. Return of used filter elements

If you wish to return used filter elements to Sartorius Stedim Biotech for inspection, it is necessary to decontaminate and sterilize them prior to shipment. This procedure has to be documented in the return shipment form, available from your local Sartorius Stedim Biotech representative. This documentation is required by country law. Returned filter elements cannot be handled unless the appropriate return shipment form is completed.

12. Liability

Sartorius Stedim Biotech cannot assume liability if Virosart® HC Minisart® are subject to improper use. In the interest of product development we reserve the right to make changes.

Thank you for working with Sartorius Stedim Biotech, we appreciate your business.

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Sartorius Stedim Biotech GmbH