

High Recovery of Cationised Protein in Centrisart® I



Application
Note

#07

#08

#09

#10

#11

Abstract | Brief

Protein cationisation is a method used in biochemical research whereby a protein of interest (POI) is chemically modified to produce a protein with an excess of overall positive charge compared to its native form; sometime referred to as a 'super-charged' protein. Protein cationisation can be used for a variety of purposes; a well-documented use is of cationised BSA due to its greatly enhanced immunogenicity compared to non cationised BSA.

Here, a 38 kDa protein was isolated and purified from *E. coli* cell lysates but at a low concentration of 0.1 mg/mL. Freeze drying (lyophilisation) was used initially due to its high yield results to increase concentration, but as a dialysis step is required upstream to remove unwanted salts, protein was lost due to some aggregation. Lyophilisation is also a time consuming process involving multiple steps.

Use of a standard PES ultrafiltration device although a much faster process proved insufficient for the level of protein recovery required. The Centrisart[®] I provided the ideal solution to this issue with its highly non-specific binding CTA membrane and novel self-cleaning method of concentration, giving both a fast process time compared to lyophilisation methods and a high recovery compared to alternative centrifugal ultrafiltration devices.

Equipment

- Centrisart[®] I, 10K CTA (Sartorius, 13239-E)
- Allegra X12-R with SX4750 swing out rotor (Beckman Coulter)
- Cary-60 (Agilent)
- BioDesign dialysis tubing (Fisher Scientific)

Method

1. The isolated protein was cationised with N,N'-Dimethyl-1,3-propanediamine (DMPA) using 1-ethyl-3-(3-dimethylaminopropyl) carbodiimide (EDC) activation.
2. The cationised protein was dialysed into a 10 mM phosphate buffer (pH7).
3. 2 mL of dialysed cationised protein was loaded into each Centrisart 1[®] 10K CTA device.
4. Centrisart devices were centrifuged at 2,500 xg for 30 min.
5. Protein recovery was quantified by UV-vis.

Results and discussion

Protein recovery of 90 % was obtained using the Centrisart[®] I for the protein concentration process. At higher sample concentrations, additional downstream analysis was possible, such as MALDI mass spectrometry for confirming protein presence, dynamic light scattering and zeta potential for determining the particle size and charge, respectively.

Summary

The Centrisart[®] I device demonstrated its use as a low process time, high yield concentrator for a purified cationised protein. By using the Centrisart[®] I over traditional lyophilisation methods process time was compressed and by using CTA over PES membrane increase protein yield. The Centrisart[®] I device demonstrates that choice of the correct ultrafiltration product for a given sample is critical and if correct can dramatically improve the process workflow for a given sample in laboratory research.

Notes


Details of this protein and the exact method are proprietary information. We respect the individual's right to reserve certain information. However, additional information on cationisation and similar methods can be found in the references section.

Testimonial

"Very simple to use."

"Impressed with ease of recovery: pipettes fit into both main compartment and floater, no narrow sections to consume protein."

"Excellent sample recovery compared to other concentrators I've used."



Above: Schematic use of ultrafiltration in standard laboratory sample preparation compared to traditional method.

1. Purified sample is obtained.
2. Sample is loaded into lyophilisation equipment or Centriscart® I for centrifugation.
3. Sample is concentrated to aid downstream analytics.

References

<http://www.piercenet.com/method/carbodiimide-crosslinkerchemistry>

Sharma, K.P., Collins, A.M., Perriman, A.W., and Mann, S. (2013) Enzymatically active self-standing protein-polymer surfactant films prepared by hierarchical self-assembly. *Adv Mater.* 25, 2005-2010

Futami, J., Kitazoe, M., Maeda, T., Nukui, E., Sakaguchi, M., Kosaka, J., Miyazaki, M., Kosaka, M., Tada, H., Seno, M., Sasaki, J., Huh, N.H., Namba, M., and Yamada, H. (2005) Intracellular delivery of proteins into mammalian living cells by polyethylenimine-cationization. *Journal of bioscience and bioengineering.* 99, 95-103

Lawrence, M.S., Phillips, K.J., and Liu, D.R. (2007) Supercharging proteins can impart unusual resilience. *Journal of the American Chemical Society.* 129, 10110-10112

Apple, R J., Domen P L., Muckerheide A and Michael J G. (1988) Cationization of protein antigens. IV. Increased antigen uptake by antigen-presenting cells *The Journal of Immunology.* 140 10 3290-3295

J.Gabriel Michael. (1991)

Cationization of protein antigens: VI. Effects of cationization on the immunoregulatory properties of a bovine serum albumin peptide, a.a. 506-589

Cellular immunology. 138, 1 121-129

Sales and Service Contacts

For further contacts, visit www.sartorius.com

Europe

Germany

Sartorius Lab Instruments
GmbH & Co. KG
Weender Landstrasse 94-108
37075 Goettingen

Phone +49.551.308.0
Fax +49.551.308.3289

France & Suisse Romande

Sartorius France
2, rue Antoine Laurent de Lavoisier
ZA de la Gaudrée
91410 Dourdan

Phone +33.1.70.62.50.00
Fax +33.1.64.59.76.39

Austria

Sartorius Austria GmbH
Franzosengraben 12
1030 Vienna

Phone +43.1.7965760.0
Fax +43.1.7965760.24

Belgium

Sartorius Belgium N.V.
Leuvensesteenweg, 248/B
1800 Vilvoorde

Phone +32.2.756.06.71
Fax +32.2.253.45.95

Finland & Baltics

Sartorius Biohit Liquid Handling Oy
Laippatie 1
00880 Helsinki

Phone +358.9.755.951
Fax +358.9.755.95.200

Hungary

Sartorius Hungária Kft.
Kagyló u. 5.
2092 Budakeszi

Phone +3623.457.227
Fax +3623.457.147

Ireland

Sartorius Ireland Ltd.
Unit 41, The Business Centre
Stadium Business Park
Ballycoolin Road
Dublin 11

Phone +353.1.8089050
Fax +353.1.8089388

Italy

Sartorius Italy S.r.l.
Viale A. Casati, 4
20835 Muggiò (MB)

Phone +39.039.4659.1
Fax +39.039.4659.88

Netherlands

Sartorius Netherlands B.V.

Phone +31.30.60.53.001
Fax +31.30.60.52.917

info.netherlands@sartorius.com

Poland

Sartorius Poland sp.z o.o.
ul. Wrzesinska 70
62-025 Kostrzyn

Phone +48.61.6473830
Fax +48.61.6473839

Russian Federation

LLC "Sartorius ICR" and LLC "Biohit"
Uralskaya str. 4, Lit. B
199155, Saint-Petersburg

Phone +7.812.327.5.327
Fax +7.812.327.5.323

Spain & Portugal

Sartorius Spain, S.A.
Avda. de la Industria, 32
Edificio PAYMA
28108 Alcobendas (Madrid)

Phone Spain +34.902.123.367
Phone Portugal +351.800.855.800
Fax Spain +34.91.358.96.23
Fax Portugal +351.800.855.799

Switzerland

Sartorius Mechatronics Switzerland AG
Ringstrasse 24a
8317 Tagelswangen (ZH)

Phone +41.44.746.50.00
Fax +41.44.746.50.50

U.K.

Sartorius UK Ltd.
Longmead Business Centre
Blenheim Road, Epsom
Surrey KT19 9QQ

Phone +44.1372.737159
Fax +44.1372.726171

Ukraine

LLC "Biohit"
Post Box 440 "B"
01001 Kiev, Ukraine

Phone +380.44.411.4918
Fax +380.50.623.3162

America

USA

Sartorius Corporation
5 Orville Drive, Suite 200
Bohemia, NY 11716

Phone +1.631.254.4249
Toll-free +1.800.635.2906
Fax +1.631.254.4253

Argentina

Sartorius Argentina S.A.
Int. A. Avalos 4251
B1605ECS Munro
Buenos Aires

Phone +54.11.4721.0505
Fax +54.11.4762.2333

Brazil

Sartorius do Brasil Ltda
Avenida Senador Vergueiro 2962
São Bernardo do Campo
CEP 09600-000 - SP- Brasil

Phone +55.11.4362.8900
Fax +55.11.4362.8901

Canada

Sartorius Canada Inc.
2179 Dunwin Drive #4
Mississauga, ON L5L 1X2

Phone +1.905.569.7977
Toll-Free +1.800.668.4234
Fax +1.905.569.7021

Mexico

Sartorius de México S.A. de C.V.
Circuito Circunvalación Poniente
No. 149
Ciudad Satélite
53100, Estado de México
México

Phone +52.5555.62.1102
Fax +52.5555.62.2942

Asia | Pacific

Australia

Sartorius Australia Pty. Ltd.
Unit 5, 7-11 Rodeo Drive
Dandenong South Vic 3175

Phone +61.3.8762.1800
Fax +61.3.8762.1828

China

Sartorius Scientific
Instruments (Beijing) Co., Ltd.
33 Yu An Road, Airport Industrial Park
Zone B,
Shunyi District, Beijing 101300,
P.R.China

Phone +86.10.8042.6300
Fax +86.10.8042.6486

Hong Kong

Sartorius Hong Kong Ltd.
Unit 1012, Lu Plaza
2 Wing Yip Street
Kwun Tong
Kowloon, Hong Kong

Phone +852.2774.2678
Fax +852.2766.3526

India

Sartorius Weighing India Pvt. Ltd.
#69/2-69/3, NH 48, Jakkasandra,
Nelamangala Tq
562 123 Bangalore, India

Phone +91.80.4350.5250
Fax +91.80.4350.5253

Japan

Sartorius Japan K.K.
4th Fl., Daiwa Shinagawa North Bldg.
8-11, Kita-Shinagawa 1-chome
Shinagawa-ku, Tokyo, 140-0001 Japan

Phone +81.3.3740.5408
Fax +81.3.3740.5406

Malaysia

Sartorius Malaysia Sdn. Bhd
Lot L3-E-3B, Enterprise 4
Technology Park Malaysia
Bukit Jalil
57000 Kuala Lumpur, Malaysia

Phone +60.3.8996.0622
Fax +60.3.8996.0755

Singapore

Sartorius Singapore Pte. Ltd
1 Science Park Road,
The Capricorn, #05-08A,
Singapore Science Park II
Singapore 117528

Phone +65.6872.3966
Fax +65.6778.2494

South Korea

Sartorius Korea Ltd.
8th Floor, Solid Space B/D,
PanGyoYeok-Ro 220, BunDang-Gu
SeongNam-Si, GyeongGi-Do, 463-400

Phone +82.31.622.5700
Fax +82.31.622.5799

Thailand

Sartorius (Thailand) Co. Ltd.
129 Rama 9 Road,
Huaykwang
Bangkok 10310

Phone +66.2643.8361-6
Fax +66.2643.8367



◀ www.sartorius.com