
sartorius stedim
biotech

Operating Manual | Bedienungsanleitung | Mode d'emploi | Manuale d'uso | Manual de funcionamiento

BioWelder[®] Total Containment

From serial number BWTC#2000

1000046827

Ver2019-03

Contents

English	Page 3
Deutsch	page 47
Français	Page 93
Italiano.....	página 139
Español.....	pagina 185

Table of Contents

1 Hazard Information	5	6.3 Web browser	26
2 Conventions used in this document	7	6.3.1 Establish connection	27
3 General	7	6.3.2 Menu tree	27
3.1 Warranty	7	6.3.3 Operator controls	28
3.2 Liability disclaimer	8	6.3.4 Browser representation	29
3.3 Patent	8	6.3.5 General settings	29
4 Introduction	9	6.3.6 User login/logout	30
4.1 Features	9	6.3.7 Date and time	31
4.2 Safety	9	6.3.8 User management	31
4.3 Limitations	9	6.3.9 Log	32
4.4 Function	10	6.3.10 Info	32
4.5 Blocked device	10	6.4 Automatic Logout	32
4.6 Cleaning instruction	11	7 Sequences	33
4.7 General Cleaning	11	7.1 Initialise axes	33
4.8 Periodic Cleaning	11	7.2 Weld	33
4.9 Cleaning the Thermopile Sensor	11	7.2.1 Critical section	33
4.10 Decontamination	12	7.2.2 Abort in the critical section	34
5 Start-up	13	7.2.3 Reweld	35
5.1 Delivery kit	13	7.3 Reverse	35
5.2 Familiarizing yourself with the device	13	7.4 Temperature Validation	36
5.3 Installation	14	8 Operational problems	37
5.4 Login	14	8.1 Switching off device	37
5.5 Maintenance	14	8.2 Foreign objects in the device	37
5.6 Entering settings	14	8.3 Wrong type of tube inserted	37
5.7 Initializing axes	14	8.4 Wrong tube dimension in the tube holder	37
5.8 Blade handling	15	8.5 Memory card	37
5.9 Welding	16	8.6 Device reverted	38
6 Operation	18	8.7 Network	38
6.1 Mutual exclusion	18	8.8 Web browser settings	38
6.2 Touch screen	19	8.9 Password forgotten	38
6.2.1 Menu tree	19	9 Appendix	39
6.2.2 Operator controls	20	9.1 Copyright	39
6.2.3 Summary of symbols	21	9.2 Disposal	39
6.2.4 User login/logout	22	9.3 Technical data	39
6.2.5 Date and time	22	9.4 Technical contact	40
6.2.6 Language	23	9.5 External size of device	40
6.2.7 Type of tube	23	9.6 User groups and user interfaces	41
6.2.8 Network parameters for reset	23	9.7 Log	42
6.2.9 User management	24	9.8 Tube holder identification	43
6.2.10 Log	25	9.9 Declaration of conformity	44
6.2.11 Info	26	9.10 Error reports and Warning messages	45

Technical Safety Information for the User

This description contains the necessary information for the designated use of the product described herein. It is intended for technically qualified personnel.

Qualified personnel are persons who, because of their education, experience and training as well as their knowledge of relevant standards, specifications, accident-prevention regulations and operating conditions, have been authorized by persons responsible for system safety to execute the necessary activities, and can recognize and avoid possible dangers in such activities.

1 Hazard Information

The following information is for the personal safety of operating personnel as well as the safety of the described product.

Hand injury warning
Never grab hold of anything inside openings during operation.
Noncompliance can lead to serious injury.

In case of emergency, press emergency off button!
This will shut down the machine immediately and will result in a power failure sequence.

To disconnect the BWTC from the mains supply, remove the power cord from the machine.

Remove power cord before cleaning.
Never open up the device!

Electrostatic sensitive devices (ESD)
Ground straps, table mats and tools with conducting grips are used to ground people, to achieve equipotential welding and to prevent electrostatic discharge. This is always a requisite measure if electronics or electronic components are handled or mounted.

Make sure that rotating tubes do not pose a safety risk during the welding process.

Follow the accident-prevention and safety regulations relating to the specific application.

Ground conductor connections must function perfectly.

The device is not intended for use in/with explosive atmospheres.

A discontinuity during the critical section of the welding process can result in damage and is to be absolutely avoided. The device will attempt to limit the damage. It is imperative that you follow the operating instructions. Read the corresponding Chapter 7.2.2 Abort in the critical section

In case of damage, have the device cleaned and examined by an authorized specialist.

The device must be decontaminated before being sent to suppliers or manufacturers.

The surface in this area can have very high temperatures.
Skin contact can result in burns | injury.
Avoid direct contact with skin in this area and wear corresponding protective clothing | equipment.
Always use the blade remover tool (accessory) to remove the blade.

Before start-up check whether the rated voltage of the device coincides with the local line voltage.

The device was developed for laboratory use and must be protected from every type of fluid or vapour.

The device can be operated between +5°C and +40°C.
The welding parameters are qualified for room temperature (20°C–25°C) operation only.

Optimal relative humidity is 60%.

The device is programmed with standard parameter sets for welding C-Flex[®] 374, AdvantaFlex[®], PharMed[®] BPT and SaniPure[™] BDF[™].
These parameter sets have been validated at ambient temperature.

The degree of protection is IP20.

Make sure that the device is not switched off during the welding process.

Ensure that there is an uninterrupted power supply (UPS).

Use only the blades provided by the manufacturer and use it only once.

Use only tested control cables and accessories.

Isolate the supply voltage before replacing a fuse.

The device has an Ethernet connection. Only skilled personnel may connect the device to a network.

Only shielded Ethernet cables may be used.

Only skilled personnel may configure the device for connection to a network.

It is recommended to operate the device via the touch screen during sequences.

The device contains electric drives which can be blocked | destroyed by foreign substances. Make sure that no foreign substances reach the interior of the device.

After using the BioWelder[®] TC, you should always log out. Otherwise another user is able to manipulate the device with your account.

The auto login is used to work more efficiently with the BioWelder[®] TC.
The auto login can only be activated for one user. This user will be logged in AUTOMATICALLY without entering any username and password after each restart or after changing the user interface.

A user of the Supervisor group should NOT use the auto login because of safety reasons, otherwise any user is able to take essential changes on the BioWelder[®] TC under the Supervisor account.

Not following the manual instructions of the manufacturer may lead in to unreliability, hazardous situations and disturbances of safety.

2 Conventions used in this document

This document is intended for users of the BioWelder[®] TC from the Supervisor and Operator user group.

Screenshots from a Web browser were created using Internet Explorer 8.

Warning: Situations in which you or others could be injured, as well as, information about process safety

Information: Tip or additional information

▷ This symbol indicates other noteworthy aspects or information.

1., 2., ... Numbered paragraphs indicate operating steps that should be performed in the order given.

3 General

Please read this operating manual carefully. It allows you to start up and operate the device safely and successfully.

3.1 Warranty

We shall, in principle, be liable for defects for a period of one year after delivery. This warranty excludes the spring contact and cap, the crosshead tube holder disengaging assembly and the shaft blade fixing.

In the event of a defect we shall, at our own discretion, be entitled by way of subsequent performance to effect a replacement delivery or to rectify the defect. In the event of a defect in title, we shall be entitled to remedy such defect by modifying the delivered item in a manner reasonably acceptable for the customer. The customer shall only be entitled to claim damages in lieu of performance in the event of intentional or grossly negligent breach of duty on the part of Sartorius Stedim Biotech. If in the event of a defect in title modification is not feasible on economically reasonable terms or within a reasonable period, both the customer and we shall be entitled to withdraw from the contract. Where Sartorius Stedim Biotech has provided a warranty as to durability, the customer shall be required to furnish evidence that the material defect existed at the time of transfer of risk, unless the customer cannot reasonably be expected to do so in individual cases. In all other respects, the provisions contained below under "Liability" shall also apply to the warranty.

3.2 Liability disclaimer

Proper use

Sartorius Stedim Biotech will not accept any liability for faults and damage that occur through improper handling and because of inappropriate use of the device. Improper handling takes place in the case of noncompliance with the operating manual, especially instructions for start-up. Return the device only in the original packaging or in appropriately secure packaging. We accept no liability for transit damage. The BioWelder® TC was designed for industrial use and may not be used for human medical applications. We accept no liability for errors made when connecting the device to a network. The system is intended only for the purpose mentioned in the manual. Another use above and beyond this or retrofitting of the device without a written agreement with the Sartorius Stedim Biotech is considered improper use. Sartorius Stedim Biotech is not liable for any damage resulting from this. The operator alone bears the risk.

Liability

In the event of intent or gross negligence on the part of legal representatives or employees at executive level, and in the event of culpable breach of material contractual duties by the same or other employees, Sartorius Stedim Biotech shall be liable in accordance with the legal provisions. With regard to any breach of minor contractual duties by other employees, Sartorius Stedim Biotech shall not be liable, nor in the event of slight negligence of legal representatives or of employees at executive level. Sartorius Stedim Biotech shall not be liable for losses that are not typically associated with this type of contract or that are scarcely foreseeable. This shall apply, in particular, without being limited to, claims due to missed profit and consequential damages. With regard to claims under Product Liability Law, such Law shall apply without restriction. Claims for damages in tort shall only apply if they are caused intentionally or by gross negligence. This shall also apply with regard to actions of vicarious agents of Sartorius Stedim Biotech. Sartorius Stedim Biotech shall be liable in accordance with the statutory provisions in the event of culpable physical injury and hazards to life and health.

If determination of fault is required for any claim of the customer, the burden of evidence shall rest with the customer unless otherwise prescribed by law and the customer cannot reasonably be expected to bear the burden of evidence.

3.3 Patent

BioWelder® Total Containment is a registered trademark of Sartorius Stedim Biotech GmbH.

The manufacturer reserves the right to make any changes and improvements in the product and process described in this manual at any time without notice.

All registered trademarks and patents are the property of their respective owners.

4 Introduction

The BioWelder[®] TC works by means of a thermal welding technique that uses no laminar flow. In this welding technique the blade is kept sterile and pyrogen-free during the welding process. This occurs according to the reaction equation or Arrhenius's law by means of a temperature and time constant (Wallhäuser: The D-value resistance indicator for bacillus subtilis is less than 1 sec at 200°C of dry heat). The inserted pre-sterilized tubing (gamma irradiated and/or autoclaved) tubes are separated by the heated blade and thermally welded to each other by the device. The welding is completely automated and is released only after a finished cycle. BioWelder[®] TC was developed and programmed to weld C-Flex[®] 374, Advantaflex, Pharmed BPT and SaniPure[™] BDF[™]. These parameters are validated at room temperature.

4.1 Features

- ▷ Fully automatic mode
- ▷ Welds thermoplastic tubes from different manufacturers up to 1 inch
- ▷ Liquid filled tubes
- ▷ Welds waste tubes without spilling fluid
- ▷ Suitable for the taking of samples
- ▷ Pre-installed standard settings for the most common tubes
- ▷ Reliable thermal welding
- ▷ Aseptic connections without laminar-flow bank
- ▷ Data logging and output to memory card or via Web browser to the computer
- ▷ Easy to validate with temperature validation set (option)
- ▷ Compact and portable design
- ▷ Easily expandable for new tube materials and tube geometries
- ▷ Simple operation via a touch screen
- ▷ Simple operation via a Web browser without additional installation of software
- ▷ Pre-installed languages, additional languages upon request

4.2 Safety

As the blade can reach temperatures of over 400°C during the welding process and drives move with high forces, the cover must not be opened during the process to avoid injury. The device releases the cover automatically as soon as the blade has cooled off adequately. Use only tubes that are compatible with the device. Use only BioWelder[®] TC single use blades from Sartorius. Use the blades only once to obtain optimum results. Use the blade remover tool (accessory) all the time to remove the blade from the device.

4.3 Limitations

The welds made with the BioWelder[®] TC are guaranteed to be of sufficient quality and strength for sterile fluid transfer. However, it is ultimately the responsibility of the end consumer to check and validate the device for suitability. Under no circumstances can data from this operating manual be presented to an official body as reference to manufacture a clinical product without the written consent of the Sartorius Stedim Biotech. The BioWelder[®] TC was designed for industrial use and may not be used for medical applications. Please refer to Chapter 1 Hazard Information for additional limitations. The connection created by the BioWelder[®] TC cannot meet the requirements for medical applications. It is the responsibility of the user to decide whether the device is suitable. Examine welds thoroughly before using them.

4.4 Function

The BioWelder[®] TC is used for the sterile welding of two tube pieces. Every piece of tube is connected on one end to a container, bag or process unit. The other end is sealed with a clamp, plug or another seal. To create a sterile connection, first place the two tubes in the tube holder of the BioWelder[®] TC. Then a new blade is inserted into the blade holders of the device. Closing the cover and pressing the start button activates the steps in the welding programme. In the first step the blade is heated to over 400°C and held for 5 seconds at this temperature. This step guarantees that the blade is sterile and pyrogen-free. Then the blade is cooled to the required cutting and welding temperature which depends on the type of tube. Afterwards the heated blade simultaneously cuts through the two pieces of tube. The blade temperature is controlled during the cutting process. Once the tubes have been cut through, the device turns the tubes so that they line up. The blade is then pulled back, the tubes pressed together and the weld cooled off. The release of the cover signals the end of the welding and the device can be opened. The welded tube is taken out and is now ready for use. The used blade and the waste tube are removed. The BioWelder[®] TC is a fully automatic unit and displays the successful conclusion of the welding if all steps were followed properly. The process is completely reproducible.

The device is programmed with standard parameter sets for welding C-Flex[®] 374, AdvantaFlex[®], PharMed[®] BPT and SaniPure[™] BDF[™]. These parameter sets have been validated at room temperature.

The use of a new blade after each welding prevents cross-contamination and to ensure proper welding quality. To simplify the logging of welding data for the user, the BioWelder[®] TC has a memory card to which data can be written as needed. In addition, data can be loaded to a computer through a Web browser. To simplify the validation or to facilitate the validation process for the customer, a validation set matching the device can be acquired for the BioWelder[®] TC as an option.

4.5 Blocked device

Check all equipment inside the BioWelder[®] TC device continuously for operational reliability. In addition, monitor every step during a sequence. If equipment such as a drive or the heating gives out, or the execution of a step in a sequence fails, then an error is present. Such a fault leads to a blocked device. When the device is blocked, all on-going processes are stopped and the operator control displays this message: "All runs blocked (X)"

A blocked device is present if the operator control displays the following text: "All runs blocked (X)"

If a sequence of a process cannot be finished for some reason, this does not necessarily indicate a defective device. It could be that the user has chosen a wrong type of tube, for example.

A blocked device does NOT necessarily contain a defective component. It is possible that a sequence could not be completed successfully.

If there is a blocked device because of a failed sequence, in most cases the device can be reversed. In such a case the blockage can usually be resolved by switching the device on and off.

If there is no equipment defect, in most cases a blockage in the BioWelder[®] TC can be remedied by switching it off and then switching it on again.

If there is an equipment defect, the device must be repaired. In this case, contact your Sartorius Stedim Biotech Service Representative.

Before you contact your Sartorius Stedim Biotech Service Representative, describe the fault in detail and save the "All" log to the SD card as described in Chapter 6.2.10 Log. Send these files together with the fault description to your Sartorius Stedim Biotech Service Representative.

4.6 Cleaning instruction

Before cleaning the BioWelder[®] TC the cleaning cover for the temperature sensor must be installed.

Before cleaning the cleaning agent must be tested on compatibility with the respective materials.

The BioWelder[®] TC may not be cleaned with a spray- or steam procedure. Spray- or steam cleaning can cause a reduced life time of single components or lead to damage of single components.

The manufacturer disclaims any liability if components are damaged by improper cleaning or in the case of non-observance of the cleaning instructions or when persons are being injured by not observing the proper cleaning instructions.

4.7 General Cleaning

Immerse a cleaning tissue with a standard disinfectant, e. g. microcide, and clean the contaminated parts.

4.8 Periodic Cleaning

To increase process security it is recommended to periodically clean the contact springs of the Guide Block - Turning squeeze axis, the contact points for the tubing holders and the tubing holders in the pinching area. For cleaning the contact springs and the contact points, it is recommended to use a cotton bud which is immersed with a standard disinfectant, e. g. microcide.

4.9 Cleaning the Thermopile Sensor

The cleaning of the temperature sensor should only be performed in case the temperatures are outside of the range of tolerance.

After cleaning, it must be ensured that the visibility field of the temperature sensor does not bear any residues.

For cleaning the temperature sensor will be removed from the BioWelder[®] TC. Use an appropriate cotton bud with a max. diameter of 4 mm and immerse it with a standard disinfectant, e. g. microcide, and clean the visibility field of the temperature sensor.

4.10 Decontamination

The machine must be decontaminated with either:

- Isopropanol 70% (e.g. Klericide 70/30)
- Sodium hypochlorite 0.75% (e.g. Klericide-CR)
- Propanol 45%, Isopropanol 25%, Ethanol 4.7% (e.g. Bacillol AF)

The customer is responsible for correct decontamination!

After correct decontamination, the customer must fill out the "Decontamination form" that the machine is free from biological, chemical or radioactive contaminants. This form can be get from the local Sartorius representative.

5 Start-up


5.1 Delivery kit

This Chapter is a step-by-step guide for the start-up of the BioWelder® TC.

- ▷ BioWelder® TC device
- ▷ Blade remover tool
- ▷ Tube holder set
- ▷ Validation set, optional
- ▷ Power cord
- ▷ QC samples made with your unit
- ▷ SD memory card inserted into slot, optional
- ▷ Quick start guide
- ▷ Operating manual
- ▷ BioWelder® TC certificate
- ▷ Acceptance report
- ▷ EC-declaration of conformity

5.2 Familiarizing yourself with the device

Item	Description
1	Blade remover tool
2	Tube holder set (by pair)
3	Validation sensor with holder
4	SD memory card slot
5	Ethernet cable terminal
6	Device socket with main switch and fuse carrier
7	Touch screen
8	Closing handle for cover
9	Carrying handle, front
10	Carrying handle, rear
11	Nameplate with serial number
12	Emergency off button

5.3 Installation

The BioWelder® TC is portable and can be used in any location. Use a secure and horizontal base for the installation. Make sure to place the machine in a way that the access to the on/off switch and the emergency off button is well accessible! Connect the power cord with safety earthing to the device socket and switch it on. The power supply must be according to the technical data. After a brief start-up section the device is ready to operate.

The raiser is responsible for the safety of the system and its integration.

5.4 Login

The device is set upon delivery so that a user with the name "Supervisor" from the "Supervisor" user group or the name "Example" from the "Operator" user group can login without a password. In this case, the instructions needed to execute initialisation are displayed on the touch screen.

An automatic login can be activated. See Chapter 6.2.9 User management

5.5 Maintenance

There is no spare part list for BioWelder® TC available for the customer. Sartorius recommends to have a yearly intervention and maintenance executed by your local Sartorius Service representative that includes temperature sensor validation. Please contact your local Service representative from Sartorius to establish a maintenance contract.

This information has to be confirmed with the OK button on the screen. Contact your Sartorius Stedim Biotech Service Representative, if maintenance is needed.

5.6 Entering settings

Different settings such as the following can be entered from a Supervisor of the "Supervisor" user group via the menu: language, date, time and type of tube.

Read the corresponding Chapter 6.2 Touch screen

5.7 Initializing axes

After the device starts up, an initialisation must be executed. Follow the instructions on the touch screen.

**The tube holder and the blade are removed from the device.
The user then closes the cover.**

5.8 Blade handling

To insert a blade into the device, proceed as follows:

Insert the blade by handling it with fingers avoiding touching the centre.

– Ensure the blade is completely pushed to the end of the blade holder.

Pay attention! The blade must be pushed completely to the end to avoid risk of burning and damaging the device.

Follow the instruction displayed on the screen for the steps until the weld is completed and the tubing removed.

1. Take the blade remover tool in your hand.

2. Push the blade remover tool over the blade.

3. Pull it back together with blade.

5.9 Welding

The device is programmed with standard parameter sets for welding C-Flex[®] 374, AdvantaFlex[®], PharMed[®] BPT and SaniPure[™] BDF[™].

These parameter sets have been validated at room temperature.

To execute welding, follow the instructions on the touch screen. If the device has not been previously initialised, no welding can be executed at this point.

Proceed as follows:

Tube holder and tubes are available, and the right type of tube is inserted.
See Chapter 6.2.7 Type of tube

1. Insert blade (see Chapter 5.8 Blade handling for details)

Make sure that the blade is correctly inserted up to the stop in the holder provided for it. See Chapter 5.8 Blade handling

2. Insert tube holder set

The tube holders are inserted so that the clamp bars are facing each other.

3. Insert the tubes into the tube holder and close cover.

Clean the exterior of the tubes using a suitable disinfectant (see Chapter 4.6 Cleaning instruction).
Insert the tubes carefully on the holders up to the stop.

Make sure that the diameters of the tube holder and of the tubes match.
Use a minimum of 450 mm tubing



For reasons of safety, seal off the tubes on all sides with a clamp, plug or another seal before every welding. Attach the clamps from a distance between 2 and 10 cm close to the device but outside the closed cover.


Make sure that the tubes on the left side of the device can rotate freely. Otherwise, this can lead to a discontinuation of welding.

No tension or stress shall be applied on tubing during the entire welding process.

4. Start welding sequence

Check on the touch screen whether the type of tube entered matches with the type of tube inserted. Check inside and outside diameter as well. If all specifications match, press the start button to start the welding sequence.

The welding is finished when the cover opens. Afterwards proceed as follows:


1. Remove blade with the blade remover tool (ref to Chapter 5.8 Blade handling)

The blade cooling remains switched on until the blade is removed by the user or the after-cooling time of two minutes have elapsed. Each blade may only be used once.

2. Remove tubes from the holder

If the welding seam sticks together, pinch it gently to open up.

3. Check quality of welding

For an optimal connection, let the tube cool for about one minute before checking the welding and before the thermal weld is mechanically loaded.

4. Dispose of waste tube appropriately

The welding is now finished and a new welding sequence can be started.

6 Operation

The BioWelder® TC comes with two user interfaces which cover some of the common functionalities. The device can either be operated via the touch screen integrated into the device or via Ethernet through a Web browser.

Wherever possible the user should operate the device using the local touch screen. Switching between touch screen and Web browser is not recommended.

The following chapter describes the differences between the two user interfaces.

6.1 Mutual exclusion

The following rules are implemented so that it becomes clear which user interface has control of the device:

You can only work with one user interface at a time. The other user interface remains inactive. When the device is started, neither user interface is active or inactive. Web browser and touch screen become active by pressing an operational control, i.e. by pressing any button.

Remote control
Touch to take
control!

Display in an inactive touch screen user interface:

Display in an inactive Web browser user interface:

Rules for activating/switching over to a user interface:


The touch screen becomes immediately active if it is touched. The Web browser only becomes active if no user has logged in on the device or, after start-up, if an operational control in the Web browser is activated. If you wish to activate the Web browser user interface, you must log out of the touch screen. What they have in common is that switching over from one user interface to the other logs out the current user. If a user has activated automatic login, this user is logged in, otherwise the login dialogue appears.

6.2 Touch screen

The touch screen operates through a touch-sensitive screen. The operator controls displayed can be pressed directly, in the same way as buttons. If the touch screen will be operated with more than 1 finger it can lead to unexpected results.




6.2.1 Menu tree

The following diagram shows the menu tree of the touch screen for the Supervisor user group.

6.2.2 Operator controls

The touch screen has only a few different standardized operator controls which are also common in other user interfaces.

Operator control	Example	Meaning
Button		Start: Start process
Button		Stop: Discontinue process
Button		Buttons have a symbol but they can also be labelled with a text. See Chapter 6.2.3 Summary of symbols
Spinner		A spinner offers a selection which you can scroll through using the arrow button.
Spinner with input		User name input is possible with this spinner which allows for text input. To do so, touch the text of the selection and switch to the input window.
Text box		Changing the characters in a text box is treated in the same way as the spinner with input. A text box has a frame. A label indicates the meaning; in this case, password. Touching the area inside the frame opens the input window.
Input window		Alphanumeric inputs without special characters but with upper and lower case are possible. Switch from number input to upper and lower-case input using the "1>A>a" button. The keyboard changes accordingly. Pressing Del deletes the whole input. ← deletes the last number.
Menu		The user can navigate to different functions from the Menu. Navigate to the next menu using the arrow button. An Operator has only one menu page available. Exit the menu through the Cancel button.

6.2.3 Summary of symbols

The following tables show the meaning of the symbols found on the touch screen.

General symbols

Symbol	Meaning
	Cancel: Exit the current screen. Changes are discarded.
	OK: Changes are accepted and saved by the user.
	Meaning depends on the context. In a selection box the next object is selected, in a menu the next screen of the same level is selected.

Menu window 1/2

	User login, refer to Chapter 6.2.4 User login/logout
	User logout, refer to Chapter 6.2.4 User login/logout
	General information, see Chapter 6.2.11 Info
	Save log data to memory card, refer to Chapter 6.2.10 Log
	User management, refer to Chapter 6.2.9 User management
	Selection of tube type, refer to Chapter 6.2.7 Type of tube

Menu window 2/2

	Validation, see Chapter 7.4 Temperature Validation
	Language of user interfaces, refer to Chapter 6.2.6 Language
	Set date, refer to Chapter 6.2.5 Date and time
	Set time, refer to Chapter 6.2.5 Date and time
	Network parameters for reset, see Chapter 6.2.8 Network parameters for reset

6.2.4 User login/logout

Login: If no user has the automatic login attribute, the user must log in with a user name and password

Screen

Meaning


Enter the user name and the corresponding password in the User Login dialogue. Confirm the entries by pressing the OK button.

If a user has the automatic login attribute, then this user is automatically logged in every time the device is started or after switching modes between the touch screen and Web browser.

Logout: After finishing work on the device, the user logs out to ensure that no other user can tamper with the device under a wrong user name.

Screen

Meaning

A user logs out via the User Logout menu item by pressing the OK button on the machine.

6.2.5 Date and time

The date and the time are used as a time stamp for logs. The setting on the touch screen differs from that on the Web browser.

Screen

Meaning

Date:
The day, month and year can be adjusted using the arrow button. Pressing the number once increases it by one count. Holding down arrow button increases the number by 5 counts until the button is released.

Time:

The hour, minute and second can be adjusted using the arrow button. Pressing the number once increases it by one count. Holding down arrow button increases the number by 5 counts until the button is released.

6.2.6 Language

Screen

Meaning

Language:
The language is changed via the arrow button and confirmed with OK. The language is updated in the operator control when the machine is switched on and switched off.
Available language: English, German, French.

6.2.7 Type of tube

Screen

Meaning

Type of tube:
The type of tube is changed via the arrow button and confirmed with OK.

6.2.8 Network parameters for reset

Screen

Meaning

Network parameters for reset:
Using the OK button the network parameters are set to the following values:
DHCP: off
Static IP address: 192.168.1.199
Static subnet mask: 255.255.255.0
Static default gateway address: 192.168.1.1
Static Server address: 192.168.1.1

In case of a network conflict, the network parameters can be reset to the factory settings through this menu.

MAC address will be reset to default value.

6.2.9 User management

User management contains all users. A user can edit their user data or data from the same or a deeper user group. The Operator can only change their own password.

User management is divided into three different screens on the touch screen. For the input, please read Chapter 6.2.2 Operator controls. You can find the different screens and their parameters in the following table.

Parameter

Meaning

Screen 1:
Screen change with the down arrow

All users are offered for selection which the logged in user can change. The user can select a user for deletion or add a new user to the spinner with Input. Users cannot delete themselves.

Every user must belong to a user group. The user group defines the rights. The current user cannot change his or her user group. Chapter 9.5 User groups and user interfaces displays all users and their rights.

A password must be assigned. The password must be entered once.

Screen 2:
Screen change with the down arrow

All users are offered for selection which the logged in user can delete. Users cannot delete themselves.

Screen 3:
Screen change with the down arrow

Automatic Login can be assigned exclusively to one single user. When switching on the device and when changing the user interface, the login dialogue in this case is skipped and the user with Automatic Login is logged in. A user loses Automatic Login if Automatic Login has been assigned to another user.

6.2.10 Log

For general information about logs, please read the Chapter 9.7 Log.
Note: The BioWelder® TC is not CFR21 part 11 compliant.

The log output via the touch screen is always saved to the SD card. The sole purpose of the SD card is to save logs.

The device does not support every type of memory card. Therefore, a compatible card is supplied and it is imperative that this card be used. Incompatible memory cards can cause the machine to behave erroneously.

Logs of several devices can be saved to one card. These logs are saved in different folders. The folder name contains the serial number of the device and is hence unique. The following table shows an example of how the content of a memory card can look.

Folder name	Log files in the folder
BWTC-101	Raw.log State.xml Weld.csv Weld.html
BWTC-156	Weld.csv Weld.html

As you can see the content of the folders in this example is different, the reason being that the user has not always made the same selection in the output. Details about the selection criteria are listed in the following table.

Selection criterion	Meaning
Weld Log	All weld logs from the selected day to the last welding done.
All	Weld logs and Raw logs from the selected day to the last entry. The current status

Selecting the date

By selecting a date, the user can specify how many logs will be downloaded to the SD card. The date entered will set the start of the period for downloading. The end date is always the current day. If entered 01/01/2001, the whole memory content will be downloaded. This may take some time!

The saving process is shown by a progress display and as a sequence can be stopped.

When outputting logs to the SD card, all data saved on the device are outputted. When outputting logs via Web browser, only a portion of the saved data is outputted. See Chapter 9.7 Log

6.2.11 Info

In this menu the current status is displayed according to following table:

Page 1

Item	Meaning
IP address	Currently used IP address
DHCP	Yes or no, depending on whether DHCP is used
Serial number	Serial number of the device

Page 2

Item	Meaning
Firmware	Firmware version
File system	File system version
Weld Cycles	Number of welding cycles which the device has executed so far
Tube Holder Cycles	9999 (fixed number) - (no tubeholder inserted)
Wear	Equipment abrasion indicator. If the wearing down of one or more parts has reached $\geq 100\%$, information will appear on the operator screen. Refer to Chapter 5.5 Maintenance.

6.3 Web browser

The Operator and Supervisor users can perform their tasks equally with the touch screen or the Web browser.

Advantages of working with the Web browser:

No additional software needs to be installed.

Several BioWelder[®] TC devices can be operated simultaneously over a Web browser.

Operation via mobile devices such as laptop, tablet PCs, etc. is possible.

All known browsers on different platforms are supported.

The prerequisites for working via Web browser are:

A computer with Web browser and a network cable are available.

The network configuration in Chapter 6.3.5 General settings was correctly done.

The following Web browsers are recommended:

Internet Explorer from version 8

Firefox from version 7

Google Chrome from version 16


Greyed out menu items cannot be chosen as the prerequisites have not been met.

6.3.1 Establish connection

The current IP address is requested from the touch screen in the Info menu and it must be entered in the URL box of the browser as shown in the diagram below. See also Chapter 6.2.11 Info

Example: The IP address of the device is 192.168.1.104. Input in the URL box: "192.168.1.104" followed by the Enter key. The homepage of the device then appears as follows:





Depending on the browser used, shortcuts or buttons can be set up to operate several BioWelder® TC devices efficiently from one browser.






If the device is initialised with the Web browser at the time of welding, the Initialise Axes menu item is greyed out. If not, the Weld and Validation menu items are greyed out.

6.3.2 Menu tree

Menu	Meaning
	User login or logout, see Chapter 6.3.6 User login/logout
	Initialise axes, see Chapter 7.1 Initialise axes
	Weld, see Chapter 7.2 Weld
	Validation, see Chapter 7.4 Temperature Validation
	General settings, see Chapter 6.3.5 General settings
	User management, see Chapter 6.3.8 User management
	Info, see Chapter 6.3.10 Info

6.3.3 Operator controls

Element	Example	Meaning
Button		OK: Accept changes
Button		Cancel: Reject changes
Button		Login: User login
Button		Logout: User logout
Button		Start: Start process
Button		Stop: Discontinue process
Spinner		A spinner offers a selection
Spinner with input		One more input is possible with this Spinner with Input. You select the text "Type..." and then enter the corresponding text via the keyboard.
Check box		A functionality can be switched on or off with the Check box.
Text box		The input in the text box must be meaningful. The input will be checked in terms of its value range.

6.3.4 Browser representation

The user interface consists of a heading with product name and logo. The selectable menu items are located on the left side of the browser window. Device status or information is displayed in the main window depending on the selected menu. The following screen section shows an example of an initialised machine. Axis Initialisation is greyed out and hence cannot be selected. The Login/Logout menu is selected. A user with the Supervisor user name from the Supervisor group is logged in.

The screenshot shows the BioWelder TC web interface. The top header includes the product name 'BioWelder TC' and the Sartorius Stedim logo. On the left side, there is a vertical menu with buttons for 'Login/-out', 'Initialize Axes', 'Weld', 'Temperature Validation', 'General Parameters', 'User Management', and 'Info'. The 'Login/-out' button is selected. The main content area is titled 'Login/-out' and displays the message 'User 'Supervis' of group 'Supervisor' is logged in'. Below this, it prompts the user to 'Type or choose user name and password'. A form contains a 'User name' dropdown menu with 'Supervis' selected, a 'Password' input field, and 'Login' and 'Logout' buttons.

6.3.5 General settings

The general settings are used to connect the device to a network. The language and the standard tube type can also be selected here. General settings can be found under Web Browser -> General Parameters:

The screenshot shows the BioWelder TC web interface with the 'General Parameters' menu selected. The left sidebar menu is the same as in the previous screenshot. The main content area is titled 'General Parameters' and contains a table of settings:

Parameter	Action	Value
Tube type	Choose tube type	C-Flex 374
Auto logout time	in [minutes]	11
Language	Choose language	English
DHCP	Click to enable	<input checked="" type="checkbox"/>
MAC address	Type MAC address	00:0E:3A:10:05:F5
Static IP address	Type address	192.168.1.199
Static subnet mask	Type subnet mask	255.255.255.0
Static default gateway address	Type address	192.168.1.1
Static server address	Type address	192.168.70.1
Weld Log Screen	Click to enable	<input type="checkbox"/>

At the bottom of the settings table are 'OK' and 'Cancel' buttons.

General parameters	Description
Tube type	Preset standard tube type
Auto logout time	After the specified time, the BWTC logs out automatically if no user activity happens
Language	Operator control language The language is updated in the operator control when the machine is switched on and switched off.
DHCP	If the device is connected to a DHCP server and the user would like to work with a dynamically assigned IP, the check box must be activated.
MAC address	This address only needs to be adjusted if two devices in the same network have the same MAC address. The MAC addresses are assigned by the manufacturer and every device has its own address.
Static IP address	If the DHCP is off, then the user with this address can log in on the machine via Web browser
Static subnet mask	Subnet mask of the local network
Static default gateway address	Address of a possible gateway
Static server address	Address of the DHCP server Is used if DHCP is activated.
Weld Log Screen	If check box is active, a weld log screen will be shown after each weld with according informations.

Only a specialist can modify the network settings. Improper use can cause the machine to stop responding via the Web browser. New network parameters are accepted when the machine is switched on and off.

6.3.6 User login/logout

Login:

If no user has the automatic login attribute, the user must log in with a user name and password. To do this, select the user name in the Login/Logout dialogue and enter the corresponding password. Confirm the entries by pressing the login button.

Logout:

A user logs out by pressing the Logout button on the machine.

Login/-out

User 'Supervis' of group 'Supervisor' is logged in

Type or choose user name and password

User name	Password		
Supervis ▼	<input type="text"/>	Login	Logout

If a user has the automatic login attribute, then this user is automatically logged in every time the device is started or after switching modes between touch screen and Web browser. Automatic logout is than deactivated.

6.3.7 Date and time

Synchronization at the touch of a button with the computer on which the browser is running is available under the Info menu:
Synch BioWelder® TC with this Host.

```
Date and Time
On BioWelder TC: Dienstag, 11. Oktober 2011 15:14:30
On this host:      Dienstag, 11. Oktober 2011 15:14:33
Synch BioWelder TC with this Host
```

Date and time on the computer must be correctly set.

6.3.8 User management

The user management contains all users. A user can edit their user data or data from the same or a deeper user group. The Operator can only change their own password. The following table shows all parameter values of user management in the same order as they are listed in the Web browser:

Parameter	Meaning
Select or delete user name	All users are offered for selection which the logged in user can change. The user can also delete a user. Users can not delete themselves.
Add user name	A new user can be added here. Please note that a maximum of 9 characters can be used for a user name
Password	A password must be assigned. The password must be entered twice. Please note that a maximum of 9 characters can be used for a password.
User group	Every user must belong to a user group. The user group defines the rights. The current user cannot change his or her user group. Chapter 9.5 External size of device User groups and user interfaces shows all users and their rights.
Automatic login	Automatic Login can be assigned exclusively to one single user. When switching on the device and when changing the user interface (see Chapter 6.1 Mutual exclusion), the login dialogue in this case is skipped and the user with Automatic Login is logged in.

Only one user can have automatic login. If a user has automatic login, this user loses it if the automatic login is assigned to another user.

If automatic login is activated, automatic log out is deactivated.

Changing between automatic login and automatic logout requires a reboot.

6.3.9 Log

For general information about logs, please read Chapter 9.7 Log

Through the Info menu logs can be downloaded to the computer.

Selection	Meaning
Weld (CSV)	Weld log as CSV file
Weld (HTML)	Weld log as HTML file
All	Raw log data

When outputting logs via Web browser, not all data saved in the device is outputted. The most current log is output up to a file size of around 500 kB.

6.3.10 Info

Detailed information can be downloaded in this menu. The following table lists this information:

Function	Meaning
State	The current status of the device can be downloaded and viewed via the Download button.
Download log files	A log file can be downloaded and opened by pressing the corresponding selection button. Further details can be found in Chapter 6.3.9 Log.
Maintenance	Informs you about the abrasion condition of maintaining parts. If the wearing down of one or more parts has reached $\geq 100\%$, an information will occur on the operator screen. Refer to Chapter 5.5 Maintenance.
Date and Time	see Chapter 6.3.7 Date and time.

6.4 Automatic Logout

Login and automatic login is already described in chapter 6.2.9 User management and 6.3.8 User management.

The automatic logout feature is based on CFR21 Part 11. After a certain time without user interaction, the machine logs out automatically. Changes not saved will be discarded. However, there are certain exceptions.

No user interaction is required during a process (e.g. welding). The time starts counting after the weld is completed. Weld log screen or screen for removing the blade do not count to the process. If the timer runs out during one of those screens, the login screen will appear after the last screen (remove blade).

The same situation is by the validation or initialisation. No user interaction is needed as long as the BWTC is doing some work (moving motors or heating). After the work is done (waiting for the user to do something) the logout timer is running. The login screen will appear immediately after time run out or after some ending up user interaction (e.g. removing the blade).

The timer can be adjusted in the web interface under general parameters. If a user has the automatic login, the automatic logout is deactivated.

7 Sequences

Sequences are fully automatic processes in which drives are moved and/or the blade is heated when the cover is closed. A sequence is terminated if the cover of the BioWelder[®] TC opens. Only sequences for which the user has authorization and which the user can execute with the current machine status are available on the user interface. For example, an Operator cannot start a validation or the welding is possible only after an initialisation is successfully run. The following chapter explains the special features of the various sequences.

7.1 Initialise axes

All axes are initialised in this sequence. During initialisation all component systems are tested to ensure greater process safety in the welding process.

The device decides when an initialisation is needed and prompts the user to do it (e.g. after the machine is switched on).

Closing the cover starts the initialisation, and it can be stopped at any time. The initialisation sequence is aborted. Initialisation can be restarted at any time.

7.2 Weld

The welding process is described in Chapter 4.4 Function. Read Chapter 5.9 Welding to find out how a welding process is started.

A welding process can be aborted by the user through the Stop button. It is impossible to stop the sequence during the critical section. During this section the Stop button in the operator control disappears.

If the user aborts a welding process through the Stop button, the device automatically reverts to the starting position. Read more about it in Chapter 7.3 Reverse.

7.2.1 Critical section

The critical section of a welding process is that range in which the automatic sequence of the welding process should in no way be discontinued. The critical section begins when the tubes are separated by the blade and ends after they are welded.

The user cannot abort the critical section of the welding process through the Stop button. An abort during this section is possible only when there is network failure. Such a network failure must be avoided. Therefore, never switch off the device. We recommend the use of an uninterrupted power supply (UPS).

A network failure during the critical section of the welding process must absolutely be avoided.

If the welding process was nonetheless discontinued during this section because of a network failure or another problem, try to minimize the damage by means of a subsequent rewelding or reverse. For more about this, see Chapter 7.2.2 Abort in the critical section.


7.2.2 Abort in the critical section

Basically, there are two ways to abort the welding process in the critical section. If the operator control shows the "Power has failed" message, a network interruption has discontinued the welding process in the critical section. If the operator control shows the "Failed in critical section" message, the welding process was aborted in the critical section because of a drive blockage.

Sartorius Stedim Biotech recommends that a rewelding be selected to minimize the risk of failure. When in the critical section while using liquid filled tubing welding, select "Reweld" only to avoid liquid on the device.

Follow the operating instructions. Proceed according to the following diagram:

The blade is brought to the starting position with „Device reverted from the critical section“. Both welding axes remain closed so that fluid can not flow out.

Clamp the tubes very close to the housing but outside the cover. This minimizes damage in case of failure.

7.2.3 Reweld

Rewelding can be executed only after the welding process is aborted in the critical section. It is the attempt to limit the ensuing damage and should prevent the medium from leaking out the tubes. If rewelding is selected, the user is given the possibility before welding to clamp the tubes and to choose the right type of tube in the settings.

The user must check the quality of welding after a welding. The welding must be repeated if needed. A risk assessment regarding the sterility of the weld has to be performed by the user.

Make sure that the rewelding is not discontinued. If this occurs, no reweld is possible anymore and the device has to be fixed by a technician.

7.3 Reverse

During a reverse all drives are brought back to the starting position and the blade cooled off to the release temperature.

This sequence is started if the user has pressed the Stop button, or after a power failure, or if a drive or the heating fails to reach the default setting.

After the reverse of the device due to an abort of the welding process in the critical section, there is risk of weld failure and damage to the device. A risk assessment regarding the sterility of the weld has to be performed by the user.

7.4 Temperature Validation

To validate is to check the built-in temperature sensor using an external measuring system. The validation set consisting of a type K thermocouple and a holder was designed for this purpose. Validation helps boost the reproducibility of welds. The validation process consists of several stages. These stages are:

1. Choose temperature Validation menu

To be able to execution a validation, an initialised device is required. The Validate sequence is called up via the corresponding symbol in the menu. Only the Supervisor user group has access to the validation icon in the menu.

2. Install blade

Close the cover after the blade is inserted. The device now runs with the drives to the positions for validation. The cover then opens automatically.

3. Click in validation set

The validator sensor is clicked into the holder provided for it, see Chapter 5.2 Familiarizing yourself with the device. Once the sensor is locked into place, close the cover and press the Start button.

**The cable of the temperature sensor must be free to move.
Use only a calibrated gauge for the temperature measurement.**

An automatic sequence controls three predefined blade temperatures for 10 seconds each. The predefined temperatures are displayed on the touch screen as soon as they are reached. The validation is finished when the cover opens.

4. Remove validator set

The sensor can now be removed from the holder. Close the cover. The drives are brought back to the starting position.

5. Remove blade

To finish the process of validation, you must remove the blade with the blade remover tool and then exit the Validate menu by pressing the Cancel button.

If the user executes an abort, please follow the instructions on the operator control.

8 Operational problems

8.1 Switching off device

The user must not switch off the BioWelder® TC device at any time. Make sure that the device remains switched on during an on-going process.

Before the user switches off the device, the on-going process must be aborted via the operator control. For more about this, see Chapter 7.2.1 Critical section

8.2 Foreign objects in the device

The BioWelder® TC is a fully automated device with several electromechanical drives which sometimes generate very high forces.

It is imperative that you ensure that no foreign objects accidentally reach the mechanical parts. This can lead to a blockage or even a defect in the device.

8.3 Wrong type of tube inserted

The type of tube installed is not automatically recognized by the device and must be set through the operator control. Only the Supervisor user group can change the type of tube. Therefore, there is always the possibility that a user could insert the wrong type of tube.

The consequence of this can be that the welding quality is insufficient or the welding is aborted.

Make sure that you only insert the type of tube displayed on the operator control into the tube holder.

There are special welding parameters for pre-treated tubes (irradiated, autoclaved, etc.).

If the wrong type of tube is chosen, a failure may occur.

8.4 Wrong tube dimension in the tube holder

The tube dimensions are not automatically recognized by the device. Therefore, there is always the possibility that a user could insert a tube with the wrong dimensions.

The consequence of this can be that the welding quality is insufficient or the welding is aborted.

The tube holders are labelled with the inside and outside diameter of the tube dimension. Only tube types with identical dimensions must be inserted. Dimensions are indicated in inches.

8.5 Memory card

The SD memory card is not recognized by the device. Possible causes: The card is defective, or was not FAT formatted.

8.6 Device reverted

See Chapter 7.3 Reverse.

8.7 Network

The current IP address and the setting of the DHCP client of the BioWelder® TC can be queried from the Info menu item. See Chapter 6.2.11 Info.

The network configuration of the BioWelder® TC can be reset to the factory settings via the Reset Network Parameters menu. Read more about this in Chapter 6.2.8 Network parameters for reset.

DHCP activated

Use DHCP only if there is a DHCP server in your network. Make sure that the network configuration was done correctly.

If no connection can be established with the Web server of the device via the network and the settings are correct, it is possible that your network has a device with the same MAC address as the BioWelder® TC. To resolve this conflict, change the MAC address of the BioWelder® TC so that it is unique.

DHCP deactivated

Make sure that the network to which you connect the BioWelder® TC has no second device with the same IP address. The BioWelder® TC device is delivered by Sartorius Stedim Biotech with the static IP address of 192.168.1.199. There are various solution approaches to resolve this conflict:

Network Parameter Reset	The network parameters of the device can be set to the original state via the Reset Network Parameters menu. If the network already has a device with the IP address 192.168.1.199, the IP address of one of the two clients must be changed.
Isolated network	The network settings of the BioWelder® TC can be changed only via the operator control via Web browser. You might need to switch the device in an isolated network to a unique IP address.

The BioWelder® TC device can be sent a ping to see if it is present in the network.

8.8 Web browser settings

To work with the Web browser, Java Scripts and cookies must be allowed. Make sure as well that the Automatic or Western Europe setting is chosen under Encoding, otherwise characters may be incorrectly displayed.

8.9 Password forgotten

There are two ways to receive a forgotten password. A user can ask another user from the Supervisor user group to reset the password or the emergency input is used. The technician in charge must request the login data for the emergency input.

Users from the Operator user group can change only their own password.

9 Appendix

9.1 Copyright

Parts of the device are patented. No part this device may be reproduced, distributed, translated or transferred, in any form or through any means, electronically or mechanically, including photocopying, recording or saving to an information storage and retrieval system, without the previous written approval of Sartorius Stedim Biotech GmbH.

9.2 Disposal

The manufacturer disposes of the BioWelder® TC.

The device must be decontaminated before being sent back to Sartorius Stedim Biotech.

9.3 Technical data


Type designation	BioWelder® TC, BWTC
Power connection	100 VAC – 240 VAC
Input frequency	50/60 Hz
Power input	300 VA
In and out connections	Device plug C14 max. 250VAC Ethernet jack type RJ45
Power connection of fuse	2 × 3.15 A T (Type FST)
Battery	CR2032
Operating temperature	+5°C – +40°C *
Place of use	Indoor (Laboratory)
Transient overvoltage	Overvoltage category II
Pollution degree	2
Altitude	up to 2000 m
Humidity	80% up to 31°C, linearly diminishing to 50% relative humidity at 40°C, not condensing
Degree of protection	IP20
Weight	16.4 kg
External size (L × W × H)	555 mm × 261 mm × 269 mm
Power cord	According to local regulations minimum 3 × AWG18 or 3 × 0.75 mm ² minimum local mains supply voltage

* The device is programmed with standard parameter sets for welding C-Flex® 374, AdvantaFlex®, PharMed® BPT and SaniPure™ BDF™. These parameter sets have been validated at room temperature.

9.4 Technical contact

Sartorius Stedim Switzerland AG
Ringstrasse 24a
CH-8317 Tagelswangen
+41 52 354 36 36
biotech.switzerland@sartorius-stedim.com

9.5 External size of device

9.6 User groups and user interfaces

The following table shows the authorizations of user groups and the user interfaces.

The two user groups are indicated in the User Group column. Operator group has lower-level authorizations and Supervisor group has higher-level authorizations. The Operator Control column indicates which user interface supports which function.

Whether a function is usable depends on the machine status. Welding, for example, is impossible if no initialisation is done or if the device is blocked. These limitations are not visible here.

Example: A user from the Supervisor group can validate, but a user from the Operator group cannot.

User group	Operator control	Function		Comments	
		Group or menu item	Detail		
Operator	Touch screen and Web browser	Initialise axis			
		Weld			
		Log in			
		Log out			
	Parameter	User management	Only the user's own password, no assignment by auto login		
	Touch screen	Info			
		SD card			
Supervisor	Touch screen and Web browser	Validate			
		Parameter	Type of tube		
			Language		
			User management	No restriction for the processing of users from the user's own group or those of a lower level	
			DHCP		
			MAC address		
	Web browser	Static IP address			
		Static subnet mask			
		Static default gateway address			
		Static server address			
		Info			
		Touch screen	Time		
			Date		
Network reset					

9.7 Log

BioWelder[®] TC records operations for its on-going operation and saves them. A circular buffer is provided for this purpose. If the memory is full and a new log appears, the oldest log is deleted permanently and the new log is saved. Depending on setting and operating mode, it is possible that recordings are never overwritten even if the device has executed thousands of welds.

These recordings can be copied to a storage medium through the Web browser or through the touch screen. See Chapter 6.2.10 Log for touch screen or Chapter 6.3.9 Log for the Web browser.

The output format of the weld log is CSV or html. To read a csv file, various tabular calculation programs such as Microsoft Excel or OpenOffice Calc can open these files without problem or import functions for databases can be written. To open a html file, use a web browser such as Firefox, Chrome or Internet Explorer.

The html file has a self-validation function built in. This function prevents accidentally altering of the data stored inside the html log file. At the end of the file it shows either "Valid" in green or "not Valid" in red.

A criminal action cannot be prevented with this self-validation function!

The internal memory of the BioWelder[®] TC has a capacity of around 50'000 Weld Logs. If the memory is full, the oldest entries will be overwritten without warning.

The customer is responsible to archive the downloaded weld logs.

The following table lists all logs and their content which the device outputs:

Log	Content
Weld log	<p>The weld log is always output in the language that is currently chosen.</p> <p>The weld log is a table with the following columns: Date (dd.mm.yyyy) Time (hh:mm:ss) User (Name of the user who performed the welding) Tube type (for example: C-Flex 374) Tube description (for example: OD 1in / ID 3/4in) Tube holder identification Cycle (A counter that is incremented after every welding.) Depyrogeneration temperature Welding temperature Welding time (Cutting in to the tubes – tubes are pressed together) Process time (Pressing start – Fan starting) Event (Important event during the welding such as abort by user, blocked condition, etc.)</p>
Raw log	Please always send this file with a problem report/complaint.
Status	Please always send this file with a problem report/complaint.

9.8 Tube holder identification


The following list shows all currently supported tube holder identifications and their tube geometries.

Tube holder identification color code	Tube geometry			
	Inside diameter		Outside diameter	
	[mm]	[inches]	[mm]	[inches]
1 blue	19.0	3/4	25.4	1
2 green	15.9	5/8	22.2	7/8
3 grey	12.7	1/2	19.0	3/4
4 white	9.5	3/8	15.9	5/8
12 red	6.4	1/4	11.1	7/16
9 orange	6.4	1/4	9.5	3/8
11 yellow	3.2	1/8	6.4	1/4

9.9 Declaration of conformity

The BioWelder® TC is compliant with the following certifications:

- CE
- TÜV SÜD NRTL
- REACH WEEE
- RoHS

Note: The BioWelder® TC is not compliant with the CFR21 part 11.

9.10 Error reports and Warning messages

The following tables list all error reports and warning messages which can appear in the operator control.

Error report	Possible cause	Remedy
All runs blocked (X)	No. Error	see Chapter 4.5 Blocked device
(X) = Error number	0 Turning axis	
	1 Turning squeeze axis	
	2 Fixed squeeze axis	
	3 Cutting axis	
	4 Forward axis	
	5 Heating*	
	6 Fan	
	7 Placeholder	
	8 Placeholder	
	9 Cover lock	
	Run stopped by power fail	Follow the instructions on the user interface.

* The (5) Heating Error is triggered when the blade temperature has not reached the depyrogenisation temperature of < 400°C

Warning messages	Possible cause	Elimination
Please remove tube holder	Tube holder inserted	Tube holder removed
No tube holder inserted	1. No tube holder inserted 2. Contact surface of tube holder or the spring contacts are dirty.	1. Insert tube holder set 2. Clean the contact surfaces of the tube holders and those of the spring contacts with a suitable cleansing agent
No matching tube was found	No welding parameters available for the selected type of tube and the inserted tube holder set.	Was the right type of tube chosen? Request new welding parameters from Sartorius Stedim Biotech.
Only one tube holder inserted	1. Only one tube holder is inserted 2. Only one tube holder is recognized.	1. Insert both tube holders. 2. Clean the contact surfaces of the tube holders and those of the spring contacts with a suitable cleansing agent
Tube holders do not match	Tube holders with different color marking were inserted.	Insert tube holders with the same color marking.
Remove blade	Blade has not yet been removed up to now	Remove blade
Blade not inserted	Blade not completely inserted to back of holder	Insert blade

Warning messages	Cause	Elimination
There are no valid parameters on this machine. Please load parameters.	There are no parameters on the device.	Switch machine on and off. If there is still no current data, contact your Sartorius Stedim Biotech Service Representative.
SD missing. Please insert SD and retry!	No SD card inserted	Insert card in SD slot and repeat the action again.
SD not valid	Wrong SD card inserted	Use the SD card supplied by Sartorius Stedim Biotech.
Login failed	Wrong user name or password entered	Enter correct user name with the associated password
Insert validation holder	<ol style="list-style-type: none"> 1. Validation holder not inserted 2. Validation holder inserted, identification defective 	<ol style="list-style-type: none"> 1. Insert validation holder 2. Contact Sartorius Stedim Biotech Service Representative.
Remove validation holder	<ol style="list-style-type: none"> 1. Validation holder was not removed 2. Validation holder was removed, identification defective 	<ol style="list-style-type: none"> 1. Remove validation holder 2. Contact your Sartorius Stedim Biotech Service Representative.
Run stopped by user	User interrupts a sequence	Confirm message with OK
Log transfer stopped by user	Abort by the user while data being written to the SD card	Confirm message with OK

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The screenshot shows the BioWelder TC web interface. The top header contains the product name 'BioWelder TC' and the Sartorius Stedim logo. On the left side, there is a vertical menu with buttons for 'Login/-out', 'Initialize Axes', 'Weld', 'Temperature Validation', 'General Parameters', 'User Management', and 'Info'. The 'Login/-out' button is selected. The main content area is titled 'Login/-out' and displays the message 'User 'Supervis' of group 'Supervisor' is logged in'. Below this, it prompts the user to 'Type or choose user name and password'. There is a form with a 'User name' dropdown menu set to 'Supervis', an empty 'Password' input field, and 'Login' and 'Logout' buttons.

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Log-In/Out

User 'Superuser' of group 'Superuser' is logged in

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User name	Password		
Pitsch	<input type="password"/>	Login	Logout

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The screenshot shows the BioWelder TC web interface. The header includes the product name 'BioWelder TC' and the Sartorius Stedim logo. On the left, a vertical menu contains buttons for 'Login/-out', 'Initialize Axes', 'Weld', 'Temperature Validation', 'General Parameters', 'User Management', and 'Info'. The 'Login/-out' button is selected. The main content area is titled 'Login/-out' and displays the message 'User 'Supervis' of group 'Supervisor' is logged in'. Below this, it prompts the user to 'Type or choose user name and password'. A form contains a 'User name' dropdown menu with 'Supervis' selected, a 'Password' input field, and 'Login' and 'Logout' buttons.

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The screenshot shows the BioWelder TC web interface with the 'General Parameters' menu selected. The left sidebar menu is the same as in the previous screenshot, but 'General Parameters' is highlighted. The main content area is titled 'General Parameters' and contains a table of settings:

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The screenshot shows the BioWelder TC web interface. The header includes the product name "BioWelder TC" and the Sartorius Stedim logo. On the left side, there is a vertical menu with buttons for "Login/-out", "Initialize Axes", "Weld", "Temperature Validation", "General Parameters", "User Management", and "Info". The "Login/-out" button is selected and highlighted. The main content area is titled "Login/-out" and displays the message "User 'Supervis' of group 'Supervisor' is logged in". Below this message, it says "Type or choose user name and password". There is a form with two input fields: "User name" (containing "Supervis" with a dropdown arrow) and "Password" (empty). To the right of these fields are "Login" and "Logout" buttons.

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The screenshot shows the BioWelder TC web interface with the "General Parameters" menu selected. The header is the same as in the previous screenshot. The left menu has "General Parameters" highlighted. The main content area is titled "General Parameters" and contains a table of settings:

Tube type	Choose tube type	C-Flex 374
Auto logout time	in [minutes]	11
Language	Choose language	English
DHCP	Click to enable	<input checked="" type="checkbox"/>
MAC address	Type MAC address	00:0E:3A:10:05:F5
Static IP address	Type address	192.168.1.199
Static subnet mask	Type subnet mask	255.255.255.0
Static default gateway address	Type address	192.168.1.1
Static server address	Type address	192.168.70.1
Weld Log Screen	Click to enable	<input type="checkbox"/>

At the bottom of the settings table are "OK" and "Cancel" buttons.

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Phone: +49.551.308.0
www.sartorius.com

The information and figures contained in these instructions correspond to the version date specified below.

Sartorius reserves the right to make changes to the technology, features, specifications and design of the equipment without notice.

Masculine or feminine forms are used to facilitate legibility in these instructions and always simultaneously denote the other gender as well.

Copyright notice:

This instruction manual, including all of its components, is protected by copyright.

Any use beyond the limits of the copyright law is not permitted without our approval.

This applies in particular to reprinting, translation and editing irrespective of the type of media used.

© Sartorius Germany

Last updated:
10 | 2019

C-Flex[®], SaniPure[™] and PharMed[®] are registered trademarks of Saint-Gobain Performance Plastics Corporation.