Virosart® CPV Minisart®

Directions for use

1. Introduction

The purpose of this document is to provide the appropriate method for identifying, installing and using Sartorius Stedim Biotech Virosart® CPV Minisart®. The procedure should be strictly followed to ensure operator safety and to maintain optimum performance of Virosart® CPV Minisart®.

This document pertains to Sartorius Stedim Biotech Virosart® CPV filter elements with the following part number:

Virosart® CPV Minisart®:

- 545VM-----B (box of 5)
- 545VM-----ITA (box of 4)

Virosart® CPV Minisart® has a female luer lock inlet & male luer slip outlet connection.

2. Labeling of Virosart® CPV Minisart®

The device consists of a hydrophilic polyethersulphone (PES) double layer membrane with 5 cm² effective filtration surface within a polycarbonate housing.

The type, lot number and serial number are imprinted on the outer support of the housing of the Virosart® CPV Minisart®.


The Virosart® CPV Minisart® tested for integrity by diffusion testing in house at Sartorius Stedim Biotech are single packed into a blister package which is labeled with type, order code, lot number and use before date. In addition, these specifications are printed on the labels of the cardboard box.

3. Installing of Virosart® CPV Minisart®

⚠ The Virosart® CPV Minisart® tested for integrity in house at SSB (545VM-----ITA) is sterilized in a validated production step by ETO sterilization according to EN ISO 11135. Sterility, however, is only guaranteed when this primary packaging of the Minisart® is intact before using the Virosart® CPV Minisart® as a sterile product. If the packaging is damaged in any way the Virosart® CPV Minisart® cannot be handled as a sterile product!

If the Virosart® CPV Minisart® is to be used as a sterile product, care must be taken to prevent contamination when connecting the Minisart® to an appropriate filtration system (e.g. connection within a clean bench).

The Virosart® CPV Minisart® must be connected in the direction of filtration at the female luer lock (printed side) to an appropriate filtration set-up (see picture below). An application specialist from Sartorius Stedim Biotech can provide assistance in installing the system at your facility to make sure that requirements for a correct scale down trial are met. Virosart® CPV Minisart® can be connected to the liquid reservoir by using a sterilized disposable 2-way valve. The order code for such valves is 7ZV--K0112, the order code for the stainless steel liquid reservoir is 16249-----6.

Picture 1: Typical filtration set-up for a Virosart® CPV Minisart®

4. Wetting of Virosart® CPV Minisart®

(In the direction of filtration)

Before operation the Virosart® CPV Minisart® as well as the entire system needs to be rinsed with buffer or WFI to ensure that it is particle-free as well as for complete wetting of the membrane and proper use of the whole filter capacity. The Virosart® CPV Minisart® must be wetted with buffer or WFI prior use or autoclaving. Please connect the Minisart® to the filtration system. Sartorius Stedim Biotech strongly recommends using constant pressure mode for operating the Virosart® CPV Minisart®. Carefully increase the inlet pressure to 2 bar | 29 psi and flush the filter approx. 15 minutes with buffer or WFI. Please make sure that no air is entrapped within the Minisart®.

Air bubbles within the Minisart® may decrease the effective filtration area of the device. Typical WFI flow rates may vary between 0.8–1.4 ml/min at 2 bar pressure. If the Virosart® CPV Minisart® has been flushed according to the instructions it is thoroughly wetted. Now carefully remove the buffer or WFI from the reservoir and replace it with product. Please ensure there is no air trapped upstream of the Minisart® and start the filtration at the desired filtration pressure.

5. Autoclaving Virosart® CPV Minisart®

Sterilization temperature: up to 121 °C, 30 min.

Prior to autoclaving, the Minisart® must be wetted according to the directions in the paragraph "Wetting of Virosart® CPV Minisart®". For optimum wetting after autoclaving and to avoid air trapped inside the housing the Minisart® is placed into a beaker of water and autoclaved submerged in water. A liquid autoclaving cycle should be used as the Minisart® device is filled with WFI or buffer.

Once the sterilization temperature has been attained, the Minisart® can be autoclaved for 30 minutes at a maximum sterilization temperature of 121 °C, 1 bar | 14.5 psi. If the Minisart® is not placed into a beaker of water the inlet and outlet of the Virosart® CPV Minisart® should be wrapped with steam permeable autoclave paper or aluminum foil. It is important not to close the Minisart® with plugs (or a pre-filter) to let the water expand during heating up. During the autoclave cycle, please make sure that any connecting piece on the Minisart® (e.g. made of stainless steel) are not attached too tightly and do not have a tractive effect on the Minisart® housing which could cause deformation.

The cooling phase has to be conducted without deliberately decreasing the pressure (opening the autoclave at temperatures above 85°C), because pressure drops at elevated temperature may lead to bumping and device damages.

Autoclaved Minisart® must be flushed with buffer for 15 min @ 2 bar | 29 psi and used for filtration within 4 h after cooling down to room temperature. With longer storage times after autoclaving there is a possibility of trapping air in the housing. Air in the housing would lead to less effective filtration area and with that misleading filtration results.

Virosart® CPV Minisart® can be sterilized by autoclaving up to 2 cycles.

6. Maximal Allowable Differential Pressure

Virosart® CPV Minisart®

In the direction of filtration:
At 20 °C max. 5.0 bar | 73 psi

7. Filtration with Virosart® CPV Minisart®

Virus filtration with Virosart® CPV should be realized at the end of the purification process prior to the last buffer exchange | form and fill. The protein solution which is subject to the virus filtration should be as much as possible free of aggregates, host cell protein & DNA to avoid premature blockage of Virosart® CPV filter. Even if these contaminants should be removed either during the initial recovery process or latest during the polishing process of your target molecule, small amounts might be sufficient to cause premature blockage of the final virus filter. To prevent this, an efficient adsorptive pre-filtration step, such as the Virosart® Max, might be required as protection for the Virosart® CPV. The optimum pre-filter – final filter ratio has to be identified during development of this process step. It is recommended to use fresh and representative material for these studies.

To make full use of the filtration area, make sure that no air is entrapped within the Virosart® CPV Minisart®. The pre-filter should be used in-line with the final virus filter under the same conditions. The filtration can be realized at any pressure up to 5 bar | 73 psi pressure. The filter qualification was performed at an operating pressure of 2 bar | 29 psi to meet customer standard requirements. The flow rate achieved through the Virosart® CPV Minisart® during the filtration is a function of the purity of the feed stream, the target molecule size and concentration. The flow rate will decline with the progress of filtration and indicate the respective blockage of the filter.

After filtration the Virosart® CPV Minisart® can be flushed with buffer to maximize protein transmission. Therefore please flush the Minisart® with >3 ml of buffer. Make sure that no air is entrapped within the Virosart® CPV Minisart® at this step.

8. Post Use Integrity Test for Virosart® CPV Minisart®

After filtration, an integrity test can be performed. Therefore, please submerge the Minisart® still connected to the filtration system, into a beaker of WFI or WFI equivalent. Now please pressurize the Minisart® carefully up to 4.5 bar | 65 psi.

The specification for a passed integrity test is that no air bubbles may be detected at the Minisart® outlet for 30 seconds.

9. Changing out the Virosart® CPV Minisart®

The Virosart® CPV Minisart® must be changed, at the latest, when the maximum allowable differential pressure is reached.

10. Applications Support by Sartorius Stedim Biotech

If required, please contact your local Sartorius Stedim Biotech representative to obtain further information on technical data, product integrity testing or general information concerning specific applications. Our application specialists from Purification Technologies will gladly assist in your filtration trial work if needed. Moreover our validation service and support team are always available to assist with questions concerning product and process validation, as well as training of staff and optimization of filtration processes.

11. Return of used filter elements

If you wish to return used filter elements to Sartorius Stedim Biotech for inspection, it is necessary to decontaminate and sterilize them prior to shipment. This procedure has to be documented in the return shipment form, available from your local Sartorius Stedim Biotech representative. This documentation is required by country law. Returned filter elements cannot be handled unless the appropriate return shipment form is completed.

12. Liability

Sartorius Stedim Biotech cannot assume liability if Virosart® CPV Minisart® are subject to improper use. In the interest of product development we reserve the right to make changes.

Thank you for working with Sartorius Stedim Biotech, we appreciate your business.

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Sartorius Stedim Biotech GmbH