Virosart® HC 539... Standard Filter Cartridges & T-Style MaxiCaps®

Directions for Use

1. Summary

The purpose of this document is to provide the appropriate method for identifying, installing and using a Sartorius Stedim Biotech Virosart® HC standard filter cartridge and Virosart® HC T-Style MaxiCap®. The procedure should be strictly followed to ensure operator safety and to maintain optimum performance of the Virosart® HC filter elements.

This document pertains to Sartorius Stedim Biotech Virosart® HC filter elements with the following part numbers:

Virosart® HC Filter Cartridges:

- 5392528V1
- 5392528V2
- 5392528V3

Virosart® HC T-Style MaxiCaps®:

- 5398328V1-SS
- 5398328V2-SS
- 5398328V3-SS

All Virosart® HC T-Style MaxiCaps® have sanitary inlet & outlet connections.

2. Labeling

The type, lot number and serial number are imprinted on the outer support of the cartridge. The Virosart® HC T-Style MaxiCaps® are labelled on the banderole of their housings with the same data in addition to the order number. These specifications together with the order number are printed on the labels of the plastic bag and the cardboard box. For Virosart® HC T-Style MaxiCaps® the direction of filtration is indicated by an arrow on the banderole of the capsule.

3. Installing Virosart® HC cartridges | T-Style MaxiCaps®

Virosart® HC filter cartridges should be placed into a filter housing and connected with piping. Cut the plastic bag at the bottom end of the cartridge (O-ring). Leave the open plastic bag on the cartridge for protection until the cartridge has been properly installed. Wet the O-rings with sterile water or with the liquid to be filtered. Holding onto the base of the cartridge, insert the filter cartridge by carefully turning it into place. Once the cartridge is securely in place, remove the plastic bag and close the housing. Virosart® HC T-Style MaxiCaps® are connected to the processing line at the inlet and outlet via triclamp (please do not forget the gaskets!) connections.

4. Wetting of Virosart® HC cartridges | T-Style MaxiCaps®

(In the direction of filtration)

Before every sterilization (steaming and autoclaving) and prior to every integrity test the entire system needs to be rinsed to ensure that it is particle-free. To rinse the filter element, open the vent valve on the housing (T-Style MaxiCap® housing) and fill the housing with WFI until liquid comes out of the vent valve. To ensure that all air has been removed from the housing, keep the vent valve slightly opened for the first minute of flushing. Then close the vent valve and rinse for > 10 minutes at a pressure of 2.0 bar | 29 psi. The housing (T-Style MaxiCap® housing) can be drained by opening the drain valve and using compressed air or nitrogen at a maximal pressure of 0.5 bar | 7 psi. As the membrane of the Virosart® HC is much tighter than a 0.2 µm membrane, it is recommended to drain the housing by pressure to speed up this procedure. After the housing is empty, close all the valves. After the Virosart® HC has been flushed according to the instructions, it is thoroughly wetted.

5. Sterilization

5.1 In-line Steam Sterilization

Only Virosart® HC standard filter cartridges are dedicated for steam sterilization! Virosart® HC T-Style MaxiCaps® cannot be in-line steam sterilized!

Sterilization temperature: up to 130°C, 30 min
 Open the vent, drain and if necessary the test | sampling valves. Slowly start a flow of steam so that the differential pressure does not exceed 0.2 bar | 2.9 psi during the pressure build-up process. Throttle the valves as soon as steam escapes from them. Set the required sterilization steam pressure and temperature at the outlet by adjusting the various valves (steam, vent, test | sampling). After reaching the sterilization temperature, continue to steam the system for 30 minutes, keeping conditions as constant as possible. The maximum steaming temperature of 130°C and the maximum differential pressure of 0.2 bar | 2.9 psi should not be exceeded during the whole steaming process. Upon completion of sterilization, turn off the steam, wait until steam no longer emerges from the outlet, and both inlet and outlet pressure have dropped to zero. If necessary, close the valves on the sterile side of the housing. Now an integrity test can be performed with closed inlet valves.

5.2 Autoclaving

Sterilization temperature: up to 121°C, 30 min
 Once the sterilization temperature has been attained the Virosart® HC filter cartridge | T-Style MaxiCap® can be autoclaved for 30 minutes. The inlet and outlet as well as any open valves of the Virosart® HC T-Style housing should be wrapped with steam permeable autoclave paper or aluminium foil. During the autoclave cycle of the Virosart® HC T-Style MaxiCaps®, make sure that connecting pieces on the Virosart® HC T-Style MaxiCaps® (e.g. made of stainless steel) are not attached too tightly and do not have a torque effect on the Virosart® HC T-Style MaxiCaps® housing which could cause deformation.

6. Integrity Testing

The integrity test of a Virosart® HC filter must be done using an automated integrity test unit, preferably the Sartocheck® 3 or Sartocheck® 4 integrity test unit. When using a Sartorius Stedim Biotech Sartocheck® Integrity Test System, the diffusion tests can be automatically performed, recording the corresponding integrity test data. For additional information on how to perform an automatic integrity test please see our Installation and Operating Instructions for the Sartocheck® 3 or Sartocheck® 4.

Integrity Test Data for Virosart® HC standard filter cartridges and Virosart® HC T-Style MaxiCaps®.

Integrity test	Filter type	Test pressure	Limit
Water based Diffusion test	5392528V1	4.5 bar 65.25 psi	24 ml/min
Water based Diffusion test	5392528V2	4.5 bar 65.25 psi	48 ml/min
Water based Diffusion test	5392528V3	4.5 bar 65.25 psi	72 ml/min
Water based Diffusion test	5398328V1--SS	4.5 bar 65.25 psi	24 ml/min
Water based Diffusion test	5398328V2--SS	4.5 bar 65.25 psi	48 ml/min
Water based Diffusion test	5398328V3--SS	4.5 bar 65.25 psi	72 ml/min

Remark:

All values given above were obtained at room temperature and using water as wetting agent.

7. Maximal Allowable Differential Pressure

Virosart® HC standard filter cartridges and Virosart® HC T-Style MaxiCaps®

In the direction of filtration

At 20°C max. 5.0 bar | 73 psi

At 121°C max. 0.2bar | 2.92 psi (T-Style MaxiCaps®)

At 130°C max. 0.2 bar | 2.92 psi (Cartridges)

In the reversed direction of filtration

At 20°C max. 0.2 bar | 2.92 psi

8. Filtration with Virosart® HC

Virus filtration with Virosart® HC should be realized at the end of the purification process prior to the last buffer exchange | form and fill. The protein solution which is subject to the virus filtration must be free of aggregates, host cell protein & DNA to avoid premature blockage of Virosart® HC filter.

These contaminants should be removed either during the initial recovery process or latest during the polishing process of your target molecule. Standard prefiltration prior the virus filtration employs a 0.2 | 0.1 µm Sartopore® 2 membrane filter.

To make full use of the filtration area, the housing must be vented completely free of air. To do this, briefly open the vent valve until the liquid to be filtered comes out. Then apply the desired filtration pressure. The virus filtration should be realized at 2 bar | 29 psi pressure.

The flow rate achieved through the Virosart® HC filter capsule during the filtration is a function of the purity of the feed stream, the target molecule size and concentration. The flow rate will decline with the progress of filtration and indicate the respective blockage of the filter.

9. Changing the Virosart® HC filter

The Virosart® HC filter must be changed, at the latest, when the maximum allowable differential pressure is reached or when the validated process volume has been processed over the filter.

10. Applications Support

If required, please contact your local Sartorius Stedim Biotech representative to obtain further information on technical data, product integrity testing or general information concerning specific applications. Our application specialists from Purification Technologies will gladly assist in your filtration trial work if needed. Moreover our Validation Service and support team are always available to assist in questions concerning product and process validation, as well as training of staff and optimization of filtration processes.

11. Return of used filter elements

If you wish to return used filter elements to Sartorius Stedim Biotech for inspection, it is necessary to decontaminate and sterilize them prior to shipment. This procedure has to be documented in the Return Shipment Form, available from your local Sartorius Stedim Biotech representative. This documentation is required by country law. Returned filter elements cannot be handled unless the appropriate Return Shipment Form is completed.

12. Liability

Sartorius Stedim Biotech can not assume liability if Virosart® HC filter cartridges and Virosart® HC T-Style Maxi-Caps® are subjected to improper use.

In the interest of product development we reserve the right to make changes.

Thank you for working with Sartorius Stedim Biotech, we appreciate your business.

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Sartorius Stedim Biotech GmbH