



## Virosart® Max\*

Maximize your virus filter performance with adsorptive pre-filtration

### Product Information

Virosart® Max is a specifically optimized virus pre-filter significantly increasing downstream virus filter performance. This filter combines size exclusion mechanism with efficient adsorptive capacities to increase the robustness of the following virus filter. As a result, Virosart® Max ensures the highest protection of your final virus-retentive membrane, significantly increasing its robustness and capacity.



### Description

Choose your perfect fit from the Sartorius virus clearance strategy summarizing orthogonal technologies, manufacturing solutions, validation support and consultancy. The orthogonal technologies from Sartorius consisting of virus adsorption by chromatography, virus inactivation and virus filtration. The Virosart® product range includes four different virus retentive membranes, in order to provide the best solution for every application.

The performance of the final virus filter mainly depends on the feed stream properties. Therefore, appropriate feed stream conditioning i.e. through efficient pre-filtration is a major aspect within the development of a virus filtration step.

### Application & Positioning

The main application for Virosart® Max is the protection of the final virus retentive membrane. For that purpose this filter is used at the end of the purification process in-line with the final virus filtration step of the biopharmaceutical product. The optimized pre-filter – final-filter ratio should be determined during development of the virus filtration step.

### Product Benefits

Virosart® Max provides highest protection to the virus retentive membrane with significant increase in robustness and capacity of the final virus filter. The most challenging molecules for the final virus retentive membrane are aggregates and | or small hydrophobic molecules. Virosart® Max binds aggregates very efficiently through hydrophobic interactions, independently of process conditions such as conductivity. The clean triple layer membrane material provides highest adsorptive capacities optimized for this process step.




With the T-Style design the Virosart® Max MaxiCaps® is ideal for easy installation of multiple filters in series or parallel.

### Integrity Testing

Virosart® Max filters are tested for integrity using a water-based integrity test with the Sartocheck® technology of Sartorius Stedim Biotech.

\* The patented technology (DE 10 2011 105 525 B4) binds aggregates efficiently through hydrophobic interactions with polyamide, independently of process conditions such as conductivity from biological feed streams (mAbs, plasma derivatives or recombinant proteins).

## Technical Data

	Minisart® & Sartolab	Capsule & MidiCaps®	MaxiCaps® & Cartridge
			
<b>Nominal filtration area</b>	5 & 21 cm <sup>2</sup>	220 cm <sup>2</sup> 0.24 ft <sup>2</sup>	0.21 m <sup>2</sup> 2.3 ft <sup>2</sup>
<b>To be used for</b>	<ul style="list-style-type: none"> <li>• Scale-down work</li> <li>• Filtration work for capacity studies</li> <li>• Optimizing the pre-filter final-filter-ratio of the final virus filtration step</li> </ul>	<ul style="list-style-type: none"> <li>• Scale-up studies</li> </ul>	<ul style="list-style-type: none"> <li>• Large scale manufacturing</li> </ul>
<b>Typical filtration volume</b>	< 200 mL	< 80 L	> 80 L
<b>Available connectors</b>	<ul style="list-style-type: none"> <li>• Female Luer Lock inlet &amp; male Luer Lock outlet</li> </ul>	<ul style="list-style-type: none"> <li>• <b>MidiCaps® &amp; Capsule:</b> 3/4" triclamp (sanitary) connector inlet &amp; outlet</li> <li>• <b>Capsule:</b> 1/2" hose barb connector inlet &amp; outlet</li> </ul>	<ul style="list-style-type: none"> <li>• <b>MaxiCaps®:</b> Sanitary inlet &amp; outlet adapter</li> <li>• <b>Cartridge:</b> S-adapter top, 2 flange bayonet adapter with double o-ring bottom</li> </ul>
<b>Sterilization</b>	<ul style="list-style-type: none"> <li>• Autoclaving: 121°C @ 2.0 bar   29 psi for 30 min up to 2 cycles</li> </ul> <p>⚠ No inline steaming of Minisart®</p>	<ul style="list-style-type: none"> <li>• Autoclaving: 121°C @ 2.0 bar   29 psi for 30 min up to 2 cycles</li> </ul> <p>⚠ No inline steaming of Capsule &amp; MidiCaps®</p>	<ul style="list-style-type: none"> <li>• <b>MaxiCaps®:</b> Autoclaving 121°C @ 2.0 bar   29 psi for 30 min up to 2 cycles</li> </ul> <p>⚠ No inline steaming of MaxiCaps®</p> <ul style="list-style-type: none"> <li>• <b>Cartridges:</b> Steaming 121°C @ 2.0 bar   29 psi for 30 min up to 2 cycles</li> </ul>
<b>Operating parameters</b>	<ul style="list-style-type: none"> <li>• In the direction of filtration: max. 5.0 bar   72.5 psi at 20°C; max. 0.5 bar   7.3 psi at 121°C</li> <li>• In the reversed direction of filtration: max. 0.5 bar   7.3 psi at 20°C</li> </ul>		
<b>Water based diffusion test at 2.0 bar   29 psi</b>	N   A	<ul style="list-style-type: none"> <li>• 2 mL/min (220 cm<sup>2</sup>)</li> <li>• 10 mL/min (0.27 m<sup>2</sup>)</li> </ul>	<ul style="list-style-type: none"> <li>• 8 mL/min (0.7 m<sup>2</sup>)</li> <li>• 16 mL/min (1.4 m<sup>2</sup>)</li> <li>• 24 mL/min (2.1 m<sup>2</sup>)</li> </ul>

## Materials

### Device

#### Cartridge, MaxiCaps®, MidiCaps® & Capsule

**Supportive Fleece**  
Polyester

**Core (not capsule)**  
Polypropylene

**End Caps**  
Polypropylene

**Capsule Housing**  
Polypropylene

#### Minisart®

**Housing**  
Polypropylene

#### Sartolab


**Housing**  
Polycarbonate

#### Membrane

**Material**  
Optimized polyamide

**Pore size**  
0.1 µm (nominal)

**Format**  
Triple layer




Construction of Viroart® Max cartridge and capsule with zoom on cross section and membrane

# Performance


## Characteristic Water Flow Rates

The flow rate achieved through the Virosart® Max, in-line with the final virus filter, is determined by the flow rate of the final virus filter!


### Virosart® Max 30" Cartridges & 30" T-Style MaxiCaps® (2.1 m<sup>2</sup> | 22.6 ft<sup>2</sup>)




### Virosart® Max 20" Cartridges & 20" T-Style MaxiCaps® (1.4 m<sup>2</sup> | 15 ft<sup>2</sup>)




### Virosart® Max 10" Cartridges & 10" T-Style MaxiCaps® (0.7 m<sup>2</sup> | 7.5 ft<sup>2</sup>)



### Virosart® Max MidiCaps® (220 cm<sup>2</sup> | 0.24 ft<sup>2</sup>)



### Virosart® Max Capsules (180 cm<sup>2</sup> | 0.19 ft<sup>2</sup>)



## Regulatory Compliance

- Each individual filter is tested for integrity (except 54AMI-----B and 54SA-----V)
- Designed, developed and manufactured in accordance with an ISO 9001 certified Quality Management System
- Meet or exceed the requirements for WFI quality standards set by the current USP
- Non pyrogenic according to USP Bacterial Endotoxins
- USP Plastic Class Test VI

## Technical References

Validation Guide	SPK5802-e
Extractables Guide	SPK5803-e

## Ordering Information



### Cartridge

54A 25 58

#### Adapter

25: S-adapter top,  
1 flange bayonet adapter  
with double O-ring bottom

#### Filter size

N1: 10" 0.7 m<sup>2</sup> | 7.5 ft<sup>2</sup>  
N2: 20" 1.4 m<sup>2</sup> | 15 ft<sup>2</sup>  
N3: 30" 2.1 m<sup>2</sup> | 22.6 ft<sup>2</sup>



### T-Style MaxiCaps®

54A 83 58 -- SS

#### Filter size

N1: 10" 0.7 m<sup>2</sup> | 7.5 ft<sup>2</sup>  
N2: 20" 1.4 m<sup>2</sup> | 15 ft<sup>2</sup>  
N3: 30" 2.1 m<sup>2</sup> | 22.6 ft<sup>2</sup>

#### Adapter

SS: Sanitary inlet &  
outlet adapter



### MidiCaps® & Capsules

54A 58 -- FF --

13: Capsules  
53: MidiCaps®

Filter size  
N4: 220 cm<sup>2</sup> | 0.24 ft<sup>2</sup>  
N9: 0.21 m<sup>2</sup> | 2.3 ft<sup>2</sup>

#### Connector

FF: 1½" Tri-Clamp (sanitary)  
inlet & outlet

#### Units per package

B: Five pieces  
V: Two pieces



### Minisart® & Sartolab

54A -- -- -- -- --

Filter size  
MI: Minisart®  
SA: Sartolab

Units per package  
B: Five pieces  
V: Two pieces

## Accessories & Services

### Virus Filtration

Virosart® HF is a high speed virus filter for mAbs and recombinant proteins. The filter combines highest virus safety with superior capacities.



### Integrity Testing

Fully automated Virosart® integrity testing to guarantee intactness of the Virosart® filter pre- and post diffusion test.



### Ready-to use Filter Transfer Sets

Simplify your daily routine work by using modular filter assembly.




### Single-use Systems

Flexible processing with FlexAct® VR system for production from pilot plants up to commercial processing.



### Customized Systems

High level of automation and individual requirements can be relegalized by customized single-use or hybrid solutions.




### Testing Service

Your partner to assure virus safety for your process by MCB | WCB characterization, bulk harvest testing and spiking studies.



### Services Worldwide

Trust our comprehensive range of services for your virus filtration processes: We gladly assist you with tasks like process validation, process optimization and many more.



Sartorius Stedim Biotech GmbH  
August-Spindler-Strasse 11  
37079 Goettingen, Germany  
Phone +49.551.308.0  
[www.sartorius-stedim.com](http://www.sartorius-stedim.com)

USA Toll-Free +1.800.368.7178  
Argentina +54.11.4721.0505  
Brazil +55.11.4362.8900  
Mexico +52.55.5562.1102  
UK +44.1372.737159  
France +33.442.845600  
Italy +39.055.63.40.41  
Spain +34.913.586.098  
Russian Federation +7.812.327.53.27  
Republic of South Africa +27.11.315.54  
Japan +81.3.4331.4300  
China +86.21.6878.2300

Specifications subject to change without notice.  
Copyright Sartorius Stedim Biotech GmbH. Printed  
in the EU on paper bleached without chlorine.  
Publication No.: SPK2181-e190603  
Order No.: 85037-561-70  
Ver. 06 | 2019