Centrisart® I


Directions for Use

Centrisart® I 100K device for *in vitro* diagnostic use


Centrisart® I 5K, 10K, 20K, and 300K devices for research use only; not for use in diagnostic procedures

Storage conditions | shelf life

Centrisart® I ultrafiltration spin columns should be stored at 15–30°C. The devices should be used before the expiry date printed on the box.

1. Remove the closing cap, allow the floater (inner tube) to slide out and stand it on a clean surface with open end down, avoiding contact with the membrane. Pour the sample (0.1–2.5 ml) into the centrifuge tube (outer tube, Fig. 1).

2. Re-insert the floater, membrane side side down, and let it slide down onto the surface of the sample. Allow the centrifuge tube to stand (Fig. 2) for approx. 5 minutes to ensure complete wetting of the membrane. Exception: to filter noncoagulated blood, precentrifuge the tube for about 2 min. before inserting floater in order to avoid hemolysis. Do not re-insert closing cap for centrifugation.

3. Centrifuge in a laboratory centrifuge accepting up to 17 × 100 mm tubes (Fig. 3). Maximum relative centrifugal force: 2,500 g with swing head rotors, 2,000 g with angular head rotors (see graph). When filtering concentrated protein solutions, for example serum or blood, start with 5 minutes centrifugation at half of the maximum rpm. When filtering solutions containing lipids such as milk, only use an angular head rotor.

4. Pipet out the filtrate (Fig. 4) as soon as possible after completion of centrifugation (filtrate will slowly diffuse back into the concentrate).

5. Or remove the floater with opened forceps to have access to the concentrate (Fig. 5).

Note:
Not for use with organic solvents.

The ultrafilters contain glycerol. Due to its hygroscopicity, transparent spots of darker appearance may develop on the filters at higher levels of humidity. Such spots, however, do not have any adverse effects on the retention rate or on the flow rate of these filters.


If removal of glycerol (approx. 2 mg) is required, it can be effected by a 5-minute dip in water or other suitable liquid, or by precentrifugation, immediately prior to use. Do not allow membrane to dry out. Replace closing cap if losses by evaporation are to be avoided.

For the purpose of product development we reserve the right to make modifications.

CE

The Centrisart® I product line includes 5 different cutoffs (Molecular Weight Cutoff, MWCO):
– Centrisart® I 5K device: 5,000 MWCO
– Centrisart® I 10K device: 10,000 MWCO
– Centrisart® I 20K device: 20,000 MWCO
– Centrisart® I 100K device: 100,000 MWCO
– Centrisart® I 300K device: 300,000 MWCO

Centrisart® I 100K filter devices are for *in vitro* diagnostic use and can be used to concentrate serum, urine, cerebrospinal fluid, and other body fluids prior to analysis. Centrisart® I 5K, 10K, 20K and 300K filter devices are for research use only and not for use in diagnostic procedures. The Centrisart® I devices are supplied non-sterile and are for single use only.

Maximum rpm for each rotor radius (r) between 10 and 20 cm.
Upper curve A: for swing head rotors (2,500 g)
Lower curve B: for angular head rotors (2,000 g)

In Vitro Diagnostic Product Labeling
The following table defines the symbols found on Centrisart® I 100K device labels.

Symbol	Definition
	In vitro diagnostic medical device
	Catalogue number
	Do not reuse
	Use by
	Batch code

Symbol	Definition
	Date of manufacture
	Manufacturer
	Temperature limitation
	Non-sterile product
	CE conformity marking

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Specifications subject to change without notice. Printed and copyrighted by Sartorius Stedim Biotech GmbH
Sartorius Stedim Lab Ltd
Sartorius Stedim North America Inc.