Virosart® HC Capsule

Directions for use

1. Summary

The purpose of this document is to provide the appropriate method for identifying, installing and using Sartorius Stedim Biotech Virosart® HC capsules. The procedure should be strictly followed to ensure operator safety and to maintain optimum performance of Virosart® HC capsules.

This document pertains to Sartorius Stedim Biotech Virosart® HC filter elements with the following part number:

Virosart® HC Capsules (box of 5):
5391328V4--FF--B

All Virosart® HC capsules have sanitary inlet & outlet connections.

2. Labeling of Virosart® HC Capsules

The type, lot number and serial number are imprinted on the outer support of the housing of the Virosart® HC capsule. This information is identical to the information found on the label of both the plastic bag and the cardboard box where in addition to type and lot number the order number is indicated. The direction of filtration is indicated by an arrow on the capsule housing.

3. Installing Virosart® HC Capsules

The plastic bag of the Virosart® HC capsule should be cut open at one end. Leave the open plastic bag on the filter for protection as long as possible. Virosart® HC capsules are connected to the processing line at the inlet and outlet via triclamp (please do not forget the gaskets!) connections.

4. Wetting of Virosart® HC Capsules

(In the direction of filtration)

Before every sterilization (autoclaving) and prior to every integrity test the entire system needs to be rinsed to ensure that it is particle-free. To rinse the filter element, open the vent valve on the capsule housing and fill the housing with WFI until liquid comes out of the vent valve. To ensure that all air has been removed from the housing, keep the vent valve slightly opened for the first minute of flushing.

Then close the vent valve and rinse for >10 minutes at a pressure of 2.0 bar | 29 psi.

The housing of the Virosart® HC capsule can be drained by opening the drain valve and using compressed air or nitrogen with a maximal pressure of 0.5 bar | 7 psi in the direction of filtration.

As the membrane of the Virosart® HC is much tighter than a 0.2 µm membrane, it is recommended to drain the housing by pressure to speed up this procedure. After the housing is empty, close all the valves. After the Virosart® HC has been flushed according to the instructions, it is thoroughly wetted.

5. Sterilization of Virosart® HC Capsules

5.1 Autoclaving

Sterilization temperature: up to 121 °C, 30 min

The inlet and outlet as well as any open valves of the Virosart® HC capsule housing should be wrapped with steam permeable autoclave paper or aluminium foil. Once the sterilization temperature has been attained the Virosart® HC capsule can be autoclaved for 30 minutes. During the autoclave cycle of the Virosart® HC capsule, make sure that connecting pieces on the Virosart® HC capsule (e.g. made of stainless steel) are not attached too tightly and do not have a torque effect on the Virosart® HC capsule housing which could cause deformation.

5.2 In-line Steam Sterilization

No in-line steam sterilization of Virosart® HC capsule!

6. Integrity Testing

It is important to follow the wetting procedure for Virosart® HC to make sure that the filter is wetted out completely prior to the integrity test. The integrity test of a Virosart® HC capsule must be done using an automated integrity test unit, preferable the Sartocheck® integrity test unit. When using a Sartorius Stedim Biotech Sartocheck® integrity test system, the diffusion tests can be automatically performed, recording the corresponding integrity test data. To ensure reliable integrity testing, it is recommended to use 10 minutes stabilisation time and 5 minutes measuring time. For additional information on how to perform an automatic integrity test please see our installation and operating instructions for the Sartocheck® integrity test unit.

Integrity Test Data for Virosart® HC Capsules

Integrity test	Filter type	Test pressure	Limit
Water based diffusion test	5391328V4	4.5 bar 65.2 psi	2 ml/min

Remark

All values given above were obtained at room temperature and using water as wetting agent.

7. Maximum Allowable Differential Pressure

Virosart® HC Capsules

In the direction of filtration

At 20 °C max. 5.0 bar | 73 psi

At 121 °C max. 0.2 bar | 2.92 psi

In the reversed direction of filtration

At 20 °C max. 0.2 bar | 2.92 psi

8. Filtration with Virosart® HC Capsules

Virus filtration with Virosart® HC should be realised at the end of the purification process prior to the last buffer exchange | form and fill. The protein solution which is subject to the virus filtration must be free of aggregates, host cell protein & DNA to avoid premature blockage of Virosart® HC capsules.

These contaminants should be removed either within the initial recovery process or latest by the polishing process of your target molecule. Standard pre-filtration prior the virus filtration employs a 0.2/0.1µm Sartopore 2 membrane filter.

To make full use of the filtration area, the housing must be vented completely free of air. To do this, briefly open the vent valve until the liquid to be filtered comes out. Then apply the desired filtration pressure. The virus filtration should be realised at 2 bar | 29 psi pressure.

The flow rate achieved through the Virosart® HC capsule during the filtration is a function of the purity of the feed stream, the target molecule size and concentration. The flow rate will decline with the progress of filtration and indicate the respective blockage of the filter.

9. Changing out the Virosart® HC Capsules

The Virosart® HC capsule must be changed, at the latest, when the maximum allowable differential pressure is reached or when the validated process volume has been processed over the filter.

10. Applications Support

If required, please contact your local Sartorius Stedim Biotech representative to obtain further information on technical data, product integrity testing or general information concerning specific applications. Our application specialists from Purification Technologies will gladly assist in your filtration trial work if needed. Moreover our validation service and support team are always available to assist in questions concerning product and process validation, as well as training of staff and optimization of filtration processes.

11. Return of used Filter Elements

If you wish to return used filter elements to Sartorius Stedim Biotech for inspection, it is necessary to decontaminate and sterilize them prior to shipment. This procedure has to be documented in the return shipment form, available from your local Sartorius Stedim Biotech representative. This documentation is required by country law. Returned filter elements cannot be handled unless the appropriate return shipment form is completed.

12. Liability

Sartorius Stedim Biotech cannot assume liability if Virosart® HC capsules are subjected to improper use.

In the interest of product development we reserve the right to make changes.

Thank you for working with Sartorius Stedim Biotech, we appreciate your business.

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