



# Sartobind<sup>®</sup> Cassettes, 4 and 8 mm bed height

## Process-Scale Membrane Adsorbers

### Benefits

- Simple pod-like modular cassette system up to 100 L bed volume
- Flexible and modular system fits exactly to needed capacity
- Disposable or batch-reusable
- Directly scalable from a small capsule
- Designed for capture chromatography of virus vaccines, virus like particles (VLP), large proteins such as blood factors and conjugated proteins, glycoproteins, antibody drug conjugates (ADCs) as well as for flow-through polishing at largest scale
- Q anion exchanger also available in validated gamma-sterile version




### Product Information

Membrane chromatography was limited by size to 5 L membrane volume for almost a decade ignoring the needs for bind and elute chromatography. As capture of large proteins (such as blood factors, conjugated proteins or viruses and virus like particles) plays an increasing role, a new modular pod-like format is developed. The design is derived from the 4 | 8 mm Sartobind capsules to keep the flow path, bed heights, void volumes and last but not least the membrane types constant. A multiple number of cassettes can be quickly set up in a holder to get the exact volume needed. As a result the cassettes can be directly scaled from a smaller capsule size and pressure flow relation, shape of breakthrough curves and binding capacities per volume are independent of the number of cassettes applied.


The same Sartobind membranes are incorporated in the cassettes as in the capsules to be used for large scale capturing and impurity removal at high flow rates. Sartobind strong anion (Q) and cation (S) ligands are covalently attached on a flexible hydrogel onto the stabilized reinforced cellulose. For Sartobind STIC PA (primary amine) salt tolerant anion exchanger and Sartobind Phenyl hydrophobic interaction chromatography (HIC) adsorber, ligands are directly attached to the matrix. The membrane pore size of >3  $\mu\text{m}$  allows large proteins, bioparticles and viruses to enter the macroporous structure and achieves high binding capacity without size exclusion effects. The membrane is stacked to a bed height of 4 mm (0.8 L) or 8 mm (1.6 L) in a pod-like cassette setup.

If capture or polishing needs to be run under sterile conditions, Sartobind Q cassettes are the device of choice, as they can be delivered also gamma-irradiated and with validated sterility. Two applications that could especially benefit are virus purification processes that cannot be sterile filtered prior to filling due to the virus size, and continuous manufacturing operations that run over a prolonged period of time. Sterile Q cassettes fulfill all the same specifications as the non-sterile versions.


## Cassette Design



The cassette design is directly derived from the Sartobind capsules. Cutting and unrolling the membrane stacks results in two membrane stacks of 4 or 8 mm bed height. To achieve a small footprint, the stacks are unrolled. The feed enters between the membrane stack separated by a central spacer, flows through the membrane stacks and is collected at both sides of the outlet.





The side view of the cassette displays the upstream distribution channel where the feed enters between the two membrane stacks. The lower downstream channel collects the outlet of each cassette.




The direct comparison of capsule (upper cutaway) and cassette design (bottom) shows same construction principles, bed height and flow scheme.


## Breakthrough curve and pressure flow comparison to capsules



Sartobind devices 8 mm bed height from nano to Jumbo capsules are compared to 3 x and 10 x cassettes and display same breakthrough behavior. Sample: 2 g/L bovine serum albumin in Tris HCl, pH 7.2, 20 mM NaCl. Flow rate: 5 MV/min.



Break through curves of Sartobind devices with 4 mm and 8 mm bed height from nano to Jumbo capsules and cassette display same breakthrough curve. Sample: 2 g/L bovine serum albumin in Tris HCl, pH 7.2, 20 mM NaCl. Flow rate: 10 membrane volumes (MV) per minute in 4 mm and 5 MV/min in 8 mm devices



Pilot filter holder equipped with 13 Sartobind cassettes

### Capture

- Viruses (e.g. adeno- or lentivirus)
- Virus like particles (VLP)
- Large proteins such as blood factors
- Q can be delivered in a sterile cassette version if terminal sterile filtration is not possible

### Polishing

- Host cell proteins
- DNA and viruses
- Aggregates
- Endotoxins

## Technical Data

### Specifications

<b>Membrane materials</b>		
Matrix	Stabilized reinforced cellulose	
Membrane thickness   membrane volume = membrane area	275 µm   1 mL = 36.4 cm <sup>2</sup>	
Nominal pore size	> 3 µm	
Ion exchanger ligand Q	Strong anion Q: quaternary ammonium (R-CH <sub>2</sub> -N <sup>+</sup> (CH <sub>3</sub> ) <sub>3</sub> )	
Ion exchanger ligand STIC PA	Weak anion STIC PA: primary amine (R-NH <sub>2</sub> )	
Ion exchanger ligand S	Strong cation S: sulfonic acid (R-CH <sub>2</sub> -SO <sub>3</sub> <sup>-</sup> )	
Hydrophobic interaction ligand	HIC: Phenyl (R-NH-C <sub>6</sub> H <sub>5</sub> )	
<b>Cassette materials</b>		
Outer cage, seal, nonwoven, fleece	ABS, silicone, polyethylene	
<b>Operation</b>		
Depyrogenation	1 N NaOH for 30 minutes at 20°C	
Integrity testing	By the diffusion test method with Sartocheck® 4 Plus	
<b>Typical dynamic binding capacity at 10% breakthrough</b>		
Q (bovine serum albumin, 20 mM Tris/HCl, pH 7.5)	29 mg/mL (0.8 mg/cm <sup>2</sup> )	
STIC PA (BSA, 20 mM Tris HC, 150 mM NaCl, pH 7.5)	50 mg/mL (1.4 mg/cm <sup>2</sup> )	
STIC PA (salmon sperm DNA, 20 mM Tris/HCl, 150 mM NaCl pH 7.2)	10.9 mg/mL (0.3 mg/cm <sup>2</sup> )	
S (lysozyme, 10 mM potassium phosphate, pH 7.0)	25 mg/mL (0.7 mg/cm <sup>2</sup> )	
Phenyl (polyclonal bovine IgG, 50 mM potassium phosphat, 0.9 M (NH <sub>4</sub> ) <sub>2</sub> SO <sub>4</sub> , pH 7.5)	14.6 mg/mL (0.4 mg/cm <sup>2</sup> )	
<b>Ligand density</b>		
Q	2 – 5 µeq/cm <sup>2</sup>	
STIC PA	18 – 22 µeq/cm <sup>2</sup>	
S	2 – 5 µeq/cm <sup>2</sup>	
Phenyl	3 µeq/cm <sup>2</sup>	
<b>pH stabilities</b>		
	Short term	Long term
Q	pH 1 – 14	pH 2 – 12
STIC PA	pH 2 – 14	not defined
S	pH 3 – 14	pH 4 – 13
Phenyl	pH 2 – 14	pH 3 – 13
<b>Gamma irradiated products</b>		
	Gamma dose	Sterile validated
Sterile Sartobind Q cassettes 4 mm and 8 mm	≥ 25 kiloGrays (kGy) maximum dose: 50 kGy	Yes, minimum dose for sterility is 7 – 8 kGy
Gamma irradiated manifold	See above	No

## Dimensions and Connections

Membrane volume   4 mm	1 mL	75 mL	200 mL	400 mL	600 mL	2,500 mL	800 mL
Membrane volume   8 mm	3 mL	150 mL	400 mL	800 mL	1,200 mL	5,000 mL	1,600 mL
Size	Nano	5"	10"	20"	30"	Jumbo	Cassette



Membrane area cm <sup>2</sup>   4 mm	36.4	2,700	7,300	14,600	22,000	91,000	29,000*
Membrane area cm <sup>2</sup>   8 mm	110	5,500	14,600	29,000	44,000	182,000	58,000*
Void volume mL MV**   4 mm	3.5   3.5	2.7   200	2.7   540	2.7   1,080	2.7   1,600	2.8   7,000	3.1   2500*
Void volume mL MV**   8 mm	1.3   4	1.3   200	1.4   540	1.4   1,080	1.3   1,600	1.4   7,000	1.8   2900*
Dimensions one device (height × diameter) mm	37 × 33 h × Ø	190 × 77 h × Ø	350 × 100 h × Ø	570 × 100 h × Ø	810 × 100 h × Ø	850 × 302 h × Ø	634 × 387 × 49 w × h × t
Connectors	Luer female	Sanitary 3/4" 25 mm outer, 14 mm inner diameter	Sanitary 1 1/2" 50.5 mm outer, 36 mm inner diameter	Sanitary 1 1/2" 50.5 mm outer, 36 mm inner diameter	Sanitary 1 1/2" 50.5 mm outer, 36 mm inner diameter	Sanitary 1 1/2" 50.5 mm outer, 36 mm inner diameter	Sanitary 1 1/2" 50.5 mm outer, 36 mm inner diameter via manifolds (accessory)
Hose barb version (height × diameter) mm	n.a.	203 × 77	n.a.	n.a.	n.a.	n.a.	n.a.
Approximate weight	10 g	400 g	760 g	1.3 kg	1.9 kg	16 kg 20 kg wet 23 kg filled	4.9 kg* 6.0 kg* wet

n.a. = not available, \*Multiply with the number of used cassettes \*\*MV = membrane volume (including the porosity of the membrane which is 80%)

## Size | Membrane volume options for cassettes in the Pilot Filter Holder

Number of cassettes	1	2	3	4	5	6	7	8	9	10	11	12	13
Membrane volume L   4 mm	0.8	1.6	2.4	3.2	4	4.8	5.6	6.4	7.2	8	8.8	9.6	10.4
Membrane volume L   8 mm	1.6	3.2	4.8	6.4	8	9.6	11.2	12.8	14.4	16	17.6	19.2	20.8
Membrane area 4 mm. m <sup>2</sup>	2.9	5.8	8.7	11.6	14.5	17.4	20.3	23.2	26.1	29.0	31.9	34.8	37.7
Membrane area 8 mm. m <sup>2</sup>	5.8	11.6	17.4	23.2	29.0	34.8	40.6	46.4	52.2	58.0	63.8	69.6	75.4
Void volume 4 mm. L	2.5	5.0	7.5	10.0	12.5	15.0	17.5	20.0	22.5	25.0	27.5	30.0	32.5
Void volume 8 mm. L	2.9	5.8	8.7	11.6	14.5	17.4	20.3	23.2	26.1	29.0	31.9	34.8	37.7
Approximate weight kg	4.9	9.8	14.7	19.6	24.5	29.4	34.3	39.2	44.1	49	53.9	58.8	63.7
Approximate weight wet kg	6	12	18	24	30	36	42	48	54	60	66	72	78
Weight of manifold set. 1 in   1 out kg	11.2	11.2	11.2	11.2	11.2	11.2	11.2	11.2	11.2	11.2	11.2	11.2	11.2
Weight Pilot Filter Holder kg	160	160	160	160	160	160	160	160	160	160	160	160	160
Total weight kg	176	181	186	191	196	201	206	210	215	220	225	230	235
Total weight wet kg	182	193	204	215	226	237	248	258	269	280	291	302	313

Size of the Pilot Filter Holder for all configurations: 1327 × 955 × 790 mm (w × h × d)





## Ordering Information






### Sartobind Cassettes

Order number	Description	Quantity	Bed height [mm]	Protein binding capacity [g]	Recommended flow rate [L/min]	Maximum pressure [MPa] (bar   psig)
<b>Sartobind Q</b>						
98IEXQ42D-L	Sartobind Q Cassette 0.8 L, 4 mm, 1½" sanitary clamp via manifold set (accessory), manual, certificate	1	4	23.2	16	0.2 (2   29)
98IEXQ42E-P	Sartobind Q Cassette 1.6 L, 8 mm, 1½" sanitary clamp via manifold set (accessory), manual, certificate	1	8	46.4	8	0.2 (2   29)
98IEXQ42DGL	Sartobind Q Cassette 0.8 L, 4 mm, sterile by gamma irradiation at ≥ 25 kGy, max.50 kGy, 1½" sanitary clamp via gamma manifold set (accessory), manual, certificate	1	4	23.2	16	0.2 (2   29)
98IEXQ42EGP	Sartobind Q Cassette 1.6 L, 8 mm, sterile by gamma irradiation at ≥ 25 kGy, max.50 kGy, 1½" sanitary clamp via gamma manifold set (accessory), manual, certificate	1	8	46.4	8	0.2 (2   29)

Order number	Description	Quantity	Bed height [mm]	Protein binding capacity [g]	Recommended flow rate [L/min]	Maximum pressure [MPa] (bar   psig)
<b>Sartobind S</b>						
98IEXS42D-L	Sartobind S Cassette 0.8 L, 4 mm, 1½" sanitary clamp via manifold set (accessory), manual, certificate	1	4	20.3	16	0.2 (2   29)
98IEXS42E-P	Sartobind S Cassette 1.6 L, 8 mm, 1½" sanitary clamp via manifold set (accessory), manual, certificate	1	8	40.6	8	0.2 (2   29)
<b>Sartobind STIC PA</b>						
98STPA42D-L	Sartobind STIC PA Cassette 0.8 L, 1½" sanitary clamp via manifold set (accessory), manual, certificate	1	4	40.6	16	0.2 (2   29)
<b>Sartobind Phenyl</b>						
98HICP42E-P	Sartobind Phenyl Cassette 1.6 L, 8 mm, 1½" sanitary clamp via manifold set (accessory), manual, certificate	1	8	23.2	8	0.2 (2   29)

## Accessories

Order number	Description	Quantity
29Z-S00001	 Manifold set, 1½" sanitary clamp, inlet and outlet adapter plate	2
29Z-S00003	 Gamma irradiated manifold set, 1½" sanitary clamp, inlet and outlet adapter plate, ≥ 25 kGy, max.50 kGy	2
2ZGL--0005	 Pilot Filter Holder	1
2ZGL--0006	 Process filter holder	1

Order number		Description	Quantity
2ZGL--0007		Double process filter	1
2ZGL--0008		Drip pan for Pilot Filter holder	1
2ZGL--0015		Drip pan for Process and double Process Filter Holder	1
2ZGL--0013		Separation Plate for serial or parallel fluid distribution	1
2ZGL--0014		Pressure Safety Device, to prevent excessive pressure during compressed air draining of the cassettes	1
16290		Sartocheck® 3 Plus Integrity Tester	1
26288---FT		Sartocheck® 4 Plus Filter Integrity Tester	1

Sartorius Stedim Biotech GmbH  
August-Spindler-Strasse 11  
37079 Goettingen, Germany  
Phone +49.551.308.0  
www.sartorius-stedim.com

USA Toll-Free +1.800.368.7178  
Argentina +54.11.4721.0505  
Brazil +55.11.4362.8900  
Mexico +52.55.5562.1102  
UK +44.1372.737159  
France +33.442.845600  
Italy +39.055.63.40.41  
Spain +34.913.586.098  
Russian Federation +7.812.327.53.27  
Japan +81.3.4331.4300  
China +86.21.6878.2300

Specifications subject to change without notice.  
Copyright Sartorius Stedim Biotech GmbH.  
Printed in the EU on paper bleached without chlorine.  
Publication No.: SL-2089-e180402  
Order No.: 85037-557-43  
Ver. 04 | 2018

Chromatography as easy as filtration



◀ Sartobind App  
[www.sartorius-stedim.com/apps](http://www.sartorius-stedim.com/apps)