Virosart® Max MidiCaps®

Directions for Use

1. Summary

The purpose of this document is to provide the appropriate method for identifying, installing and using Sartorius Stedim Biotech Virosart® Max MidiCaps®. The procedure should be strictly followed to ensure operator safety and to maintain optimum performance of Virosart® Max MidiCaps®.

This document pertains to Sartorius Stedim Biotech Virosart® Max filter elements with the following part number:

Virosart® Max MidiCaps® (box of 2):
54A5358N9--FF--V

All Virosart® Max MidiCaps® have sanitary inlet & outlet connections.

2. Labeling of Virosart® Max MidiCaps®

The type, pore size, lot number and serial number are imprinted on the outer support of the housing of the Virosart® Max MidiCaps®. This information is identical to the information found on the label of both the plastic bag and the cardboard box where in addition to type and lot number the order number is indicated. The direction of filtration is indicated by an arrow on the MidiCap® housing.

3. Installing of Virosart® Max MidiCaps®

The plastic bag of the Virosart® Max MidiCaps® should be cut open at one end. Leave the open plastic bag on the filter for protection as long as possible.

⚠ In a validated production step Virosart® Max MidiCaps® are heat sterilized in gas permeable protection plastic bags according to DIN | EN 554. Sterility, however, is only guaranteed when this primary packaging is intact before using the MidiCaps® as a sterile product. If the packaging is damaged in any way the MidiCaps® must be sterilized by the end user according to the procedure described under S5.1 or 5.3 before use.

If the MidiCap® is to be used as a sterile product care must be taken to prevent contamination when connecting the capsule to the filtration system. Virosart® Max MidiCaps® are connected to the processing line at the inlet and outlet via triclamp (please do not forget the gaskets!) connections.

4. Wetting of Virosart® Max MidiCaps®

(In the direction of filtration)

Before every sterilization (autoclaving) and prior to every integrity test the entire system needs to be rinsed to ensure that it is particle-free. To rinse the filter element, open the vent valve on the MidiCap® housing and fill the housing with WFI until liquid comes out of the vent valve. To ensure that all the air has been removed from the housing, keep the vent valve slightly opened for the first minute of flushing. Then close the vent valve and rinse for approx. 3-4 min at a differential pressure of 2.0 bar | 29 psi or alternative 4-5 min at a differential pressure of 1.0 bar | 14.5 psi.

The housing of the Virosart® Max MidiCaps® can be drained by opening the drain valve and using compressed air or nitrogen with a maximal pressure of 0.5 bar | 7 psi in the direction of filtration. It can also be drained by simply opening drain and vent valves, allowing the solution to drain out. After the housing is empty, close all the valves. After the Virosart® Max has been flushed according to the instructions, it is thoroughly wetted.

5. Sterilization

5.1 Autoclaving of Virosart® Max MidiCaps®

Sterilization temperature: up to 121°C, 30 min

The inlet and outlet as well as any open valves of the Virosart® Max MidiCaps® housing should be wrapped with steam permeable autoclave paper or aluminium foil. Once the sterilization temperature has been attained the Virosart® Max MidiCaps® can be autoclaved for 30 minutes. During the autoclave cycle of the Virosart® Max MidiCaps®, make sure that connecting pieces on the Virosart® Max MidiCaps® (e.g. made of stainless steel) are not attached too tightly and do not have a torque effect on the Virosart® Max MidiCaps® housing which could cause deformation.

5.2 In-line Steam Sterilization

No in-line steam sterilization of Virosart® Max MidiCaps®!

6. Integrity Testing of Virosart® Max MidiCaps®

It is important to follow the wetting procedure for Virosart® Max to make sure that the filter is wetted out completely prior to the integrity test. The integrity test of a Virosart® Max MidiCaps® must be done using an automated integrity test unit, preferably the Sartocheck® integrity test unit. When using a Sartorius Stedim Biotech Sartocheck® integrity test system, the diffusion test can be performed automatically, recording the corresponding integrity test data. For additional information on how to perform an automatic integrity test please see our installation and operating Instructions for the Sartocheck® integrity test unit.

Integrity Test Data for Virosart® Max MidiCaps®

Integrity test	Filter type	Test pressure	Limit
Water based diffusion test	54A5358N9	2 bar 29 psi	3 ml/min

Remark:

All values given above were obtained at room temperature and using water as wetting agent.

7. Maximum Allowable Differential Pressure of Virosart® Max MidiCaps®

In the direction of filtration
At 20°C max. 5.0 bar | 73 psi
At 121°C max. 0.5 bar | 7.3 psi

In the reversed direction of filtration
At 20°C max. 0.5 bar | 7.3 psi

8. Filtration with Virosart® Max MidiCaps®

Virosart® Max is especially designed as a pre-filter for the final virus filtration step. Virus filtration in general should be realized at the end of the purification process prior to the last buffer exchange | form and fill. The protein solution which is subject to the virus filtration must be free of aggregates, host cell protein & DNA to avoid premature blockage of the final virus filter. Even if these contaminants should be removed either during the initial recovery process or latest during the polishing process of the target molecule, small amounts are sufficient to cause premature blockage of the final virus filter. To prevent this, an efficient pre-filtration step, like the Virosart® Max, is required as protection for the final virus filter. The optimum pre-filter – final filter ratio has to be identified during development of this process step. To make full use of the filtration area, the housing must be vented completely free of air. To do this, briefly open the vent valve until the liquid to be filtered comes out. Then apply the desired filtration pressure. The pre-filter should be used in-line with the final virus filter under the same conditions regarding to pressure, in general 2 bar | 29 psi inlet pressure. The flow rate achieved through the Virosart® Max, in-line with the final virus filter, is determined by the flow rate of the final virus filter.

9. Changing out the Virosart® Max MidiCaps®

The Virosart® Max MidiCaps® must be changed, at the latest, when the maximum allowable differential pressure is reached or when the validated process volume has been processed over the filter.

10. Applications Support

If required, please contact your local Sartorius Stedim Biotech representative to obtain further information on technical data, product integrity testing or general information concerning specific applications. Our application specialists from Purification Technologies will gladly assist in your filtration trial work if needed. Moreover our Validation Service and support team are always available to assist in questions concerning product and process validation, as well as training of staff and optimization of filtration processes.

11. Return of used filter elements

If you wish to return used filter elements to Sartorius Stedim Biotech for inspection, it is necessary to decontaminate and sterilize them prior to shipment. This procedure has to be documented in the Return Shipment Form, available from your local Sartorius Stedim Biotech representative. This documentation is required by country law. Returned filter elements cannot be handled unless the appropriate Return Shipment Form is completed.

12. Liability

Sartorius Stedim Biotech cannot assume liability if Virosart® Max MidiCaps® are subjected to improper use. In the interest of product development we reserve the right to make changes.

Thank you for working with Sartorius Stedim Biotech, we appreciate your business.

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Sartorius Stedim Biotech GmbH