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Contents | Inhalt

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
Table of Contents

Sterisart® NF gamma – Overview	4
Product description	5
Packaging	5
Intended application	6
Specification	6
1. Preparation for use in the sterilization area of an isolator	7
2. Before sample filtration	8
3. Sampling	9
4. Taking of Washing disolution solution	10
5. After filtration of the samples and the washing solution	11
6. Fill containers with liquid culture media	11
7. Incubation of the containers according to the respective regulation	12

Sterisart® NF gamma – Overview

Instructions for Use
in Sterility Testing According to
International Pharmacopeias

Pos.	Designation
1	Pre-installed tube clamp
2	Vent filter
3	Sterisart® container
4	Wing nut plug
5	Rubber cap
6	Adapter with cap for medical containers Bags with male luer lock connector
7	Double-needle metal spike with vent filter for removing diluent from closed containers
8	Tubing
9	Additional plastic bag

Scope of the instructions for use

These instructions are applicable to the following products:

16478 GBD

Sterisart® NF gamma, system for medical containers with luer lock male fittings (pack of 10, double-packed, gamma irradiation sterilized)

Product description
(see page 4)

The Sterisart® NF systems listed are ready-to-use assemblies, each consisting of two clear plastic containers that are connected by flexible tubing to an adapter on the Y-distributor for transfer of equal sample volumes into the containers.

Each container has a 50-mm diameter, 0.45 µm pore size validated membrane filter sealed at its base by a clamping technique that prevents diffusion of liquid into the filter-container boundary area. In addition, a validated sterilizing grade vent filter is attached to the top of each container. Graduated marks on each container indicate approximately 50, 75 and 100-ml levels.

The combination adapter for sampling consists of two Y-distributors connected to each other by a clued adapter piece. One Y-distributor is connected to the Sterisart® containers by a clued double hose. The other Y-distributor is connected with two single hoses of different length, each fitted with a white clamp. The shorter hose is connected to a luer lock femal connector which enables samples to be removed from medical containers|bags with male luer lock fittings. The second longer hose is fitted with a dual needle metal spike with a 2 cm needle for removal of washing-|rinsing solution or media from spikable, closed containers.

Color-coded clamps (two white, two yellow) are pre-installed on the tubing for convenient closing as required.

All accessories needed (two rubber caps for the vent filters and two plastic wing nut plugs for the container outlets) are included with each system.

Packaging

Sterisart® NF gamma:

A single box consists of a Sterisart® system packed with appropriate accessories in a gas-impermeable, clear plastic tray which is sealed on top by a gas-impermeable cover sheet.

A label on the cover gives details of the enclosed system (brief description, order number, expiration date, sterilization date and lot number). Sterilization is carried out with gamma-irradiation. The gas-impermeable boxes enable the outer surfaces to be sterilized by usual isolator procedures without passage of sterilant (e.g., H₂O₂) through to the box and the Sterisart® NF system.

A change from orange to red of the sterilization indicator located on the outer standard box as well as on the accessory package inside the individual system shows that the product has been reliably sterilized by gamma-irradiation.

A standard pack of Sterisart® gamma contains ten sterile single packed systems which are altogether additionally packed in a plastic bag. Each gamma system has a colour sterilization indicator.

A turquoise dot on the end of the box additionally identifies type 16478.

Intended application

Sterisart® NF gamma systems have been specifically designed for carrying out the sterility test by the membrane filter method described in various pharmacopeias, in combination with a dedicated peristaltic pump. The 16478 version permits sterility testing of liquids from medical containers that have male luer lock fittings. Sterisart® NF gamma is the preferred type for sterility testing in isolators.

Important notes on safe and correct usage

Sterisart® NF systems are only to be used with peristaltic pumps designed for sterility testing with such systems, such as the Sterisart® Universal pumps 16429 and 16420.

All necessary precautions must be taken to avoid secondary contamination. Handle the parts only when wearing sterile gloves.

The dual needle metal spike is very sharp. Take great caution to avoid piercing gloves or injuring hands when handling them.

Check the color indicator and the expiration date on the box before using a system.

To avoid media leakage during incubation, securely attach the plastic wing nut plugs to the outlets of the Sterisart® NF containers by turning them clockwise.

Never "glue" the rubber caps onto the Sterisart® NF vent filters! The caps are intended to serve not only as closures, but also as safety valves. Should exceptionally high pressure build up in the containers during pumping (high pump speed, clogged filter), normally fitted caps will blow off before the burst pressure of other system parts is reached.

Specifications

Max. operating pressure: 3 bar (44 psi) at 20°C


Max. operating temperature: 50°C

Chemical compatibility: Details available on request


Please follow the instructions for use starting on the next page.

1. Preparation for use in the sterilization area of an isolator


Open the standard package and the plastic bag to remove the quantity of individual packs required for the number of sterility tests to be performed.

Place the individual packs, still enclosed in the gas-impermeable box, in the sterilization area.


Sterilize the outer surface of the box according to the isolator validation specifications, for example, with H₂O₂.


Open the box or, if appropriate, pass the closed Sterisart® NF boxes through the air-lock port on the isolator.


2. Before sample filtration

Cut open the box and take out the Sterisart® NF system. Position the two Sterisart® NF containers in the pump holder.


Thread the Sterisart® NF tubing in the pump head, following the instructions given for the pump used.


Tightly seal the two vent filters using the enclosed rubber caps.

3. Sampling

Sampling from closed collapsible containers with a **male** luer lock fitting.

Turn the plug counterclockwise to detach it from the female luer lock connector.


Turn the female luer lock connector on the Sterisart® NF system counterclockwise to connect it to the sample plastic bag with a male luer lock fitting.


Switch on the pump to transfer equal volumes of the sample into the two Sterisart® NF containers; at the same time, the plastic bag deflates as vacuum builds up inside. The white clamp fitted on the hose with dual needle metal spike is open.

Switch off the pump. Follow the same procedure with further containers of sample.

4. Taking of Washing|disolution solution


Turn the female luer lock connector on the Sterisart® NF system counterclockwise to remove from the plastic back.

Alternatively the Sterisart® System can stay connected with the plastic bag. In this case the clamp on the hose connected with the plastic bag has to be closed during the filling of washing-|disolution solution with the dual needle metal spike.

The washing solution taking with the luer lock adapter from containers with male luer lock connectors is also possible. In this case the dual needle metal spike is not necessary (clamp on this hose is closed).


To transfer washing or disolution solution from closed , unvented containers with the dual needle metal spike, the clamp on the hose connected to the dual needle metal spike is closed and the dual needle metal spike is inserted sufficiently far into the port of the closed, unvented container.


Switch on the pump to transfer equal volumes into the two Sterisart® NF containers; at the same time, sterile air is drawn in through the spike.

Turn of the pump and draw out the spike. Follow the same procedure with washing solution containers and to wash out inhibitors.

5. After filtration of the samples and the washing solution


Switch off the pump. Tightly seal the outlet of each Sterisart® NF container with a plastic wing nut plug, twisting it about half a turn clockwise.


Remove the rubber caps from the vent filters.

6. Fill containers with liquid culture media


Culture medium 1 (thioglycollate broth):
Close the white tube clamp at the outlet of the Y-distributor.


Insert the dual needle metal spike into the pierceable closure of a closed container of the broth.

Switch on the pump. Switch off the pump when a sufficient volume has been transferred to the Sterisart® container. Take out the dual needle metal spike.


Culture medium 2 (tryptic soybean broth):
Open the white tube clamp, and close the adjacent yellow clamp. Transfer the medium as described above for thioglycollate broth.

Seal off the two tubes just above the Sterisart® NF containers by closing the pre-installed tube clamps.
Cut off the tubes above the clamps.

7. Incubation of the containers according to the respective regulations

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
Copyright by
Sartorius Stedim Biotech GmbH,
Goettingen, Germany.
All rights reserved. No part of
this publication may be reprinted
or translated in any form or
by any means without the prior
written permission of Sartorius
Stedim Biotech GmbH.
The status of the information,
specifications and illustrations
in this manual is indicated by
the date given below.
Sartorius Stedim Biotech GmbH
reserves the right to make changes
to the technology, features,
specifications and design of the
equipment without notice.

Status:
January 2010,
Sartorius Stedim Biotech GmbH,
Goettingen, Germany