



## Sartobind<sup>®</sup> Q and S, 4 and 8 mm bed height

### Process-Scale Strong Ion Exchange Membrane Adsorbers

#### Benefits

- Designed for high speed flow-through polishing as well as bind and elute chromatography
- Recommended flow rate of 5 to 30 membrane volumes per minute
- Easy and simple handling (like a filter)
- Ten times higher dynamic binding capacity for adeno-, lenti- or other virus like particles (VLPs) compared to resins
- Reusable or disposable
- Q anion exchanger also available in validated gamma-sterile version



#### Product Information

The Sartobind Q and S capsules and cassettes with 4 | 8 mm bed height are membrane chromatography devices for large scale capturing and impurity removal at high flow rates. Sartobind strong anion (Q) and cation (S) ligands are covalently attached as a flexible hydrogel onto the stabilized reinforced cellulose. The membrane pore size of  $>3 \mu\text{m}$  allows large proteins, bioparticles and viruses to enter the macroporous structure and achieve high binding capacity without size exclusion effects. In the capsules the membrane is rolled up to form a cylinder with a bed height of 4 | 8 mm around a central core. The cassettes contain two membrane stacks of 4 or 8 mm bed height with a central spacer inbetween.

The optimized design of the fluid channels reduces the void volume significantly, resulting in sharp breakthrough curves with minimal elution volumes (see application note: order no.: 85034-536-64). The recommended flow rate for devices with 4 mm bed height is 10 to 30 membrane volumes (MV) per minute and for 8 mm devices 5 MV/min. The internal support structures and the outer shell of the capsules are made from polypropylene and from acrylonitrile butadiene styrene (ABS) for the cassettes. The membrane geometry and the flow path is identical in all devices. Sartobind nano 1 or 3 mL is the recommended scale down model for the larger 4 or 8 mm devices. If sample amount is limited, Sartobind Q or S pico units with 0.08 mL bed volume (4 mm) are available. For buffer and conditions screening, additionally 96 well plates with a bed height of 0.8 mm can be used.

If capture or polishing needs to be run under sterile conditions, Sartobind Q cassettes are the devices of choice, as they can be delivered also gamma-irradiated and with validated sterility. Two applications that could especially benefit are virus purification processes that cannot be sterile filtered prior to filling due to the virus size, and continuous manufacturing operations that run over a prolonged period of time. Sterile Q cassettes fulfill all the same specifications as the non-sterile versions.

### Applications

Target molecules to be purified: monoclonal antibodies, recombinant proteins, conjugated proteins, glycoproteins, antibody drug conjugates (ADCs), virus particles, virus like particles (VLPs)

### Polishing


- Virus clearance >6 log
- DNA below detection limit
- Endotoxins >5 log
- Host cell proteins >99%
- Leached ligand materials
- Aggregates

### Capture

- Large proteins
- Viruses and virus like particles (VLP)
- Vaccines
- Proteins from large feed streams
- Q can be delivered in a sterile cassette version if terminal sterile filtration is not possible



96 well plates and pico devices for screening of conditions and initial tests



Cutaway view and flow path of a capsule (upper figures) and a cassette (lower figures).

## Technical Data

<b>Membrane materials</b>		
Matrix	Stabilized reinforced cellulose	
Membrane thickness   membrane volume = membrane area	275 µm   1 mL = 36.4 cm <sup>2</sup>	
Nominal pore size	> 3 µm	
Ion exchanger ligand Q	Strong anion Q: quaternary ammonium (R-CH <sub>2</sub> -N <sup>+</sup> (CH <sub>3</sub> ) <sub>3</sub> )	
Ion exchanger ligand S	Strong cation S: sulfonic acid (R-CH <sub>2</sub> -SO <sub>3</sub> <sup>-</sup> )	
<b>Capsule materials</b>		
Outer cage, inner core, end caps, capsule housing, nonwoven, fleece	Polypropylene	
O-ring in vent valve (except nano)	EPDM (ethylene propylene diene monomer)	
<b>Cassette materials</b>		
Outer cage, seal, nonwoven, fleece	ABS, silicone, polyethylene, stable to gamma irradiation	
<b>Operation</b>		
Depyrogenation	1 N NaOH for 30 minutes at 20°C	
Integrity testing	By the diffusion test method with Sartocheck® 4 Plus	
<b>Typical dynamic binding capacity at 10% breakthrough</b>		
Q (bovine serum albumin, 20 mM Tris/HCl, pH 7.5)	29 mg/mL (0.8 mg/cm <sup>2</sup> )	
S (lysozyme, 10 mM potassium phosphate, pH 7.0)	25 mg/mL (0.7 mg/cm <sup>2</sup> )	
<b>Ligand density</b>		
Q	2 – 5 µeq/cm <sup>2</sup>	
STIC PA	18 – 22 µeq/cm <sup>2</sup>	
S	2 – 5 µeq/cm <sup>2</sup>	
Phenyl	3 µeq/cm <sup>2</sup>	
<b>pH stabilities</b>		
Q	Short term pH 1 – 14	Long term pH 2 – 12
S	Short term pH 3 – 14	Long term pH 4 – 13
<b>Gamma irradiated products</b>		
Sterile Sartobind Q cassettes 4 mm and 8 mm	Gamma dose ≥ 25 kiloGrays (kGy) maximum dose: 50 kGy	Sterile validated Yes, minimum dose for sterility is 7 – 8 kGy
Gamma irradiated manifold	See above	No

## Dimensions and Connections

Membrane volume   4 mm	1 mL	75 mL	200 mL	400 mL	600 mL	2,500 mL	800 mL
Membrane volume   8 mm	3 mL	150 mL	400 mL	800 mL	1,200 mL	5,000 mL	1,600 mL
Size	Nano	5"	10"	20"	30"	Jumbo	Cassette



Membrane area cm <sup>2</sup>   4 mm	36.4	2,700	7,300	14,600	22,000	91,000	29,000*
Membrane area cm <sup>2</sup>   8 mm	110	5,500	14,600	29,000	44,000	182,000	58,000*
Void volume mL MV**   4 mm	3.5   3.5	2.7   200	2.7   540	2.7   1,080	2.7   1,600	2.8   7,000	3.1   2500*
Void volume mL MV**   8 mm	1.3   4	1.3   200	1.4   540	1.4   1,080	1.3   1,600	1.4   7,000	1.8   2900*
Dimensions one device (height × diameter) mm	37 × 33 h × Ø	190 × 77 h × Ø	350 × 100 h × Ø	570 × 100 h × Ø	810 × 100 h × Ø	850 × 302 h × Ø	634 × 387 × 49 w × h × t
Connectors	Luer female	Sanitary 3/4" 25 mm outer, 14 mm inner diameter	Sanitary 1 1/2" 50.5 mm outer, 36 mm inner diameter	Sanitary 1 1/2" 50.5 mm outer, 36 mm inner diameter	Sanitary 1 1/2" 50.5 mm outer, 36 mm inner diameter	Sanitary 1 1/2" 50.5 mm outer, 36 mm inner diameter	Sanitary 1 1/2" 50.5 mm outer, 36 mm inner diameter via manifolds (accessory)
Hose barb version (height × diameter) mm	n.a.	203 × 77	n.a.	n.a.	n.a.	n.a.	n.a.
Approximate weight	10 g	400 g	760 g	1.3 kg	1.9 kg	16 kg 20 kg wet 23 kg filled	4.9 kg* 6.0 kg* wet

n.a. = not available, \*Multiply with the number of used cassettes \*\*MV = membrane volume (including the porosity of the membrane which is 80%)

## Ordering Information

### Sartobind Q Capsules

Order number	Description	Quantity	Bed height [mm]	Protein binding capacity [g]	Recommended flow rate [L/min]	Maximum pressure [MPa] (bar/psig)
96IEXQ42DN-11	Sartobind Q nano 1 mL, 4 mm, Luer female connectors, 2 Luer male PEEK adapters to UNF 10-32 female, manual, certificate	1	4	0.029	0,02	0.4 (4/58)
96IEXQ42DN-11--A	Sartobind Q nano 1 mL, 4 mm, Luer female connectors, 2 Luer male PEEK adapters to UNF 10-32 female, manual, certificate	4	4	0.029	0,02	0.4 (4/58)
96IEXQ42D9M00--A	Sartobind Q 75 mL, 4 mm, hose barb connectors, manual, certificate	4	4	2.16	1.5	0.4 (4/58)
96IEXQ42D9MFF--A	Sartobind Q 75 mL, 4 mm, ¾" sanitary clamp, manual, certificate	4	4	2.16	1.5	0.4 (4/58)
96IEXQ42D1GSS	Sartobind Q 200 mL, 4 mm, 1½" sanitary clamp, manual, certificate	1	4	5.8	4	0.4 (4/58)
96IEXQ42D2HSS	Sartobind Q 400 mL, 4 mm, 1½" sanitary clamp, manual, certificate	1	4	11.7	8	0.4 (4/58)
96IEXQ42D3KSS	Sartobind Q 600 mL, 4 mm, 1½" sanitary clamp, manual, certificate	1	4	17.6	12	0.4 (4/58)
96IEXQ42D3NSS	Sartobind Q Jumbo 2.5 L, 4 mm, 1½" sanitary clamp, 2 protective caps, manual, autoclaving instructions, certificate	1	4	73	50	0.3 (3/43.5)
96IEXQ42EUC11--A	Sartobind Q nano 3 mL, 8 mm, Luer female connectors, 2 Luer male PEEK adapters to UNF 10-32 female, manual, certificate	4	8	0.088	0.015	0.4 (4/58)
96IEXQ42E2LSS	Sartobind Q 800 mL, 8 mm, 1½" sanitary clamp, manual, certificate	1	8	23.3	4	0.4 (4/58)
96IEXQ42E3FSS	Sartobind Q 1.2 L, 8 mm, 1½" sanitary clamp, manual, certificate	1	8	35	6	0.4 (4/58)
96IEXQ42E3ESS	Sartobind Q Jumbo 5 L, 8 mm, 1½" sanitary clamp, 2 protective caps, 15 membrane discs 30 mm, manual, autoclaving instructions, certificate	1	8	145	25	0.3 (3/43.5)

Order number	Description	Quantity	Bed height [mm]	Protein binding capacity [g]	Recommended flow rate [L/min]	Maximum pressure [MPa] (bar/psig)
98IEXQ42D-L	Sartobind Q Cassette 0.8 L, 4 mm, 1½" sanitary clamp via manifold set (accessory), manual, certificate	1	4	23.2	16	0.2 (2/29)
98IEXQ42E-P	Sartobind Q Cassette 1.6 L, 8 mm, 1½" sanitary clamp via manifold set (accessory), manual, certificate	1	8	46.4	8	0.2 (2/29)
98IEXQ42DGL	Sartobind Q Cassette 0.8 L, 4 mm, sterile by gamma irradiation at ≥ 25 kGy, max. 50 kGy, 1½" sanitary clamp via gamma manifold set (accessory), manual, certificate	1	4	23.2	16	0.2 (2   29)
98IEXQ42EGP	Sartobind Q Cassette 1.6 L, 8 mm, sterile by gamma irradiation at ≥ 25 kGy, max. 50 kGy, 1½" sanitary clamp via gamma manifold set (accessory), manual, certificate	1	8	46.4	8	0.2 (2   29)

#### Sartobind S Capsules

Order number	Description	Quantity	Bed height [mm]	Protein binding capacity [g]	Recommended flow rate [L/min]	Maximum pressure [MPa] (bar/psig)
96IEXS42DN-11	Sartobind S nano 1 mL, 4 mm, Luer female connectors, 2 Luer male PEEK adapters to UNF 10-32 female, manual, certificate	1	4	0.025	0.02	0.4 (4/58)
96IEXS42DN-11--A	Sartobind S nano 1 mL, 4 mm, Luer female connectors, 2 Luer male PEEK adapters to UNF 10-32 female, manual, certificate	4	4	0.025	0.02	0.4 (4/58)
96IEXS42D9M00--A	Sartobind S 75 mL, 4 mm, hose barb connectors, manual, certificate	4	4	1.89	1.5	0.4 (4/58)
96IEXS42D9MFF--A	Sartobind S 75 mL, 4 mm, ¾" sanitary clamp, manual, certificate	4	4	1.89	1.5	0.4 (4/58)
96IEXS42D1GSS	Sartobind S 200 mL, 4 mm, 1½" sanitary clamp, manual, certificate	1	4	5.1	4	0.4 (4/58)
96IEXS42D2HSS	Sartobind S 400 mL, 4 mm, 1½" sanitary clamp, manual, certificate	1	4	10.2	8	0.4 (4/58)

Order number	Description	Quantity	Bed height [mm]	Protein binding capacity [g]	Recommended flow rate [L/min]	Maximum pressure [MPa] (bar/psig)
96IEXS42D3KSS	Sartobind S 600 mL, 4 mm, 1½" sanitary clamp, manual, certificate	1	4	15.4	12	0.4 (4/58)
96IEXS42EUC11--A	Sartobind S nano 3 mL, 8 mm, Luer female connectors, 2 Luer male PEEK adapters to UNF 10-32 female, manual, certificate	4	8	0.077	0.015	0.4 (4/58)
96IEXS42E9BFF	Sartobind S 150 mL, 8 mm ¾" sanitary clamp, manual, certificate	1	8	3.9	0.75	0.4 (4/58)
96IEXS42E1HSS	Sartobind S 400 mL, 8 mm, 1½" sanitary clamp, manual, certificate	1	8	10.2	2	0.4 (4/58)
96IEXS42E2LSS	Sartobind S 800 mL, 8 mm, 1½" sanitary clamp, manual, certificate	1	8	20	4	0.4 (4/58)
96IEXS42E3FSS	Sartobind S 1.2 L, 8 mm, 1½" sanitary clamp, manual, certificate	1	8	31	6	0.4 (4/58)
96IEXS42E3ESS	Sartobind S 5 L, 8 mm, 1½" sanitary clamp, manual, certificate	1	8	127	25	0.3 (3/43.5)
98IEXS42D-L	Sartobind S Cassette 0.8 L, 4 mm, 1½" sanitary clamp via manifold set (accessory), manual, certificate	1	4	20.3	16	0.2 (2/29)
98IEXS42E-P	Sartobind S Cassette 1.6 L, 8 mm, 1½" sanitary clamp via manifold set (accessory), manual, certificate	1	8	40.6	8	0.2 (2/29)

## Related products

### Sartobind Q






Order number	Description	Quantity
99IEXQ42GC-----V	Sartobind Q 96-well plate	2
99IEXQ42GC-----D	Sartobind Q 96-well plate	10
92IEXQ42DD-11--D	Sartobind Q pico 0.08 mL	10
92MU0142DD-11	Sartobind Selection Kit pico 0.08 mL, Q, S and STIC PA	one each of Q, S and STIC PA

### Sartobind S







Order number	Description	Quantity
99IEXS42GC-----V	Sartobind S 96-well plate	2
99IEXS42GC-----D	Sartobind S 96-well plate	10
92IEXS42DD-11--D	Sartobind S pico 0.08 mL	10
92MU0142DD-11	Sartobind Selection Kit pico 0.08 mL, Q, S and STIC PA	one each of Q, S and STIC PA

Sartobind 96-well plate. Data sheet: order no.: 85032-543-28  
Sartobind pico 0.08 mL. Data sheet: order no.: 85032-541-71

## Accessories

Order number		Description	Quantity
1ZA---0004		Adapter Luer male to UNF 10-32 female, PEEK	1
1ZA0GV0003		Adapter UNF 10-32 female to sanitary 3/4", 25 mm, polyoxymethylene	2
5ZALB-0002		Distribution adapter for 3 × 200 to 1,200 mL (10–30") capsules, 1 × 2", 3 × 1 1/2", sanitary, stainless steel	1
7ZAL-V0013		Reducing adapter 1 1/2" (50.5 mm) to 3/4" (25 mm), sanitary	1
7ZAL-V0010		Reducing adapter 2" (64 mm) to 1 1/2" (50.5 mm), sanitary	1
9ZGL--0102		Trolley for Jumbo 2.5 or 5 L, stainless steel	1
16290		Sartocheck® 3 Plus Integrity Tester	1
26288---FT		Sartocheck® 4 Plus Filter Integrity Tester	1



Order number		Description	Quantity
29Z-S00001		Manifold set, 1½" sanitary clamp, inlet and outlet adapter plate	2
29Z-S00003		Gamma irradiated manifold set, 1½" sanitary clamp, inlet and outlet adapter plate, ≥ 25 kGy, max.50 kGy	2
2ZGL--0005		Pilot Filter Holder	1
2ZGL--0006		Process filter holder	1
2ZGL--0007		Double process filter	1
2ZGL--0008		Drip pan for Pilot Filter holder	1
2ZGL--0015		Drip pan for Process and double Process Filter Holder	1

Chromatography as easy as filtration



◀ Sartobind App  
[www.sartorius-stedim.com/apps](http://www.sartorius-stedim.com/apps)

# Sales and Service Contacts

For further contacts, visit [www.sartorius-stedim.com](http://www.sartorius-stedim.com)

## Europe

### Germany

Sartorius Stedim Biotech GmbH  
August-Spindler-Strasse 11  
37079 Goettingen  
Phone +49.551.308.0

Sartorius Stedim Systems GmbH  
Robert-Bosch-Strasse 5-7  
34302 Guxhagen  
Phone +49.5665.407.0

### France

Sartorius Stedim FMT S.A.S.  
ZI des Paluds  
Avenue de Jouques - CS 91051  
13781 Aubagne Cedex  
Phone +33.442.845600

Sartorius Stedim France SAS  
ZI des Paluds  
Avenue de Jouques - CS 71058  
13781 Aubagne Cedex  
Phone +33.442.845600

### Austria

Sartorius Stedim Austria GmbH  
Modecenterstrasse 22  
1030 Vienna  
Phone +43.1.7965763.18

### Belgium

Sartorius Stedim Belgium N.V.  
Rue Colonel Bourg 105  
1030 Bruxelles  
Phone +32.2.756.06.80

### Hungary

Sartorius Stedim Hungária Kft.  
Kagyló u. 5  
2092 Budakeszi  
Phone +36.23.457.227

### Italy

Sartorius Stedim Italy S.r.l.  
Via dell'Antella, 76/A  
50012 Antella-Bagno a Ripoli (FI)  
Phone +39.055.63.40.41

### Netherlands

Sartorius Stedim Netherlands B.V.  
Phone +31.30.60.25.080  
[filtratie.nederland@sartorius-stedim.com](mailto:filtratie.nederland@sartorius-stedim.com)

### Poland

Sartorius Stedim Poland Sp. z o.o.  
ul. Wrzesinska 70  
62-025 Kostrzyn  
Phone +48.61.647.38.40

### Russian Federation

LLC "Sartorius Stedim RUS"  
Vasilyevsky Island  
5<sup>th</sup> line 70, Lit. A  
199178 St. Petersburg  
Phone +7.812.327.53.27

### Spain

Sartorius Stedim Spain, S.A.U.  
Avda. de la Industria, 32  
Edificio PAYMA  
28108 Alcobendas (Madrid)  
Phone +34.913.586.098

### Switzerland

Sartorius Stedim Switzerland AG  
Ringstrasse 24 a  
8317 Tagelswangen  
Phone +41.52.354.36.36

### U.K.

Sartorius Stedim UK Ltd.  
Longmead Business Centre  
Blenheim Road, Epsom  
Surrey KT19 9 QQ  
Phone +44.1372.737159

### Ukraine

LLC "Sartorius Stedim RUS"  
Post Box 440 "B"  
01001 Kiev, Ukraine  
Phone +380.44.411.4918

## Americas

### USA

Sartorius Stedim North America Inc.  
5 Orville Drive, Suite 200  
Bohemia, NY 11716  
Toll-Free +1.800.368.7178

### Argentina

Sartorius Argentina S.A.  
Int. A. Ávalos 4251  
B1605ECS Munro  
Buenos Aires  
Phone +54.11.4721.0505

### Brazil

Sartorius do Brasil Ltda  
Avenida Senador Vergueiro 2962  
São Bernardo do Campo  
CEP 09600-000 - SP- Brasil  
Phone +55.11.4362.8900

### Mexico

Sartorius de México, S.A. de C.V.  
Libramiento Norte de Tepotzotlan s/n,  
Colonia Barrio Tlacateco,  
Municipio de Tepotzotlan,  
Estado de México,  
C.P. 54605  
Phone +52.55.5562.1102  
[leadsmex@sartorius.com](mailto:leadsmex@sartorius.com)

### Peru

Sartorius Peru S.A.C.  
Avenue Alberto del Campo 411  
Floor 12 - The Office  
15076 - San Isidro, Lima  
Phone +51.1.441 0158

## Asia | Pacific

### Australia

Sartorius Stedim Australia Pty. Ltd.  
Unit 5, 7-11 Rodeo Drive  
Dandenong South Vic 3175  
Phone +61.3.8762.1800

### China

Sartorius Stedim (Shanghai)  
Trading Co., Ltd.  
3rd Floor, North Wing, Tower 1  
No. 4560 Jinke Road  
Zhangjiang Hi-Tech Park  
Pudong District  
Shanghai 201210, P.R. China  
Phone +86.21.6878.2300

Sartorius Stedim (Shanghai)  
Trading Co., Ltd.  
Beijing Branch Office  
No. 33 Yu'an Road  
Airport Industrial Park Zone B  
Shunyi District, Beijing 101300  
Phone +86.10.8042.6501

Sartorius Stedim (Shanghai)  
Trading Co., Ltd.  
Guangzhou Branch Office  
Room 1105  
Xing Guang Ying Jing Building  
No. 119, Shui Yin Road  
Yue Xiu District, Guangzhou 510075  
Phone +86.20.3836.4193

### India

Sartorius Stedim India Pvt. Ltd.  
#69/2-69/3, NH 48, Jakkasandra  
Nelamangala Tq  
562 123 Bangalore, India  
Phone +91.80.4350.5250

### Japan

Sartorius Stedim Japan K.K.  
4th Fl., Daiwa Shinagawa North Bldg.  
8-11, Kita-Shinagawa 1-chome  
Shinagawa-ku, Tokyo, 140-0001 Japan  
Phone +81.3.4331.4300

### Malaysia

Sartorius Stedim Malaysia Sdn. Bhd.  
Lot L3-E-3B, Enterprise 4  
Technology Park Malaysia  
Bukit Jalil  
57000 Kuala Lumpur, Malaysia  
Phone +60.3.8996.0622

### Singapore

Sartorius Stedim Singapore Pte. Ltd.  
10 Science Park Rd  
The Alpha #02-13/14  
Singapore Science Park II  
Singapore 117684  
Phone +65.6872.3966

### South Korea

Sartorius Korea Biotech Co., Ltd.  
8th Floor, Solid Space B/D,  
PanGyoYeok-Ro 220, BunDang-Gu  
SeongNam-Si, GyeongGi-Do, 463-400  
Phone +82.31.622.5700



▶ [www.sartorius-stedim.com](http://www.sartorius-stedim.com)