## Virosart® CPV MidiCaps®

### Directions for use

#### 1. Summary

The purpose of this document is to provide the appropriate method for identifying, installing and using Sartorius Stedim Biotech Virosart® CPV MidiCaps®. The procedure should be strictly followed to ensure operator safety and to maintain optimum performance of Virosart® CPV MidiCaps®.

This document pertains to Sartorius Stedim Biotech Virosart® CPV filter capsules with the following part number:

Virosart® CPV MidiCaps® (box of 2):  
5455328V9–FF–V

All Virosart® CPV MidiCaps® have sanitary inlet & outlet connections.

#### 2. Labeling of Virosart® CPV MidiCaps®

The type, lot number and serial number are imprinted on the outer support of the housing of the Virosart® CPV MidiCaps®. This information is identical to the information found on the label of both the plastic bag and the cardboard box where in addition to type and lot number the order number is indicated. The direction of filtration is indicated by an arrow on the MidiCaps® housing.

#### 3. Installing Virosart® CPV MidiCaps®

The plastic bag of the Virosart® CPV MidiCaps® should be cut open at one end. Leave the open plastic bag on the filter for protection as long as possible. Virosart® CPV MidiCaps® are connected to the processing line at the inlet and outlet via triclamp (please do not forget the gaskets!) connections.

#### 4. Wetting of Virosart® CPV MidiCaps®

(In the direction of filtration)

Before every sterilization (autoclaving) and prior to every integrity test the entire system needs to be rinsed to ensure that it is particle-free. To rinse the filter element, open the vent valve on the MidiCaps® housing and fill the housing with WFI until liquid comes out of the vent valve. To ensure that all the air has been removed from the housing, keep the vent valve slightly opened for the first minute of flushing.

Then close the vent valve and rinse for >15 minutes at a pressure of 2.0 bar | 29 psi. The housing of the Virosart® CPV MidiCaps® can be drained by opening the drain valve and using compressed air or nitrogen with a maximal pressure of 0.5 bar | 7 psi in the direction of filtration. As the membrane of the Virosart® CPV is much tighter than a 0.2 µm membrane, it is recommended to drain the housing by pressure to speed up this procedure. After the housing is empty, close all the valves. After the Virosart® CPV has been flushed according to the instructions, it is thoroughly wetted.

#### 5. Sterilization

##### 5.1 Autoclaving of Virosart® CPV MidiCaps®

Sterilization temperature: up to 121°C, 30 min

The inlet and outlet as well as any open valves of the Virosart® CPV MidiCaps® housing should be wrapped with steam permeable autoclave paper or aluminium foil. Once the sterilization temperature has been attained the Virosart® CPV MidiCaps® can be autoclaved for 30 minutes. During the autoclave cycle of the Virosart® CPV MidiCaps®, make sure that connecting pieces on the Virosart® CPV MidiCaps® (e.g. made of stainless steel) are not attached too tightly and do not have a torque effect on the Virosart® CPV MidiCaps® housing which could cause deformation.

#### 6. Integrity Testing of Virosart® CPV MidiCaps®

It is important to follow the wetting procedure for Virosart® CPV to make sure that the filter is wetted out completely prior to the integrity test. The integrity test of a Virosart® CPV MidiCaps® must be done using an automated integrity test unit, preferably the Sartocheck® integrity test unit.

When using a Sartorius Stedim Biotech Sartocheck® integrity test system, the diffusion test can be performed automatically, recording the corresponding integrity test data.

To ensure reliable integrity testing, it is recommended to use 10 minutes stabilisation time and 5 minutes measuring time.

For additional information on how to perform an automatic integrity test please see our Installation and Operating Instructions for the Sartocheck® integrity test unit.

#### Integrity Test Data for Virosart® CPV MidiCaps®

Integrity test	Filter type	Test pressure	Limit
Water based diffusion test	5455328V9	4.5 bar   65.2 psi	10 ml/min

#### Remark:

All values given above were obtained at room temperature and using water as wetting agent.

#### 7. Maximum Allowable Differential Pressure of Virosart® CPV MidiCaps®

In the direction of filtration  
At 20°C max. 5.0 bar | 73 psi  
At 121°C max. 0.2 bar | 2.92 psi

In the reversed direction of filtration  
At 20°C max. 0.2 bar | 2.92 psi

#### 8. Filtration with Virosart® CPV MidiCaps®

Virus filtration with Virosart® CPV should be realized at the end of the purification process prior to the last buffer exchange | form and fill. The protein solution which is subject to the virus filtration must be free of aggregates, host cell protein & DNA to avoid premature blockage of Virosart® CPV MidiCaps®.

These contaminants should be removed either during the initial recovery process or latest by the polishing process of your target molecule. Standard prefiltration prior virus filtration employs a 0.2/0.1µm Sartopore 2 membrane filter.

To make full use of the filtration area, the housing must be vented completely free of air. To do this, briefly open the vent valve until the liquid to be filtered comes out. Then apply the desired filtration pressure. The virus filtration should be realized at 2 bar | 29 psi pressure. The flow rate achieved through the Virosart® CPV MidiCaps® during the filtration is a function of the purity of the feed stream, the target molecule size and concentration. The flow rate will decline with the progress of filtration and indicate the respective blockage of the filter.

#### 9. Changing out the Virosart® CPV MidiCaps®

The Virosart® CPV MidiCaps® must be changed, at the latest, when the maximum allowable differential pressure is reached or when the validated process volume has been processed over the filter.

#### 10. Applications Support

If required, please contact your local Sartorius Stedim Biotech representative to obtain further information on technical data, product integrity testing or general information concerning specific applications. Our application specialists from Purification Technologies will gladly assist in your filtration trial work if needed. Moreover our Validation Service and support team are always available to assist in questions concerning product and process validation, as well as training of staff and optimization of filtration processes.

#### 11. Return of used filter elements

If you wish to return used filter elements to Sartorius Stedim Biotech for inspection, it is necessary to decontaminate and sterilize them prior to shipment. This procedure has to be documented in the Return Shipment Form, available from your local Sartorius Stedim Biotech representative. This documentation is required by country law. Returned filter elements cannot be handled unless the appropriate Return Shipment Form is completed.

#### 12. Liability

Sartorius Stedim Biotech cannot assume liability if Virosart® CPV MidiCaps® are subjected to improper use. In the interest of product development we reserve the right to make changes.

**Thank you for working with Sartorius Stedim Biotech, we appreciate your business.**

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