



sartorius stedim  
biotech

# Celsius<sup>®</sup> CFT The Controlled-rate Freeze & Thaw Platform Technology



turning science into solutions

# Celsius® CFT – A well-controlled, scalable Freeze & Thaw Platform for Drug Substance Transfers in cGMP Manufacturing

The Celsius® CFT technology is a controlled-rate freezing and thawing system designed to minimize cryoconcentration in biopharmaceutical products. This robust and reliable platform ensures the bulk drug substance is well protected from start to finish, and gives flexibility by processing up to 100 L of solution in single-use Celsius®-Paks. A qualified shipping solution is offered as part of the Celsius® CFT platform.



Celsius® FT100 System



### ▶ Maximum Product Security

- Celsius®-Paks are secured in structural frames for robust handling – even at low temperatures – as proven through extensive functional qualification.
- Control of the entire process logistics from filling to dispensing, including frozen transport, to easily and safely manage large frozen volume.
- Fast, reproducible thermal cycling preserves product stability and mitigates the major adverse effects of large scale freezing and thawing of biologics.

### ▶ Process Reproducibility & Reliability

- Total integration with monitoring of freeze & thaw process critical quality attributes.
- Process control with CFR21 Part.11 & cGMP compliant system for production-scale equipment.
- Reproducible freeze and thaw cycles from Celsius®-Pak to Celsius®-Pak, from production batch to another production batch and from Celsius® F-T module to another Celsius® F-T module.

### ▶ Process Flexibility

- Decouple process steps to provide higher manufacturing flexibility.
- Increase capacities: extend product shelf-life, stockpiling and enable batch processing, easy transfer between DS & DP sites, facilitate management of multi-product facilities.
- Optimize overall process time with fast freeze and thaw processes (5 hours to freeze 100 L of product down to -60°C, everytime).

### ▶ Scalable from 30 mL to 100 L

- Celsius® CFT is the only fully scalable single-use technology for freeze and thaw on the market.
- Allows R&D & Process Development activities to explore & properly designed future production-scale processes.
- Celsius® S<sup>3</sup>, the scale-down system, reproduces the freeze & thaw rates observed at large scale and is a unique tool to run stability studies & formulation development.



8.3 L and 16.6 L Celsius®-Pak

# Product Security to Maximize the Product Safety

## Robustness of the Celsius®-Pak at low temperature (down to -80°C)

Because Celsius® CFT is mainly used in the highly critical application of long-term storage of the Bulk Drug Substance (BDS) and its frozen transfer or shipment to the drug product manufacturing site, an extensive functional qualification has been run to guarantee the robustness of the Celsius®-Pak under real-process conditions.

Over 100 freeze & thaw cycles using the Celsius® CFT equipment were performed to assess the design performances. The results are summarized in the Validation Guide for Celsius® bags (section 6) including connection, freeze & thaw (real and extreme conditions), liquid | frozen robustness and shipping tests.

## Manage large scale freeze & thaw processes with optimized logistics

Celsius® CFT provides safe and easy-to-use solution through a series of accessories designed to minimize manual handling of Celsius®-Pak and to mitigate the associated risks.

This gives to the end-users the possibility to build an industrial process from filling and freezing to shipping, thawing and dispensing.

## Fast freeze & thaw processes to minimize bulk-scale freeze concentration and preserve product stability

The adverse effects of slow freezing and thawing rates and its related cryoconcentration phenomenon on biologics (proteins, mAbs and cells) have been well described in the industry resulting to product loss (activity | recovery).


Please contact your local FMT Application Specialist to get access to the relevant literature highlighting the benefits of the Celsius® CFT to maintain product quality with concrete studies.

"Slow freezing rates can result in cryoconcentration. This can further lead to pH shifts and phase separation among those components resulting in protein structural damage"

Puri & Al., BioProcess International 13(1), January 2015

"Controlled-rate technologies such as Celsius® CFT prove to be much more efficient than the conventional uncontrolled-rate methods ... and also show minimal to negligible cryoconcentration"

PDA J Pharm Sci and Tech 2010, 64 290-298



## ► Product Security

- Controlled-rate Freezing and Thawing to defeat cryo-concentration, minimize product variation and prevent from product degradation
- Proven bag (same film, bag construction and quality as Flexboy®)
- 100% Leak test on bag and immediate connections
- Robust protective frames to hold securely frozen DS stored in bags with protection over tubing | connectors
- Over 150 000 Celsius®-Paks sold with proven performances in commercial drug manufacturing processes
- Complete integrated platform from filling through draining
- Filled Celsius®-Paks are never carried or directly handled throughout lifecycle

# Process Reproducibility and Reliability for Consistent Operations, Worldwide




## Reproducibility & monitoring of the COAs: Integration of the Freeze & Thaw process step into GMP operation

The entire automatic freeze, thaw and hold operations are processed through the Control Cabinet. Connected to the Temperature Control Unit and the Celsius® F-T module, the critical parameters of the different process steps are commanded, monitored and stored in an automatic manner with minimal operator intervention.



Inside to the freezing chamber of the Celsius® F-T module, the Celsius®-Pak will be held between two heat transfer plates for bi-directional and maximized heat transfer performances. Controlled-rate freezing and thawing will take place independently from different process conditions providing highly reproducible operation. Well defined, freezing and thawing processes in Celsius® CFT will be characterized at the most critical point inside to the Celsius®-Pak, the last point to freeze.



HMI panel touch screen

## ► Process Reproducibility and Reliability

- Celsius® CFT is the only fully integrated platform (Hardware & Bags) available on the market to transfer frozen DS worldwide coming from the same supplier
- Accountable platform from design to supply, performances and service
- More than 70 Production systems installed worldwide
- Excellent reproducibility from bag to bag, batch to batch and from unit to unit with strong process control for cGMP manufacturing
- Pre-qualified shippers (mechanical and thermal) with up to 100 L of BDS for International shipment



1 L and 2 L Celsius®-Pak

# Process Flexibility to Create Opportunity to Improve and Expand Manufacturing Networks

## Improve the overall manufacturing process efficiency

Besides the benefits on product stability linked to long-term frozen storage & shipment, Celsius® CFT will actively participate in maximizing productivity and flexibility of the manufacturing process. Celsius® CFT will decrease overall production costs by facilitating Drug Substance & Drug Product decoupling (one-way logistics possible) and by enabling batch processing.

Risk of microbial contamination is also mitigated as Celsius® CFT will make possible large scale freeze and thaw operation an aseptic process.



## Celsius® SSM Shipper

The Celsius® SSM shipper is a reusable shipper and allows shipment of multiple frozen Celsius®-Pak to remote location.

The SSM Shipper is qualified according to the ASTM D4169 Assurance Level 1 and the ISTA 7D extreme summer temperature profile.



Celsius® SSM Shipper

## Celsius® Single-Unit Shipper

The Celsius® Shipper is a single-use shipper and allows shipment of individual frozen Celsius®-Pak to remote locations.

The shipper is qualified according to the ASTM D4169 Assurance Level 1 and extreme summer temperature profile.



Celsius® Single-Unit Shipper

## ► Process Flexibility

- Extends Shelf- Life
- De-couple process steps, stockpile, and stop bioprocessing clock
- Improve production management and create opportunity

The fundamental purpose for Celsius® success:

- Efficient Process
- Transfers using the same system with same performances worldwide
- Leverage prior developments with benchtop data and production scale-up

Passed  
ASTM D4169 AL1  
Shipping  
Test ✓

Passed  
ISTA7D  
Shipping  
Test ✓



**Full Control of Process Logistics –**  
Optimized Accessories from Fill to Drain  
**Safe and Robust Handling –**  
Qualified Shippers and Shipping Methods


# Scalability with 30 mL to Create a 100 L Process



## Celsius® S<sup>3</sup> laboratory-scale system, the ideal tool to explore production-scale freeze & thaw operations

When selecting scale-down freeze thaw model, different parameters such as the material of construction of the container and its interaction | compatibility with the Drug Substance or the storage temperature must be based on the most representative solution to be used at large process scale. Freeze & thaw rates are known to affect product stability and cannot be scalable if using uncontrolled-rate freeze & thaw system.

Celsius® CFT is the only single-use technology on the market which provides full scalability from laboratory to commercial production scale. Risk of microbial contamination is also mitigated as Celsius® CFT will make possible large scale freeze and thaw operation an aseptic process.



## ► Scalability

- Benchtop System for quick, low-cost feasibility with 30 mL product to verify compatibility to Celsius® CFT F&T cycles
- Celsius®-Pak Range from 1 L to 16.6 L
- F-T Module is Modular (FT33 for development, add modules for clinical | commercial mfg, up to 100 L per cycle)
- Same Film: scale-up from Celsius® FFT to CFT platforms

# Sales and Service Contacts

For further contacts, visit [www.sartorius-stedim.com](http://www.sartorius-stedim.com)

## Europe

### Germany

Sartorius Stedim Biotech GmbH  
August-Spindler-Strasse 11  
37079 Goettingen

Phone +49.551.308.0  
Fax +49.551.308.3289

Sartorius Stedim Systems GmbH  
Robert-Bosch-Strasse 5 – 7  
34302 Guxhagen

Phone +49.5665.407.0  
Fax +49.5665.407.2200

### France

Sartorius Stedim FMT S.A.S.  
ZI des Paluds  
Avenue de Jouques – CS 91051  
13781 Aubagne Cedex

Phone +33.442.845600  
Fax +33.442.845619

Sartorius Stedim France SAS  
ZI des Paluds  
Avenue de Jouques – CS 71058  
13781 Aubagne Cedex

Phone +33.442.845600  
Fax +33.442.846545

### Austria

Sartorius Stedim Austria GmbH  
Modcenterstrasse 22  
1030 Vienna

Phone +43.1.7965763.18  
Fax +43.1.796576344

### Belgium

Sartorius Stedim Belgium N.V.  
Rue Colonel Bourg 105  
1030 Bruxelles

Phone +32.2.756.06.80  
Fax +32.2.756.06.81

### Hungary

Sartorius Stedim Hungária Kft.  
Kagyuló u. 5  
2092 Budakeszi

Phone +36.23.457.227  
Fax +36.23.457.147

### Italy

Sartorius Stedim Italy S.r.l.  
Via dell'Antella, 76/A  
50012 Antella-Bagno a Ripoli (FI)

Phone +39.055.63.40.41  
Fax +39.055.63.40.526

### Netherlands

Sartorius Stedim Netherlands B.V.

Phone +31.30.60.25.080  
Fax +31.30.60.25.099

[filtratie.nederland@sartorius-stedim.com](mailto:filtratie.nederland@sartorius-stedim.com)

### Poland

Sartorius Stedim Poland Sp. z o.o.  
ul. Wrzesinska 70  
62-025 Kostrzyn

Phone +48.61.647.38.40  
Fax +48.61.879.25.04

### Russian Federation

LLC "Sartorius Stedim RUS"  
Uralskaya str. 4, Lit. B  
199155 St. Petersburg

Phone +7.812.327.53.27  
Fax +7.812.327.53.23

### Spain

Sartorius Stedim Spain, S.A.U.  
Avda. de la Industria, 32  
Edificio PAYMA  
28108 Alcobendas (Madrid)

Phone +34.913.586.098  
Fax +34.913.589.623

### Switzerland

Sartorius Stedim Switzerland AG  
Ringstrasse 24 a  
8317 Tagelswangen

Phone +41.52.354.36.36  
Fax +41.52.354.36.46

### U.K.

Sartorius Stedim UK Ltd.  
Longmead Business Centre  
Blenheim Road, Epsom  
Surrey KT19 9 QQ

Phone +44.1372.737159  
Fax +44.1372.726171

### Ukraine

LLC "Sartorius Stedim RUS"  
Post Box 440 "B"  
01001 Kiev, Ukraine

Phone +380.44.411.4918  
Fax +380.50.623.3162

## Americas

### USA

Sartorius Stedim North America Inc.  
5 Orville Drive, Suite 200  
Bohemia, NY 11716

Toll-Free +1.800.368.7178  
Fax +1.631.254.4253

### Argentina

Sartorius Argentina S.A.  
Int. A. Ávalos 4251  
B1605ECS Munro  
Buenos Aires

Phone +54.11.4721.0505  
Fax +54.11.4762.2333

### Brazil

Sartorius do Brasil Ltda  
Avenida Senador Vergueiro 2962  
São Bernardo do Campo  
CEP 09600-000 - SP- Brasil

Phone +55.11.4362.8900  
Fax +55.11.4362.8901

### Mexico

Sartorius de México, S.A. de C.V.  
Libramiento Norte de Tepetzotlan s/n,  
Colonia Barrio Tlacateco,  
Municipio de Tepetzotlan,  
Estado de México,  
C.P. 54605

Phone +52.55.5562.1102  
Fax +52.55.5562.2942

[leadsmex@sartorius.com](mailto:leadsmex@sartorius.com)

### Peru

Sartorius Peru S.A.C.  
Av. Emilio Cavenecia 264 San Isidro  
15073 Lima, Perú

Phone +51.1.441 0158  
Fax +51.1.422 6100

## Asia | Pacific

### Australia

Sartorius Stedim Australia Pty. Ltd.  
Unit 5, 7-11 Rodeo Drive  
Dandenong South Vic 3175

Phone +61.3.8762.1800  
Fax +61.3.8762.1828

### China

Sartorius Stedim Biotech (Beijing) Co. Ltd.  
No. 33 Yu'an Road  
Airport Industrial Park Zone B  
Shunyi District, Beijing 101300

Phone +86.10.80426516  
Fax +86.10.80426580

Sartorius Stedim (Shanghai)  
Trading Co., Ltd.  
3rd Floor, North Wing, Tower 1  
No. 4560 Jinke Road  
Zhangjiang Hi-Tech Park  
Pudong District  
Shanghai 201210, P.R. China

Phone +86.21.6878.2300  
Fax +86.21.6878.2882

Sartorius Stedim Biotech (Beijing) Co. Ltd.  
Guangzhou Representative Office  
Unit K, Building 23  
Huihua Commerce & Trade Building  
No. 80 Xianlie Middle Road  
Guangzhou 510070

Phone +86.20.37618687 | 37618651  
Fax +86.20.37619051

### India

Sartorius Stedim India Pvt. Ltd.  
#69/2-69/3, NH 48, Jakkasandra  
Nelamangala Tq  
562 123 Bangalore, India

Phone +91.80.4350.5250  
Fax +91.80.4350.5253

### Japan

Sartorius Stedim Japan K.K.  
4th Fl., Daiwa Shinagawa North Bldg.  
8-11, Kita-Shinagawa 1-chome  
Shinagawa-ku, Tokyo, 140-0001 Japan

Phone +81.3.4331.4300  
Fax +81.3.4331.4301

### Malaysia

Sartorius Stedim Malaysia Sdn. Bhd.  
Lot L3-E-3B, Enterprise 4  
Technology Park Malaysia  
Bukit Jalil  
57000 Kuala Lumpur, Malaysia

Phone +60.3.8996.0622  
Fax +60.3.8996.0755

### Singapore

Sartorius Stedim Singapore Pte. Ltd.  
1 Science Park Road,  
The Capricorn, #05-08A,  
Singapore Science Park II  
Singapore 117528

Phone +65.6872.3966  
Fax +65.6778.2494

### South Korea

Sartorius Korea Biotech Co., Ltd.  
8th Floor, Solid Space B/D,  
PanGyoYeok-Ro 220, BunDang-Gu  
SeongNam-Si, GyeongGi-Do, 463-400

Phone +82.31.622.5700  
Fax +82.31.622.5799



▶ [www.sartorius-stedim.com](http://www.sartorius-stedim.com)