



Cultivation of Sf9 cells in the single use bioreactor BIOSTAT[®] RM 20 | 50 optical system



Application
Note

#04

#05

#06

#07

#08

Anil Kumar Rathod, Sheokant Diwakar,
Dr. Alexander Tappe*, Jacqueline Herrmann*,
Dr. Thorsten Adams*, Dr. Ashok Mundrigi,
Sartorius Stedim Biotech, Bangalore, India.

* Sartorius Stedim Biotech GmbH,
Goettingen, Germany

Introduction

Insect cells like mammalian cells are able to produce eukaryotic proteins (i.e. possess the ability of post-translational modifications, protein folding etc). The most commonly used insect cell line is Sf9 (*Spodoptera frugiperda* ovary cells). These cells can be grown in suspension culture or adherent as monolayer. A variety of recombinant proteins can be obtained by utilizing the baculovirus expression system.

The BIOSTAT® RM optical (figure 1) is the most advanced single use bioreactor using rocking motion technology. The pillow shaped cultivation chamber is rocked back and forth, creating waves which provide mixing with low shear. The liquid surface is constantly renewed, thereby enabling efficient mass transfer between head space and media. The cultivation chamber itself is a single use bag composed of a multilayer film with pharmaceutical grade ethyl vinyl acetate (EVA) as the contact layer. This single use bag is easy to use, reduces validation costs, and removes the need for cleaning, sterilizing, and provides stress free convenient culturing.

In this application note we tried to demonstrate the use of single use Cultibag RM system in cultivation of Sf9 cell line.

Material & Methods


Spodoptera frugiperda (Sf9) cell line was procured from Invitrogen and cultivated in Serum free sf-900-II SFM media (GIBCO, 10902). The Cultibag RM 2 L Optical bag (DB0002L) was filled with 800 ml of media under aseptic condition and the bag was inflated with air. Thereafter DO and pH sensors were inserted in respective ports and equilibrated for 2 hrs at 8 rocks per minutes. DO and pH sensors were calibrated as per manufacturer instructions. The pre culture was seeded in the Cultibag at seeding density of 0.7×10^6 cells/ml. Temperature was set to 28°C, DO to 60%, pH to 6.2. Rocking was varied from 18 to 25 rocks/minutes and an angle of 6°. 5 ml of sample was collected every day to determine viable cell density by trypan blue exclusion method and glucose and lactate estimation in the spent media. (Enzyme based assay kit from Chema diagnostic kit for Glucose assay and Radiant kit for Lactate assay was used).

Results

The maximum cell density obtained was 11.87×10^6 cells/ml at 9th day of cultivation (average of three consecutive trials). Cell viability decreased from 9th day of cultivation and dropped to 50% till day 12. Glucose was consumed from 10.2 gm/l to 4 gm/l during cultivation and lactate accumulation was 3 mg/ml at the end of cultivation. Temperature was maintained at 28°C constantly and pH around 6.2. Further the process was scaled up and cultivation was done in 5 L media in 10 L Cultibag RM. The peak cell density was 12.9×10^6 cells/ml similar to 2 L Cultibag RM.




Figure 1:
BIOSTAT® RM optical single use bioreactor



Graph 1:
Growth curve of Sf9 cell line in 2 L and 10 L Cultibag RM

Discussion

DO was maintained at 60% throughout the cultivation period. The variation in DO was $\pm 5\%$. DO was maintained by using nitrogen, gas flow and oxygen in cascade mode. DO was maintained primarily by air. Oxygen was required to maintain DO after the cell reaches the exponential log phase. The evaporation of media can also be curtailed by intermittent supply of oxygen instead of continuous supply of air. The pH was maintained at 6.2 by using 0.1 N alkalis (NaHCO_3) and CO_2 in cascade mode. CO_2 is not necessary at the start of the cultivation but can be used towards the end of cultivation when pH tends to rise as a result of toxic accumulation.



Graph 2:
Glucose and lactate profile of Sf9 cell line in 2 L and 10 L Cultibag RM

Rocking Speed


The rocking speed is dependent on both the culture volume and cell density. At low volumes i.e. 20–30% of the Cultibag working volume, an initial rocks of 15 is sufficient for mixing. It would be necessary to increase the number of rocks to 18–25 at 100% of Cultibag working volume or at high cell density. It is advisable to monitor the oxygen demand and adjust the rocks and rocking angle as needed. Rocking speed of 18 to 23 was suitable for cultivation of Sf9 in 2 L and 10 L Cultibag RM systems.

Rocking Angle


An initial angle of 6 degrees is sufficient. In general the angle is increased as the oxygen demand increases. When the Cultibag is at 100% of working volume, an angle of 6–10 degrees may be needed. An angle of 6° was suitable for cultivation in 2 L and 10 L Cultibag RM. The angle of rocking can be reduced when excessive foaming is observed.

Aeration Rate

DO is controlled by using nitrogen, air and oxygen in cascade mode. The controller for air flow was set maximum of 30% and minimum of 5%, corresponding to flow rates between 50 and 300 ml/min. Nitrogen and oxygen controller limits were set to 100%. The required set point of 60% is maintained throughout the cultivation run.



Graph 3:
History plot of Cultivation of Sf9 cell line in 2 L Cultibag RM



Graph 4:
History plot of Cultivation of Sf9 cell line in 10 L Cultibag RM

Conclusion

In this note, we have demonstrated that the single use bioreactor BIOSTAT® Cultibag 20 RM is suitable for the cultivation of Sf9 cell lines in serum free chemically defined media. Insect cell lines have higher oxygen demand compared to mammalian cell lines and Cultibag RM is able to support their growth. Single use bioreactors are preferred in today's cutting edge applications in R&D, process development and small scale production. Single use bioreactors remove the requirement for cleaning validation and require minimum maintenance while providing maximum operator ease of use.

Sales and Service Contacts

For further contacts, visit www.sartorius-stedim.com

Europe

Germany

Sartorius Stedim Biotech GmbH
August-Spindler-Strasse 11
37079 Goettingen

Phone +49.551.308.0
Fax +49.551.308.3289

www.sartorius-stedim.com

Sartorius Stedim Systems GmbH
Robert-Bosch-Strasse 5-7
34302 Guxhagen

Phone +49.5665.407.0
Fax +49.5665.407.2200

France

Sartorius Stedim Biotech S.A.
ZI Les Paluds
Avenue de Jouques - BP 1051
13781 Aubagne Cedex

Phone +33.442.845600
Fax +33.442.845619

Sartorius Stedim France SAS
ZI Les Paluds
Avenue de Jouques - CS 71058
13781 Aubagne Cedex

Phone +33.442.845600
Fax +33.442.846545

Austria

Sartorius Stedim Austria GmbH
Franzosengraben 12
1030 Vienna

Phone +43.1.7965763.18
Fax +43.1.796576344

Belgium

Sartorius Stedim Belgium N.V.
Leuvensesteenweg, 248/B
1800 Vilvoorde

Phone +32.2.756.06.80
Fax +32.2.756.06.81

Denmark

Sartorius Stedim Nordic A/S
Hoerskaetten 6D, 1.
2630 Taastrup

Phone +45.7023.4400
Fax +45.4630.4030

Hungary

Sartorius Stedim Hungária Kft
Kagyló u. 5
2092 Budakeszi

Phone +36.23.457.227
Fax +36.23.457.147

Italy

Sartorius Stedim Italy S.p.A.
Via dell'Antella, 76/A
50012 Antella-Bagno a Ripoli (FI)

Phone +39.055.63.40.41
Fax +39.055.63.40.526

Netherlands

Sartorius Stedim Netherlands B.V.
Edisonbaan 24
3439 MN Nieuwegein

Phone +31.30.6025080
Fax +31.30.6025099

Poland

Sartorius Stedim Poland Sp. z o.o.
ul. Wrzesinska 70
62-025 Kostrzyn

Phone +48.61.647.38.40
Fax +48.61.879.25.04

Russian Federation

000 „Sartorius ICR“ and 000 „Biohit“
Uralskaya str. 4, Lit. B
199155, Saint-Petersburg

Phone +7.812.327.5.327
Fax +7.812.327.5.323

Spain

Sartorius Stedim Spain SA
C/Isabel Colbrand 10,
Oficina 70
Polígono Industrial de Fuencarral
28050 Madrid

Phone +34.90.2110935
Fax +34.91.3589623

Switzerland

Sartorius Stedim Switzerland AG
Ringstr. 24 a
8317 Tagelswangen

Phone +41.52.354.36.36
Fax +41.52.354.36.46

U.K.

Sartorius Stedim UK Limited
Longmead Business Park
Blenheim Road, Epsom
Surrey KT19 9 QQ

Phone +44.1372.737159
Fax +44.1372.726171

America

USA

Sartorius Stedim North America Inc.
5 Orville Drive, Suite 200
Bohemia, NY 11716

Toll-Free +1.800.368.7178
Fax +1.631.254.4253

Argentina

Sartorius Argentina S.A.
Int. A. Avalos 4251
B1605ECS Munro
Buenos Aires

Phone +54.11.4721.0505
Fax +54.11.4762.2333

Brazil

Sartorius do Brasil Ltda
Av. Dom Pedro I, 241
Bairro Vila Pires
Santo André
São Paulo
Cep 09110-001

Phone +55.11.4451.6226
Fax +55.11.4451.4369

Mexico

Sartorius de México S.A. de C.V.
Circuito Circunvalación Poniente No. 149
Ciudad Satélite
53100, Estado de México
México

Phone +52.5555.62.1102
Fax +52.5555.62.2942

Asia | Pacific

Australia

Sartorius Stedim Australia Pty. Ltd.
Unit 5, 7-11 Rodeo Drive
Dandenong South Vic 3175

Phone +61.3.8762.1800
Fax +61.3.8762.1828

China

Sartorius Stedim Biotech (Beijing)
Co. Limited
Airport Industrial Zone B
No. 33 Yu'an Road
Beijing 101300, Shunyi District

Phone +86.10.80426516
Fax +86.10.80426580

Sartorius Stedim Biotech (Beijing)
Co. Limited
Shanghai Branch office
Room 618, Tower 1, German Centre,
Shanghai, PRC., 201203

Phone +86.21.28986393
Fax +86.21.28986392.11

Sartorius Stedim Biotech (Beijing)
Co. Limited
Guangzhou representative office
Room 704, Broadway Plaza,
No. 233-234 Dong Feng West Road
Guangzhou 510180

Phone +86.20.8351.7921
Fax +86.20.8351.7931

India

Sartorius Stedim India Pvt. Ltd.
#69/2-69/3, Jakkasandra
Kunigal Road, Nelamangala Tq
Bangalore - 562 123

Phone +91.80.4350.5248
Fax +91.80.4350.5253

Japan

Sartorius Stedim Japan K.K.
Kiba Park Bldg
5-11-13 Kiba
Koto-ku
Tokyo 135-0042

Phone +81.3.5639.9981
Fax +81.3.5639.9983

Malaysia

Sartorius Stedim Malaysia Sdn. Bhd.
Lot L3-E-3B, Enterprise 4
Technology Park Malaysia
Bukit Jalil
57000 Kuala Lumpur, Malaysia

Phone +60.3.8996.0622
Fax +60.3.8996.0755

Singapore

Sartorius Stedim Singapore Pte. Ltd.
1 Science Park Road,
The Capricorn, #05-08A,
Singapore Science Park 2
Singapore 117528

Phone +65.6872.3966
Fax +65.6778.2494

South Korea

Sartorius Korea Biotech Co., Ltd.
8th Floor, Solid Space B/D,
PanGyoYeok-Ro 220, BunDang-Gu
SeongNam-Si, GyeongGi-Do, 463-400

Phone +82.31.622.5700
Fax +82.31.622.5799



◀▶ www.sartorius-stedim.com