




Detecting Airborne Influenza Virus A



Application
Note

#10

#11

Collecting Viruses

#12

#13

#14

Detecting Airborne Influenza Virus A (H3N2) in a Children's Polyclinic in Greifswald, Germany, by collecting Virus Using the Water-soluble 12602 Sartorius Stedim Biotech Gelatin Membrane Filter

Introduction

Although the impinger and impaction collectors based on the impaction principle have long ago proved their reliability for the collection of microbial aerosols, filters only recently began to receive increased attention with the development of the water-soluble gelation filter technique and the wide spectrum of applications offered by such filters only for collection of bacterial aerosols; however - after Goettingen-based Sartorius Stedim Biotech GmbH had started commercial-scale manufacture of such gelatin filters.

Until now, the value of filters for detection and collection of virus aerosols has been considered limited, and the filters were thought to be useful for collecting stable, resistant viruses, if at all. There has been very little experimental data about this until recently.

Studies done by the Jaschhof team with experimentally produced static aerosols of T1 and T3 coli phages and influenza viruses, strain A|PR|8|34 (H1N1) using the Sartorius Stedim Biotech Gelatin Membrane Filter (12602) yielded a constant collection efficiency for these viruses, even under conditions involving high inlet velocities (up to 1.6 m/s max.) and an extended collection period (15 minutes max.).

The high resistance of the gelatin filter to mechanical stress and the proven stability

of the influenza virus were the reasons for testing the effectiveness of this method in practice during the seasonal rise in morbidity for acute respiratory diseases (ARD) at the beginning of 1990 (3rd - 5th week). This study builds upon results obtained from tests with experimentally produced virus and phage aerosols.

Air samples were taken over a period of 3 days in the waiting room of a Greifswald children's polyclinic and the water-soluble filter material was then used in virus cultivating experiments. Experiments for cultivating the virus isolated from nasal secretions obtained by aspiration from the children in the waiting room, who had a fresh case of ARD, were performed in parallel.


Materials and Methods

Location Studied:

Waiting room of a Greifswald children's clinic (approx 120 m³) frequented by a maximum of 57 children accompanied by a parent over a period of 2 hours (Jan. 24, 1990). Temperature: 25°C, 60% relative humidity.

Experimental | Test Materials:

Nasal secretions from children with a fresh case of ARD. Gelatin membrane filter, 12602, after 1.8 m³ of room air was filtered through it. It was then sealed in a polyethylene bag along with 5 ml of an isotonic buffer for cell cultivation at pH 7.4 with



Scanning electron micograph of the influenza virus A|PR|8|34 (H0N1). Viruses from a culture of an 11-day old embryonic chicken egg and purified with a saccharose density gradient.

Micrograph: Herrmann | Fischer

0.5% yeast extract (Difco), penicillin G (400 IU/ml) and streptomycin (400 µg/ml). All samples were processed within 4 hrs (with storage in between in a cooling chamber).

Air sampling:

3 days of air filtration with each sample consisting of 1.8 m³ of air per filter (120 l/min) during a collection period of 15 min through a Sartorius Stedim Biotech Gelatin Filter (126 02-050).

Virological Analyses:

Nasal secretions: Hank's solution (3 ml per sample) was added and the samples were centrifuged for 30 min (3,000 rpm), each given an additive of 10 µl antibiotic solution/ml of supernatant, and incubated for 1 hr at room temperature. The amnion of 10-day old incubated eggs was inoculated with these specimens and incubated for 3 days at 33°C; 3 blind subcultures were prepared; HAT was used to detect the viruses.

Air samples: The filter was liquified for 5 min at 37°C in a water bath; an antibiotic solution as described above was added to the solution which was used to inoculate incubated eggs; FL cell and ascites tumor cells were then cultured; the viruses were detected according to the usual techniques by HAT, CPE, and HAAdT.

The influenza virus isolated from the nasal secretions and air samples were typed and identified by cross titration in HAHT with rabbit immune sera against the influenza virus A|Greifswald|1|89 (H1N1) and A|Greifswald|1|90 (H3N2) and antigens from the influenza virus A|Sichuan|2|87 (H3N2) and influenza virus A|Greifswald|2|88 (H3N2).

The v-RNA and virus protein patterns were determined by PAGE.

Results and Discussion

From the experimental material taken on Jan. 24, 1990, three types of viruses isolated from the nasal secretions and an influenza virus A isolated from an air sample were successfully inoculated in the incubated embryonated eggs during the third transfer for subculturing (influenza virus A|Greifswald|2|90 [H3N2]).

HAT showed that the antigens of the isolated viruses and those of the variants A|Greifswald|2|88 and A|Sichuan|2|87 were largely related and did not exhibit any differences from a strain which was cultivated from a nasal secretion specimen taken from a patient who was in the room at the time of air sampling using filters (influenza virus A|Greifswald|1|90 [H3N2]).

The migration rates of the isolated viruses' proteins (HA, NP, NS and M) and all bands of their v-RNA following PAGE (polyacrylamide gel electrophoresis) matched to a high degree (Fig. 1). Attempts to isolate further viruses from air samples using cell cultures were not successful.

The fact that the air samples were taken in the pre-epidemic phase of the ARD morbidity dynamics demonstrates the effectiveness of the filter collection method used in detecting airborne influenza viruses.

Seventy-five percent of the children (≥ 8 years old) present on the day of the studies had a respiratory infection, but only 40% of these children at most showed acute catarrhal symptoms (coughing, sneezing, runny nose). The detection of the influenza virus A in the air of a children's clinic waiting room, during a period of only a slight increase in seasonal morbidity from acute respiratory diseases in the region, remarkably confirms* the results obtained from using experimentally produced static aerosols of coli phages and influenza virus A.

This verifies, for the first time the suitability of the gelatin filter in practice for the collection of respiratory viruses from the air.

* (see Sartorius Stedim Biotech Application Notes Publication No. SLF4028-e)

Outlook

Considering the uncomplicated method of air filtration involving comparatively less work, this principle using the gelatin membrane filter is bound to receive increased recognition in helping to answer a multitude of questions about the production, stability and spread of infectious virus aerosols (e.g., involved in genetic engineering techniques).

Supplementary issue of the condensed text version of a poster displayed during the Annual Conference of the Society of Virology, Potsdam March 19-23, 1991).

Names and addresses of the authors:

Dr. Helmut Jaschhof
Ph. D. in Natural Sciences
Gerhard-Katsch-Str. 14
17489 Greifswald, Germany

Dr. E.-J. Finke, M.D.
Bundeswehrkrankenhaus Berlin
Scharnhorststr. 13
10115 Berlin 1040, Germany

Dr. B. Herrmann
Ph. D. in Natural Sciences
und A. Reinecke
Institut für Medizinische Mikrobiologie
und Epidemiologie (Institute of Medical
Microbiology and Epidemiology)
E.-M.-Arndt-Universität Greifswald,
Martin-Luther-Straße 6,
17489 Greifswald, Germany




Figure 1: Separation of the v-RNA in polyacrylamide gel electrophoresis (PAGE) from the influenza virus A|Greifswald|1|90 (H3N2) (isolated from patient specimens) and the influenza virus A|Greifswald|2|90 (collected and isolated from the air.) 3.5% polyacrylamide gel in tris-acetate EDTA buffer.

Sales and Service Contacts

For further contacts, visit www.sartorius-stedim.com

Europe

Germany

Sartorius Stedim Biotech GmbH
August-Spindler-Strasse 11
37079 Goettingen

Phone +49.551.308.0
Fax +49.551.308.3289

www.sartorius-stedim.com

Sartorius Stedim Systems GmbH
Schwarzenberger Weg 73-79
34212 Melsungen

Phone +49.5661.71.3400
Fax +49.5661.71.3702

www.sartorius-stedim.com

France

Sartorius Stedim Biotech S.A.
ZI Les Paluds
Avenue de Jouques – BP 1051
13781 Aubagne Cedex

Phone +33.442.845600
Fax +33.442.845619

Sartorius Stedim France SAS
ZI Les Paluds
Avenue de Jouques – CS 71058
13781 Aubagne Cedex

Phone +33.442.845600
Fax +33.442.846545

Austria

Sartorius Stedim Austria GmbH
Franzosengraben 12
A-1030 Vienna

Phone +43.1.7965763.18
Fax +43.1.796576344

Belgium

Sartorius Stedim Belgium N.V.
Leuvensesteenweg, 248/B
1800 Vilvoorde

Phone +32.2.756.06.80
Fax +32.2.756.06.81

Denmark

Sartorius Stedim Nordic A/S
Hoerskaetten 6D, 1.
DK-2630 Taastrup

Phone +45.7023.4400
Fax +45.4630.4030

Italy

Sartorius Stedim Italy S.p.A.
Via dell'Antella, 76/A
50012 Antella-Bagno a Ripoli (FI)

Phone +39.055.63.40.41
Fax +39.055.63.40.526

Netherlands

Sartorius Stedim Netherlands B.V.
Edisonbaan 24
3439 MN Nieuwegein

Phone +31.30.6025080
Fax +31.30.6025099

Spain

Sartorius Stedim Spain SA
C/Isabel Colbrand 10-12,
Planta 4, Oficina 121
Polígono Industrial de Fuencarral
28050 Madrid

Phone +34.91.3586102
Fax +34.91.3588804

Switzerland

Sartorius Stedim Switzerland GmbH
Lerzenstrasse 21
8953 Dietikon

Phone +41.44.741.05.00
Fax +41.44.741.05.09

U.K.

Sartorius Stedim UK Limited
Longmead Business Park
Blenheim Road, Epsom
Surrey KT19 9 QQ

Phone +44.1372.737159
Fax +44.1372.726171

America

USA

Sartorius Stedim North America Inc.
5 Orville Drive
Bohemia, NY 11716

Toll-Free +1.800.368.7178
Fax +1.631.254.4253

Sartorius Stedim SUS Inc.
1910 Mark Court
Concord, CA 94520

Phone +1.925.689.6650
Toll Free +1.800.914.6644
Fax +1.925.689.6988

Sartorius Stedim Systems Inc.
201 South Ingram Mill Road
Springfield, MO 65802

Phone +1.417.873.9636
Fax +1.417.873.9275

Argentina

Sartorius Argentina S.A.
Int. A. Avalos 4251
B1605ECS Munro
Buenos Aires

Phone +54.11.4721.0505
Fax +54.11.4762.2333

Brazil

Sartorius do Brasil Ltda
Av. Dom Pedro I, 241
Bairro Vila Pires
Santo André
São Paulo
Cep 09110-001

Phone +55.11.4451.6226
Fax +55.11.4451.4369

Mexico

Sartorius de México S.A. de C.V.
Circuito Circunvalación Poniente No. 149
Ciudad Satélite
53100 Naucalpan, Estado de México

Phone +52.5555.62.1102
Fax +52.5555.62.2942

Asia | Pacific

China

Sartorius Stedim Beijing
Representative Office
No. 33, Yu'an Road,
Airport Industrial Zone B, Shunyi District
Beijing 101300

Phone +86.10.80426516
Fax +86.10.80426580

Sartorius Stedim Shanghai
Representative Office
Room 618, Tower 1, German Centre,
Shanghai, PRC., 201203

Phone +86.21.28986393
Fax +86.21.28986392.11

Sartorius Stedim Guangzhou Office
Room 704, Broadway Plaza,
No. 233-234 Dong Feng West Road
Guangzhou 510180

Phone +86.20.8351.7921
Fax +86.20.8351.7931

India

Sartorius Stedim India Pvt. Ltd.
10, 6th Main, 3rd Phase Peenya
KIADB Industrial Area
Bangalore – 560 058

Phone +91.80.2839.1963|0461
Fax +91.80.2839.8262

Japan

Sartorius Stedim Japan K.K.
KY Building, 8-11
Kita Shinagawa 1-chome
Shinagawa-ku
Tokyo 140-0001

Phone +81.3.3740.5407
Fax +81.3.3740.5406

Malaysia

Sartorius Stedim Malaysia Sdn. Bhd.
Lot L3-E-3B, Enterprise 4
Technology Park Malaysia
Bukit Jalil
57000 Kuala Lumpur

Phone +60.3.8996.0622
Fax +60.3.8996.0755

Singapore

Sartorius Stedim Singapore Pte. Ltd.
10, Science Park Road, The Alpha
#02-25, Singapore Science Park 2
Singapore 117684

Phone +65.6872.3966
Fax +65.6778.2494

Australia

Sartorius Stedim Australia Pty. Ltd.
Unit 5, 7-11 Rodeo Drive
Dandenong South Vic 3175

Phone +61.3.8762.1800
Fax +61.3.8762.1828