



Impact of Pressure Release and Multiple Pressure Fluctuations on Virus Retention Performance of Virosart® CPV Virus Retentive Filters



Application
Note

#02

#03

#04

#05

#06

Acknowledges to Sherri Dolan and Carl Breuning, Sartorius Stedim North America

Introduction

All biopharmaceutical products derived from human or animal origins must have a proven virus safety concept to ensure the final patients' health before the product is released to commercial manufacturing or even clinical trials. The concept must cover known and unknown potential viral contaminants. Contaminating viruses can come from external sources or be present in cell lines. Endogenous retrovirus-like particles, for example, are present in most CHO cell lines. Regulators require that manufacturers perform a risk analysis and have a strategy to remove contaminating viruses during downstream processing. The industry considers filtration to be a robust method of virus removal and hence it is a widely used method. Filtration relies on the principle of size-exclusion and can remove all types of viruses. Investigators have recently discovered that under specific conditions with some virus removal membranes, virus breakthrough can occur. Scientists believe that breakthrough is more likely to occur when the membrane is operated at high capacities, during flow-decay, under conditions of high or low process pressure or when the pressure is held and released.

Virus filtration with Virosart® virus retentive filters is an integral part of the orthogonal virus clearance technology platform of Sartorius Stedim Biotech. This orthogonal technology platform features virus clearance by filtration (size exclusion), inactivation and adsorption. The Virosart® product range includes four different virus retentive membranes, in order to provide the best solution for every application.


Objective

The objective of this study is to determine the impact of pressure fluctuation on the retention of parvovirus using bacteriophage PP7, an established model system for parvoviruses. Parvoviruses are renowned for being small viruses that are typically difficult to remove.

The well established Virosart® CPV virus retentive filter membrane, from Sartorius, was tested. Virosart® CPV is a 20 nm double layer PES membrane, which provides high flow rates and superior capacity. This filter can be steam sterilized and is ideal for the application within stainless steel systems.

Materials

- Compressed nitrogen together with the Sartorius constant pressure test system with stainless steel sample reservoir
- Virosart® CPV Minisart® 5.0 cm²
- Sartorius Balances
- Virosart® Max Minisart®, 5 cm²
- Buffer: 50 mM Tris, 10 mM NaCl, pH 8.1
- Product: mAb; c = 5 g/L



Picture 1: Set-up of Virosart® CPV Minisarts®




Method

A virus spike of PP7 with infectivity titer $>7 \log/\text{mL}$ was added to a monoclonal antibody solution. The mAb solution, $c=5 \text{ g/L}$, had been pre-filtered offline using a 5 cm^2 Virosart® Max, a $0.1 \mu\text{m}$ adsorptive pre-filter for Virosart® virus retentive filters, prior to being spiked. A pressure control device for monitoring and varying the trans-membrane pressure during the filtration experiments was connected to the pressure vessel containing the spiked feed solution.


The Virosart® CPV Minisart® was loaded to a capacity of 140 L/m^2 that was equivalent to 70 ml of mAb solution. The filtrate volume was collected in a container and measured automatically using a balance connected to data recording software.

Three studies were performed as follows:


1. Multi-step pressure reduction



2. Rapid pressure changes, low pressure and pressure hold



3. Pressure release study




Results

Study 1: Multi-step pressure reduction


Duplicate filtration runs were performed with two 5 cm^2 , Virosart® CPV Minisart®. Flow was started by pressurizing the reservoirs to $50 \text{ psi} | 3.44 \text{ bar}$ and an initial 10 ml fraction taken.

10 ml fractions were collected at each of the following pressure steps ($40 \text{ psi} | 2.75 \text{ bar}$, $32 \text{ psi} | 2.2 \text{ bar}$, $24 \text{ psi} | 1.65 \text{ bar}$, $16 \text{ psi} | 1.1 \text{ bar}$ and $8 \text{ psi} | 0.55 \text{ bar}$). The filtration was then paused for 10 minutes and the pressure returned to $40 \text{ psi} | 2.75 \text{ bar}$ before the final fraction was collected. Graph 1 shows the pressure profile during the experiment. Graph 2 shows the corresponding LRV profile of the fractions.

	Pressure (psi)	LRV (PP7) Virosart® CPV Filter 1	LRV (PP7) Virosart® CPV Filter 2
Load		7.67	7.67
Fraction 1	50	≥ 7.67	≥ 7.67
Fraction 2	40	≥ 7.67	≥ 7.67
Fraction 3	32	≥ 7.67	≥ 7.67
Fraction 4	24	≥ 7.67	≥ 7.67
Fraction 5	16	7.37	≥ 7.67
Fraction 6	8	≥ 7.67	≥ 7.67
10 minute pause	0		
Fraction 7	40	≥ 7.67	≥ 7.67



Graph 1: Pressure course during multi-step pressure reduction



Graph 2: Multi-step pressure reduction onto Virosart® CPV showing constant independent retention performance.

Summary of Study 1:


All fractions showed stable retention above $>7 \log$, allowing the conclusion that multi-step pressure reduction had no effect on the retention of Virosart® CPV Minisart®.

Study 2: Rapid pressure change, low pressure and pressure release


Duplicate filtration runs were performed with two 5 cm², Virosart[®] CPV Minisart[®]. The reservoirs were pressurized to 40 psi | 2.75 bar and the first fraction (20 ml) was taken. While the second 20 ml fraction was collected the pressure was rapidly switched off and on 5 times (fraction 2). Afterwards the filtration was allowed to proceed at a low pressure of only 8 psi | 0.55 bar during which the third 20 ml fraction was collected.

The filtration was then paused for 10 minutes and the pressure returned to 40 psi | 2.75 bar before the final fraction was collected. Graph 3 shows the pressure profile during the experiment. Graph 4 shows the corresponding LRV profile of the fractions.

	Pressure (psi)	LRV (PP7) Virosart [®] CPV Filter 1	LRV (PP7) Virosart [®] CPV Filter 2
Load		7.52	7.52
Fraction 1	40	≥ 7.52	≥ 7.52
Fraction 2	0 – 40 cycle 5 x	≥ 7.52	≥ 7.52
Fraction 3	8	≥ 7.52	≥ 7.52
10 minute pause	0		
Fraction 4	40	≥ 7.52	≥ 7.52



Graph 3: Pressure course during rapid pressure change, low pressure and pressure release



Graph 4: Rapid pressure change, low pressure and pressure release onto Virosart[®] CPV showing complete retention performance.


Summary of Study 2:

No fractions of the filtrate showed any virus breakthrough allowing the conclusion that the rapid pressure changes had no effect on the retention of Virosart[®] CPV Minisart[®].


Study 3: Pressure release study

Duplicate filtration runs were performed with a 5 cm², Virosart® CPV Minisart®. Flow was started by pressurizing the reservoir to 40 psi | 2.75 bar. Upon reaching 50% of the total loading (fraction1: 70 L/m²) the system was de-pressurized and held a 0 psi | 0 bar for 10 minutes. Following the de-pressurized hold, the system was re-pressurized to 40 psi | 2.75 bar and the remaining material collected as fraction 2 (total loading capacity 140 L/m²). Finally the system was de-pressurized, wash buffer added to the reservoir, and the wash fraction was collected. Graph 5 shows the pressure profile during the experiment. Graph 6 shows the corresponding LRV profile of the fractions.

	Pressure (psi)	LRV (PP7) Virosart® CPV Filter 1	LRV (PP7) Virosart® CPV Filter 2
Load		7.38	7.38
Fraction 1: (Load 0 – 70 L/m ²)	40	≥ 7.38	≥ 7.38
10 minute pause	0		
Fraction 2: (Load 70 – 140 L/m ²)	40	≥ 7.38	≥ 7.38
Fraction 3: Buffer flush	40	≥ 7.38	≥ 7.38



Graph 5: Pressure course during pressure release study



Graph 6: Pressure release onto Virosart® CPV showing complete retention performance.

Summary of Study 3:

No fractions of the filtrate showed any virus breakthrough allowing the conclusion that the performed pressure release during virus filtration had no effect on the retention of Virosart® CPV Minisart®.


Conclusions

This application note summarizes three filtration studies investigating the effect of pressure variations on the retention of PP7 bacteriophage by Virosart® CPV Minisarts®.

The virus retention of the Virosart® CPV Minisarts® was unaffected by pressure changes in any of the studies.

It can be concluded that the Virosart® CPV Minisarts®, provide robust retention under the conditions studied in the experiments described in this application note.

Although this membrane offers high virus retention under a wide range of processing conditions, users are recommended to evaluate the impact of pressure excursions onto membrane performance under their specific process conditions.





Sales and Service Contacts

For further contacts, visit www.sartorius-stedim.com

Europe

Germany

Sartorius Stedim Biotech GmbH
August-Spindler-Strasse 11
37079 Goettingen

Phone +49.551.308.0
Fax +49.551.308.3289

Sartorius Stedim Systems GmbH
Robert-Bosch-Strasse 5 – 7
34302 Guxhagen

Phone +49.5665.407.0
Fax +49.5665.407.2200

France

Sartorius Stedim FMT S.A.S.
ZI des Paluds
Avenue de Jouques – CS 91051
13781 Aubagne Cedex

Phone +33.442.845600
Fax +33.442.845619

Sartorius Stedim France SAS
ZI des Paluds
Avenue de Jouques – CS 71058
13781 Aubagne Cedex

Phone +33.442.845600
Fax +33.442.846545

Austria

Sartorius Stedim Austria GmbH
Modecenterstrasse 22
1030 Vienna

Phone +43.1.7965763.18
Fax +43.1.796576344

Belgium

Sartorius Stedim Belgium N.V.
Rue Colonel Bourg 105
1030 Bruxelles

Phone +32.2.756.06.80
Fax +32.2.756.06.81

Hungary

Sartorius Stedim Hungária Kft.
Kagyuló u. 5
2092 Budakeszi

Phone +36.23.457.227
Fax +36.23.457.147

Italy

Sartorius Stedim Italy S.r.l.
Via dell'Antella, 76/A
50012 Antella-Bagno a Ripoli (FI)

Phone +39.055.63.40.41
Fax +39.055.63.40.526

Netherlands

Sartorius Stedim Netherlands B.V.

Phone +31.30.60.25.080
Fax +31.30.60.25.099

filtratie.nederland@sartorius-stedim.com

Poland

Sartorius Stedim Poland Sp. z o.o.
ul. Wrzesinska 70
62-025 Kostrzyn

Phone +48.61.647.38.40
Fax +48.61.879.25.04

Russian Federation

LLC "Sartorius Stedim RUS"
Uralskaya str. 4, Lit. B
199155 St. Petersburg

Phone +7.812.327.53.27
Fax +7.812.327.53.23

Spain

Sartorius Stedim Spain, S.A.U.
Avda. de la Industria, 32
Edificio PAYMA
28108 Alcobendas (Madrid)

Phone +34.913.586.098
Fax +34.913.589.623

Switzerland

Sartorius Stedim Switzerland AG
Ringstrasse 24 a
8317 Tagelswangen

Phone +41.52.354.36.36
Fax +41.52.354.36.46

U.K.

Sartorius Stedim UK Ltd.
Longmead Business Centre
Blenheim Road, Epsom
Surrey KT19 9 QQ

Phone +44.1372.737159
Fax +44.1372.726171

Ukraine

LLC "Sartorius Stedim RUS"
Post Box 440 "B"
01001 Kiev, Ukraine

Phone +380.44.411.4918
Fax +380.50.623.3162

Americas

USA

Sartorius Stedim North America Inc.
5 Orville Drive, Suite 200
Bohemia, NY 11716

Toll-Free +1.800.368.7178
Fax +1.631.254.4253

Argentina

Sartorius Argentina S.A.
Int. A. Ávalos 4251
B1605ECS Munro
Buenos Aires

Phone +54.11.4721.0505
Fax +54.11.4762.2333

Brazil

Sartorius do Brasil Ltda
Avenida Senador Vergueiro 2962
São Bernardo do Campo
CEP 09600-000 - SP- Brasil

Phone +55.11.4362.8900
Fax +55.11.4362.8901

Mexico

Sartorius de México, S.A. de C.V.
Libramiento Norte de Tepotzotlan s/n,
Colonia Barrio Tlacateco,
Municipio de Tepotzotlan,
Estado de México,
C.P. 54605

Phone +52.55.5562.1102
Fax +52.55.5562.2942

leadsmex@sartorius.com

Peru

Sartorius Peru S.A.C.
Avenue Alberto del Campo 411
Floor 12 – The Office
15076 – San Isidro, Lima

Phone +51.1.441 0158
Fax +51.1.422 6100

Asia | Pacific

Australia

Sartorius Stedim Australia Pty. Ltd.
Unit 5, 7-11 Rodeo Drive
Dandenong South Vic 3175

Phone +61.3.8762.1800
Fax +61.3.8762.1828

China

Sartorius Stedim Biotech (Beijing) Co. Ltd.
No. 33 Yu'an Road
Airport Industrial Park Zone B
Shunyi District, Beijing 101300

Phone +86.10.80426516
Fax +86.10.80426580

Sartorius Stedim (Shanghai)
Trading Co., Ltd.

3rd Floor, North Wing, Tower 1
No. 4560 Jinke Road
Zhangjiang Hi-Tech Park
Pudong District
Shanghai 201210, P.R. China

Phone +86.21.6878.2300
Fax +86.21.6878.2882

Sartorius Stedim Biotech (Beijing) Co. Ltd.
Guangzhou Representative Office

Unit K, Building 23
Huihua Commerce & Trade Building
No. 80 Xianlie Middle Road
Guangzhou 510070

Phone +86.20.37618687 | 37618651
Fax +86.20.37619051

India

Sartorius Stedim India Pvt. Ltd.
#69/2-69/3, NH 48, Jakkasandra
Nelamangala Tq
562 123 Bangalore, India

Phone +91.80.4350.5250
Fax +91.80.4350.5253

Japan

Sartorius Stedim Japan K.K.
4th Fl., Daiwa Shinagawa North Bldg.
8-11, Kita-Shinagawa 1-chome
Shinagawa-ku, Tokyo, 140-0001 Japan

Phone +81.3.4331.4300
Fax +81.3.4331.4301

Malaysia

Sartorius Stedim Malaysia Sdn. Bhd.
Lot L3-E-3B, Enterprise 4
Technology Park Malaysia
Bukit Jalil
57000 Kuala Lumpur, Malaysia

Phone +60.3.8996.0622
Fax +60.3.8996.0755

Singapore

Sartorius Stedim Singapore Pte. Ltd.
1 Science Park Road,
The Capricorn, #05-08A,
Singapore Science Park II
Singapore 117528

Phone +65.6872.3966
Fax +65.6778.2494

South Korea

Sartorius Korea Biotech Co., Ltd.
8th Floor, Solid Space B/D,
PanGyoYeok-Ro 220, BunDang-Gu
SeongNam-Si, GyeongGi-Do, 463-400

Phone +82.31.622.5700
Fax +82.31.622.5799



▶ www.sartorius-stedim.com