



Preventive Maintenance

BioSealer® TC
Order No.: S873SSMW

Benefits

- Consistently robust sealing of thermoplastic tubes
- Optimized equipment performance
- Trouble-free operation; reduction of downtime
- Compliance with regulatory requirements



Product Information

The objective of preventive maintenance is a continuous, precise operation of your equipment, to secure your investments for the long term and to guarantee correct and proper function.

Preventive maintenance will be performed in accordance with Sartorius' fixed standards and in between the contractually stipulated time intervals between the customer and Sartorius.

Specifications

- Inspection visit at fixed intervals or as needed
- Visual inspection
- Maintenance check and cleaning according to maintenance checklist
- Temperature calibration
- Test sealing of customer-specific tube types
- Pressure test (optional)
- Issue of maintenance|test report

Prerequisites required by the customer

- Device is standing at the maintenance location
- Unobstructed access to the equipment
- Customer-specific tubes needs to be provided by the customer
- When ordering, a responsible person has to be appointed
- Access by automobile is possible | parking lot is available

Repairs, spare parts as well as travel costs and expenses are not covered in scope of services and will be invoiced separately.

Interested? Just get in touch with us for a free service consultation. You will find all important information on our website at www.sartorius.com/service.

Sartorius Stedim Biotech GmbH
August-Spindler-Strasse 11
37079 Goettingen, Germany
Phone +49.551.308.0
www.sartorius-stedim.com

USA Toll-Free +1.800.368.7178
Argentina +54.11.4721.0505
Brazil +55.11.4362.8900
Mexico +52.55.5562.1102
UK +44.1372.737159
France +33.442.845600
Italy +39.055.63.40.41

Spain +34.913.586.098
Russian Federation
+7.812.327.53.27
Republic of South Africa
+27.11.315.5444
Japan +81.3.4331.4300
China +86.21.6878.2300



www.sartorius.com