



Flexsafe[®] RM Bags

replace CultiBag[®] RM Bags



New generation of single-use RM bags

- Assurance of Supply
- Excellent Cell Growth
- Superior Robustness

Please be notified that on 31st December 2016 CultiBag[®] RM single-use bags will be discontinued. Final orders will be accepted until 30th of September 2016.

Passed
ASTM Shipping Test ✓

Filled Flexsafe[®] shipping bags in Palletank[®] passed the most stringent shipping validation test according to ASTM D4169.

Global 
SUPPLY GUARANTEE

Complete control over our manufacturing process from resin and film extrusion to final bag assembly

turning science **into solutions**



Learn more about Flexsafe[®]
www.sartorius-stedim.com/flexsafe

It's time for change now. CultiBag® RM bags will be discontinued – last order: September 2016!



Flexsafe® RM Designs

- Basic Bags
- Basic Bags with **viable biomass sensor** **NEW**
- Optical with pH and DO sensors
- Optical with pH, DO and **viable biomass sensors** **NEW**
- Perfusion with integrated membrane
- Perfusion with connection to ATF System
- Perfusion with pH, DO and **viable biomass sensors** **NEW**
- Customized bags

Ordering Information

Bag Size*	Basic Basic with ViaMass		Basic with screw cap Optical with screw cap		Optical Optical with ViaMass	
	NEW Flexsafe® Order No.	Old CultiBag® Order No.	NEW Flexsafe® Order No.	Old CultiBag® Order No.	NEW Flexsafe® Order No.	Old CultiBag® Order No.
1 L	DFB001L	DBB001L				
2 L	DFB002L	DBB002L	DFB002L----01SC DFO002L----01SC	DBB002L----01SC DBO002L----01SC	DFO002L	DBO002L
10 L	DFB010L DFB010L----VM	DBB010L	DFB010L----01SC DFO010L----01SC	DBB010L----01SC	DFO010L DFO010L----VM	DBO010L
20 L	DFB020L DFB020L----VM	DBB020L	DFB020L----01SC DFO020L----01SC	DBB020L----01SC DBO020L----01SC	DFO020L DFO020L----VM	DBO020L
50 L	DFB050L DFB050L----VM	DBB050L	DFB050L----01SC DFO050L----01SC		DFO050L DFO050L----VM	DBO050L
50 L **	DFB050L----01US	DBB050L----01US				
100 L	DFB100L	DBB100L			DFO100L	DBO100L
200 L	DFB200L	DBB200L			DFO200L	DBO200L

Bag Size*	Perfusion Perfusion with ViaMass		Perfusion ATF Perfusion ATF with ViaMass	
	NEW Flexsafe® Order No.	Old CultiBag® Order No.	NEW Flexsafe® Order No.	Old CultiBag® Order No.
2 L	DFP002L--SM	DBP002L--SM	DFP002L--AT	
10 L	DFP010L--SM DFP010L--SMVM	DBP010L--SM	DFP010L--AT DFP010L--ATVM	DBP010L--AT
20 L	DFP020L--SM DFP020L--SMVM	DBP020L--SM	DFP020L--AT DFP020L--ATVM	DBP020L--AT
50 L	DFP050L--SM DFP050L--SMVM	DBP050L--SM	DFP050L--AT DFP050L--ATVM	
100 L			DFP100L--AT	DBP100L--AT
200 L			DFP200L--AT	

* total volume of the cultivation chamber. Max. working volume = 0.5 x total volume; Min. working volume = 0.1 x total volume; (bags with sensors might require higher min. working volumes depending on rocking rate and angle.)

** alternative dimension: 740 x 720 mm; suitable for competitor rocker



◀ www.sartorius.com/flexsafe-rm-bags

Sartorius Stedim Biotech GmbH
August-Spindler-Strasse 11
37079 Goettingen

Phone +49.551.308.0
Fax +49.551.308.3289

Specifications subject to change without notice.
Printed in the EU on paper bleached without chlorine. | W
Publication No.: SBT1023-e150902
Order No.: 85037-551-71
Ver. 09 | 2015