



## Midisart® 2000

Ready-to-use filter for sterilizing gases and venting

### Product Information

The Midisart® 2000 is the best choice for small-scale sterile air | gas filtration applications. Superior flow rates and seamless integration into single-use assemblies provide Midisart® customers with a distinct advantage that directly translates into cost savings. Individual 100% factory integrity testing verifies the quality of every filter prior to use.



### Introduction

Midisart® 2000 is designed for maximum handling ease and safety. The device is lightweight only 20 g – so it will not weigh down or kink the tubing. With a diameter measuring just 64 mm, Midisart® incorporates a filter area of 20 cm<sup>2</sup>, which means that it is “packed” with high flow rate performance power.

Midisart® 2000 filters can be integrity tested (Bubble Point Test, 60% IPA, BPmin. = 1.1 bar). The Midisart® Test Manifold (Order code 1Z-LB-0002) can be used to test 10 filters in parallel which simply leads to a reduction of test time by 90%.

### Venting Applications

- Filling and Transfer Vessels
- Fermentation Carboys
- Cell Culture Chambers
- CO<sub>2</sub> incubators (6 to 120 liters)
- Holding and Storage Tanks
- Autoclave venting
- Small Bioreactors

### We prove your Trust

The single-layer hydrophobic PTFE membrane is self-dewetting and therefore provides high air flow rates even at very low differential pressures over an extended process time. All components were developed and are produced under full control of Sartorius. Furthermore the following indispensable requirements are fulfilled:

- **100% integrity tested prior to release**
- **Clear labeling** with Lot number, individual unit number and “IN” as an indicator for direction of flow. All information lasered on the top part of each filter housing.
- **Retention of  $\geq 10^7$  Brev. diminuta/cm<sup>2</sup>** in aerosol and liquid

## Technical Data

---

### Materials

---

Membrane PTFE – reinforced with polypropylene gauze  
Housing Polypropylene

### Pore Sizes

---

0.2 µm  
0.45 µm

### Filtration Area

---

20 cm<sup>2</sup> | 3 square inch

### Max. Differential Pressure for Air

---

3 bar (in-out)  
2 bar (out-in)

### Housing diameter

---

64 mm | 2,5"

### Connectors

---

Hose Barb, 1/8" NPT, Triclamp, small Hose Barb

### Sterilization

---

Max. Temperature 134°C  
Max. Autoclaving 20 Cycles  
Not suitable for gamma irradiation  
Gamma irradiatable disc filter: Midisart® Sartopore® Air

### Integrity Testing

---

Min. Bubble Point >1100 mbar | 15.95 psi  
Wetting agent: IPA/Water (60 | 40)

### Regulatory Compliance

---

- Each individual element is tested for integrity in direction of filtration
- 0.2 µm Pore size: Fully validated as sterilizing grade filters according to ASTM current F-838 guidelines
- Designed, developed and manufactured in accordance with an ISO 9001 certified Quality Management System
- USP Plastic Class VI Test
- Non fiber releasing according to 21 CFR

### Technical References

---

#### Validation Guide

SLD5702-c

## Ordering Information

Order No.	Pore Size [ $\mu\text{m}$ ]	Membrane	Connectors	Qty.	Sterile
17804-----E	0.45	PTFE	Hose Barb   Hose Barb	12	Yes
17804-----G	0.45	PTFE	Hose Barb   Hose Barb	25	Yes
17804-----NPE	0.45	PTFE	1/8"   1/8" NPT	12	Yes
17804-----NPG	0.45	PTFE	1/8"   1/8" NPT	25	Yes
17805-----E	0.2	PTFE	Hose Barb   Hose Barb	12	Yes
17805-----G	0.2	PTFE	Hose Barb   Hose Barb	25	Yes
17805-----NPE	0.2	PTFE	1/8"   1/8" NPT	12	Yes
17805-----NPG	0.2	PTFE	1/8"   1/8" NPT	25	Yes
17805-----UPN	0.2	PTFE	Hose Barb   Hose Barb	100	No
17805-----UPQ	0.2	PTFE	Hose Barb   Hose Barb	500	No
17805-----UQN	0.2	PTFE	Hose Barb   Hose Barb	100*	No
17809-----UNN	0.2	PTFE	1/8"   1/8" NPT	100	No
17812-----UNN	0.2	PTFE	1/8"   Hose Barb	100	No
17805-----TCN	0.2	PTFE	TriClamp   TriClamp	100	No
17877-----UPN	0.2	PTFE	Small Hose Barb   Small Hose Barb	100	No

\* = plastic bucket – the bucket with its smooth surface can be easily transferred into clean-room areas without risking particle release from the packaging material.



Standard Hose Barb



Small Hose Barb



1/8" NPT Thread



TriClamp

# Sales and Service Contacts

For further contacts, visit [www.sartorius-stedim.com](http://www.sartorius-stedim.com)

## Europe

### Germany

Sartorius Stedim Biotech GmbH  
August-Spindler-Strasse 11  
37079 Goettingen

Phone +49.551.308.0

Sartorius Stedim Systems GmbH  
Robert-Bosch-Strasse 5-7  
34302 Guxhagen

Phone +49.5665.407.0

### France

Sartorius Stedim FMT S.A.S.  
ZI des Paluds  
Avenue de Jouques - CS 91051  
13781 Aubagne Cedex

Phone +33.442.845600

Sartorius Stedim France SAS  
ZI des Paluds  
Avenue de Jouques - CS 71058  
13781 Aubagne Cedex

Phone +33.442.845600

### Austria

Sartorius Stedim Austria GmbH  
Modectcenterstrasse 22  
1030 Vienna

Phone +43.1.7965763.18

### Belgium

Sartorius Stedim Belgium N.V.  
Rue Colonel Bourg 105  
1030 Bruxelles

Phone +32.2.756.06.80

### Hungary

Sartorius Stedim Hungária Kft.  
Kagyló u. 5  
2092 Budakeszi

Phone +36.23.457.227

### Italy

Sartorius Stedim Italy S.r.l.  
via A. Meucci 4,  
50012 Grassano - Bagno a Ripoli (FI)

Phone +39.055.63.40.41

### Netherlands

Sartorius Stedim Netherlands B.V.

Phone +31.30.60.25.080

[filtratie.nederland@sartorius-stedim.com](mailto:filtratie.nederland@sartorius-stedim.com)

### Poland

Sartorius Stedim Poland Sp. z o.o.  
ul. Wrzesinska 70  
62-025 Kostrzyn

Phone +48.61.647.38.40

### Russian Federation

LLC "Sartorius Stedim RUS"  
Vasilyevsky Island  
5<sup>th</sup> line 70, Lit. A  
199178 St. Petersburg

Phone +7.812.327.53.27

### Spain

Sartorius Stedim Spain, S.A.U.  
Avda. de la Industria, 32  
Edificio PAYMA  
28108 Alcobendas (Madrid)

Phone +34.913.586.098

### Switzerland

Sartorius Stedim Switzerland AG  
Ringstrasse 24 a  
8317 Tagelswangen

Phone +41.52.354.36.36

### U.K.

Sartorius Stedim UK Ltd.  
Longmead Business Centre  
Blenheim Road, Epsom  
Surrey KT19 9 QQ

Phone +44.1372.737159

### Ukraine

LLC "Sartorius Stedim RUS"  
Post Box 440 "B"  
01001 Kiev, Ukraine

Phone +380.44.411.4918

## Africa

### Republic of South Africa

Sartorius South Africa (Pty) Ltd  
Unit 4, Alphen Square South  
853 16th Road  
Midrand 1685, RSA

Phone +27.11.315.5444

## Americas

### USA

Sartorius Stedim North America Inc.  
565 Johnson Avenue  
Bohemia, NY 11716

Toll-Free +1.800.368.7178

### Argentina

Sartorius Argentina S.A.  
Int. A. Ávalos 4251  
B1605ECS Munro  
Buenos Aires

Phone +54.11.4721.0505

### Brazil

Sartorius do Brasil Ltda  
Avenida Senador Vergueiro 2962  
São Bernardo do Campo  
CEP 09600-000 - SP- Brasil

Phone +55.11.4362.8900

### Mexico

Sartorius de México, S.A. de C.V.  
Libramiento Norte de Tepotzotlan s/n,  
Colonia Barrio Tlacateco,  
Municipio de Tepotzotlan,  
Estado de México,  
C.P. 54605

Phone +52.55.5562.1102

[leadsmex@sartorius.com](mailto:leadsmex@sartorius.com)

## Asia | Pacific

### Australia

Sartorius Stedim Australia Pty. Ltd.  
Unit 5, 7-11 Rodeo Drive  
Dandenong South Vic 3175

Phone +61.3.8762.1800

### China

Sartorius Stedim (Shanghai)  
Trading Co., Ltd.  
3rd Floor, North Wing, Tower 1  
No. 4560 Jinke Road  
Zhangjiang Hi-Tech Park  
Pudong District  
Shanghai 201210, P.R. China

Phone +86.21.6878.2300

Sartorius Stedim (Shanghai)  
Trading Co., Ltd.  
Beijing Branch Office  
No. 33 Yu'an Road  
Airport Industrial Park Zone B  
Shunyi District, Beijing 101300

Phone +86.10.8042.6501

Sartorius Stedim (Shanghai)  
Trading Co., Ltd.  
Guangzhou Branch Office  
Room 1105  
Xing Guang Ying Jing Building  
No. 119, Shui Yin Road  
Yue Xiu District, Guangzhou 510075

Phone +86.20.3836.4193

### India

Sartorius Stedim India Pvt. Ltd.  
#69/2-69/3, NH 48, Jakkasandra  
Nelamangala Tq  
562 123 Bangalore, India

Phone +91.80.4350.5250

### Japan

Sartorius Stedim Japan K.K.  
4th Fl., Daiwa Shinagawa North Bldg.  
8-11, Kita-Shinagawa 1-chome  
Shinagawa-ku, Tokyo, 140-0001 Japan

Phone +81.3.4331.4300

### Malaysia

Sartorius Stedim Malaysia Sdn. Bhd.  
Lot L3-E-3B, Enterprise 4  
Technology Park Malaysia  
Bukit Jalil  
57000 Kuala Lumpur, Malaysia

Phone +60.3.8996.0622

### Singapore

Sartorius Stedim Singapore Pte. Ltd.  
10 Science Park Rd  
The Alpha #02-13/14  
Singapore Science Park II  
Singapore 117684

Phone +65.6872.3966

### South Korea

Sartorius Korea Biotech Co., Ltd.  
8th Floor, Solid Space B/D,  
PanGyoYeok-Ro 220, BunDang-Gu  
SeongNam-Si, GyeongGi-Do, 463-400

Phone +82.31.622.5700



▶ [www.sartorius-stedim.com](http://www.sartorius-stedim.com)