



Installation

Configurable Single-use Systems

Benefits

- Value retention and enhancement of your system through professional service
- Smooth commissioning of your system after installation
- Friendly support in training operating personnel to work with your system



Product Information

The reliability and lifetime of your equipment strongly depend on whether they have been installed and configured according to specifications. Through professional installation and commissioning of your FlexAct® system, you will ensure that your equipment will be up and running perfectly right from the start. Sartorius Service will install and configure your equipment in compliance with regulatory requirements and quality standards.

Specifications

- Completeness check of the equipment supplied
- Inspection of the equipment for any signs of damage due to transportation
- Inspection of the manufacturer's type label for correct ratings
- Installation of system components, including the FlexAct® control unit and system components (pump, filter holder, etc.)
- Startup tests
- Brief training of operating personnel related to the system (not training in applications)
- Technical advice on equipment maintenance and maintenance cycles

The installation is available for the following single-use systems:

- FlexAct® systems for the following:
 - Preparation of buffers and media
 - Cell harvest | Depth filtration | DBF (clarification by body feed filtration)
 - Virus inactivation (pH reduction)
 - Virus filtration
 - Ultrafiltration | Diafiltration

Remarks

Installation and commissioning of further system components are not part of the scope of our services and must be ordered separately.

Customer Prerequisites

- The equipment is already at the installation site
- Sufficient space for access to the equipment is available
- Power supply is available at the place of installation
- Relatively large systems must be unpacked in the presence of a Sartorius technician and set up at the place of installation
- A helper must be available to assist with lifting and carrying the equipment
- The customer must designate the person responsible at the customer's site when placing an order
- Supply lines are available at the place of installation

Interested? Just get in touch with us for a free service consultation. You will find all important information on our website at www.sartorius.com/service.



◀ www.sartorius.com

Sartorius Stedim Biotech GmbH
August-Spindler-Strasse 11
37079 Goettingen, Germany
Phone +49.551.308.0
Fax +49.551.308.3289
www.sartorius-stedim.com

USA Toll-Free +1.800.368.7178
Argentina +54.11.4721.0505
Brazil +55.11.4451.6226
Mexico +52.5555.62.1102
UK +44.1372.737159
France +33.442.845600
Italy +39.055.63.40.41
Spain +34.913.586.098
Russian Federation +7.812.327.5.327
Japan +81.3.4331.4300
China +86.21.68782300

Specifications subject to change without notice.
Copyright Sartorius Stedim Biotech GmbH. Printed
in the EU on paper bleached without chlorine.
Publication No.: S--2219-e160901
Order No.: 85037-556-71
Ver. 09 | 2016