



Device Qualification IQ | OQ

BioSealer® TC
Order No.: S873SI00Q

Benefits

- 100% traceable documentation right from the beginning
- All required certificates provided
- Comprehensive IQ|OQ documentation
- Compliance with the requirements of GLP | GMP standards
- Strict Sartorius test standards
- Specially trained service technicians will support you during your product qualification



Product Information

Prerequisite for validation of your processes is qualification of the devices used. This is why qualification of measurement systems is obligatory in many areas of the pharmaceutical industry.

Sartorius supports you in GLP/GMPcompliant device qualification with an expert team of trained specialists. Devicespecific reports and documentation can easily be integrated in existing QM systems.

Specifications

- Device qualification and documentation are based on device-specific qualification documents
- Check and documentation of scope of delivery and provided documents
- Check of device for transport damage and documentation
- Installation and commissioning of device
- Documentation of device installation
- Functional check acc. to specifications, incl. documentation of results
- Brief instructions to operating personnel at workstation (not application training) with documentation

Prerequisites required by the customer

- Device is standing at the installation location
- Unobstructed access to the device
- Power is available at the installation location
- Tubing material needs to be provided
- When ordering, a responsible person has to be appointed

Interested? Just get in touch with us for a free service consultation. You will find all important information on our website at www.sartorius.com/service.

Sartorius Stedim Biotech GmbH
August-Spindler-Strasse 11
37079 Goettingen, Germany
Phone +49.551.308.0
www.sartorius-stedim.com

USA Toll-Free +1.800.368.7178
Argentina +54.11.4721.0505
Brazil +55.11.4362.8900
Mexico +52.55.5562.1102
UK +44.1372.737159
France +33.442.845600
Italy +39.055.63.40.41

Spain +34.913.586.098
Russian Federation
+7.812.327.53.27
Republic of South Africa
+27.11.315.5444
Japan +81.3.4331.4300
China +86.21.6878.2300



www.sartorius.com