

mLINE® Mechanical Pipette Multipacks Convenient Combo-Pack for Your Lab



mLINE® Mechanical Pipette

- Low pipetting forces
- Ergonomic finger support minimizes the grip force needed to hold the pipette
- Volume adjustment locking for preventing accidental volume changes
- Easy-to-read volume display
- Colour-coding of volume ranges for easy selection of corresponding pipette tips
- Easy user adjustment e.g. for different liquids
- Simple to clean and maintain with only three parts to disassemble
- Fully autoclavable

Sartorius Optifit Tips

- High-quality standard tips
- Excellent solution for various applications
- Versatile selection of packages and purity options

Mechanical Pipette and Tip Multipacks

Multipack Order Code	mLINE® Pipettes	Sartorius Optifit Tips	Accessories
mLINE® Pipette 3-pack 10 LH-725661	1-channel - 0,5-10 µl - 10-100 µl - 100-1 000 µl	Tip Racks - 96 tips - 0,1-10 µl - 0,5-200 µl - 10-1 000 µl	Linear Stand
mLINE® Pipette 3-pack 20 LH-725662	1-channel - 2-20 µl - 20-200 µl - 100-1 000 µl	Tip Racks - 96 tips - 0,5-200 µl - 10-1 000 µl	Linear Stand
mLINE® Pipette 4-pack LH-725663	1-channel - 0,5-10 µl - 10-100 µl - 20-200 µl - 100-1 000 µl	Tip Racks - 96 tips - 0,1-10 µl - 0,5-200 µl - 10-1 000 µl	Linear Stand
mLINE® Pipette 5-pack LH-725664	1-channel - 2-20 µl - 10-100 µl - 20-200 µl - 100-1 000 µl - 500-5 000 µl	Tip Racks - 96 tips - 0,5-200 µl - 10-1 000 µl Tip Rack - 50 tips - 100-5 000 µl	Linear Stand
mLINE® Pipette 3+1-pack LH-725665	1-channel - 0,5-10 µl - 20-200 µl - 100-1 000 µl 8-channel - 30-300 µl	Tip Racks - 96 tips - 0,1-10 µl - 0,5-200 µl - 5-350 µl - 10-1 000 µl	Linear Stand



Contact Details

Sartorius Biohit Liquid Handling Oy
Laippatie 1
00880 Helsinki
Phone +358.9.755.951
Fax +358.9.755.95.292
lhinfo.finland@sartorius.com

Sartorius Lab Instruments GmbH & Co. KG
Otto-Brenner-Strasse 20
37079 Goettingen
Phone +49.551.308.0
Fax +49.551.308.3289
www.sartorius.com