

Installation

Integrity Testers

Benefits

- Reproducible test results
- Calibrated sensors
- Trouble-free and reliable integrity testing
- Multiunits for simultaneous performance of several integrity tests
- Considerable time savings



Product Information

The reliability and lifetime of your instruments strongly depend on whether they have been installed and configured according to specifications. Through professional installation and commissioning of your Sartocheck® integrity tester, you will ensure that your equipment will be up and running perfectly right from the start. Sartorius Service will install and configure your equipment in compliance with regulatory requirements and quality standards.

Customer Prerequisites

Specifications

- Completeness check of the equipment supplied
 - Inspection of the equipment for any signs of damage due to transportation
 - Inspection of the manufacturer's type label for correct ratings
 - Installation of system components
 - Integrity test
 - Startup test
 - Brief training of operating personnel related to the system (not training in applications)
 - Technical advice on equipment maintenance and maintenance cycles
 - Only for the Sartochek® Multiunit:
 - Interfacing of the Multiunit to Sartochek® 4 plus leakage test
- The equipment is already at the installation site
 - Sufficient space for access to the equipment is available
 - Supply lines are available at the place of installation
 - Compressed air supply of at least 7 bar is available
 - The customer must designate the person responsible at the customer's site when placing an order
 - Only for Sartochek® Multiunit:
 - A compatible Sartochek® 4 plus must be provided by the customer
 - Only for Sartochek® Trolley:
 - A corresponding Sartochek® Filter Unit must be available
 - A connection to a water supply line must be available

The installation is available for the following systems:

- Sartochek® mini
- Sartochek® 3 plus | 4 plus
- Sartochek® Multiunit
- Sartochek® Trolley

Sales and Service Contacts

For further contacts, visit
www.sartorius.com

Germany

Sartorius Stedim Biotech GmbH
August-Spindler-Strasse 11
37079 Goettingen
Phone +49 551 308 0

USA

Sartorius Stedim North America Inc.
565 Johnson Avenue
Bohemia, NY 11716
Toll-Free +1 800 368 7178