

YDP50

Dot Matrix Printer
for Sartorius Balances



Description

- Acquisition and documentation of current measurement data
- Compact and robust design
- High printing speed

Scope of Delivery

- 1x YDP50 Dot Matrix Printer
- 1x Power Supply
- 10x Country-specific power plug adapter
- 1x USB connection cable (USB-B to USB-C)
- 2x Paper rolls

Technical Specifications


YDP50	
Compatible Balances	Entris II, Quintix Pro
Interfaces	USB-B, type 2.0, RS232, D-Sub, 9-pin
Dimensions (L x W x H)	188.2 x 114.6 x 86.6
Consumables	YDP50-Paper (10 Paper Rolls) 6906918 (Ink Ribbon Cartridge)

Germany

Sartorius Lab Instruments GmbH & Co. KG
Otto-Brenner-Strasse 20
37079 Goettingen
Phone +49 551 308 0

USA

Sartorius Corporation
565 Johnson Avenue
Bohemia, NY 11716
Phone +1 631 254 4249
Toll-free +1 800 635 2906

 For further contacts, visit
www.sartorius.com